DISRUPTOR GENIE™ Series

DISRUPTOR GENIE™

Digital DISRUPTOR GENIE™

Instruction Manual

EN

ES

FR

DE

IT

Change History

#	Date	Version	Descriptions
1	2026.02.24	A	New release

Table of Contents

Introduction	2
Applicable Products	2
Safety Information	2
Package Contents	3
Intended Use	4
Installation	4
Specifications	4
Control Panel	
<i>DISRUPTOR GENIE</i>	5
<i>Digital DISRUPTOR GENIE</i>	5
Applications	6
Operation Instruction	
<i>DISRUPTOR GENIE</i>	7
<i>Digital DISRUPTOR GENIE</i>	7
Maintenance	9
Technical Support Information	9
Technical Data	9
Compliance	10

Introduction

This manual contains installation, operation, and maintenance instructions for OHAUS DISRUPTOR GENIE™ cell disruptors. Please read the manual completely before using.

Applicable Products

This instruction manual is applicable to the following products:

Product	Model Name
DISRUPTOR GENIE	e-D2 (SI-xxxx*)
Digital DISRUPTOR GENIE	e-D2D (SI-xxxx*)

* SI-xxxx represents the original SI model name.

Safety Information

Safety notes are marked with signal words and warning symbols. These show safety issues and warnings. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results.

WARNING For a hazardous situation with medium risk, possibly resulting in severe injuries or death if not avoided.

CAUTION For a hazardous situation with low risk, resulting in damage to the device or the property or in loss of data, or minor or medium injuries if not avoided.

ATTENTION For important information about the product. May lead to equipment damage if not avoided.

NOTE For useful information about the product.

Warning Symbols

General hazard

Electric shock hazard

Safety Precautions

CAUTION: Read all safety warnings before installing, making connections, or servicing this equipment. Failure to comply with these warnings could result in personal injury and/or property damage. Retain all instructions for future reference.

Always operate unit on a level surface for best performance and maximum safety.

DO NOT lift unit by its attachment.

Disconnect unit from the power supply prior to maintenance and servicing.

Spills should be removed promptly.

DO NOT immerse the unit for cleaning.

DO NOT operate the unit if it shows signs of electrical or mechanical damage.

DO NOT position the equipment such that it is difficult to disconnect the power cord during use. Always place the unit on a clean and sturdy work surface.

WARNING! The protection provided by the unit may be impaired if used with accessories not provided or recommended by the manufacturer or used in a manner not specified by the manufacturer.

WARNING! DO NOT use the unit in a hazardous atmosphere or with hazardous materials for which the unit was not designed.

CAUTION: To avoid electrical shock, completely cut off power to the unit by disconnecting the power cord from the unit or unplugging from the wall outlet.

Earth Ground – Protective Conductor Terminal

Alternating Current

Package Contents

- Cell Disruptor
- Cup Head

We recommend you retain the original packaging for 90 days in case you need to return the product for any reason to your distributor or Ohaus.

Intended Use

This instrument is intended for use in laboratories, pharmacies, schools, businesses and light industry. It must only be used for processing materials as described in these operating instructions. Any other type of use and operation beyond the limits of technical specifications, without written consent from OHAUS, is considered as not intended. This instrument complies with current industry standards and the recognized safety regulations; however, it can constitute a hazard in use. If the instrument is not used according to these operating instructions, the intended protection provided by the instrument may be impaired.

Installation

Upon receiving the unit, check to ensure that no damage has occurred during shipment. It is important that any damage that occurred during transport is detected at the time of unpacking. If you do find such damage the carrier must be notified immediately.

After unpacking, place the unit on a level bench or table, away from explosive vapors. Ensure that the surface on which the unit is placed will withstand the typical heat produced by the unit. Do not position the equipment such that it is difficult to disconnect the power cord during use. Always place the unit on a clean and sturdy work surface.

The unit is supplied with a power cord that should be plugged into a properly grounded outlet. If the supplied power cord does not meet your needs, please use an approved power cord that suits local codes and electric supply. Replacement of the cord must be made by a qualified electrician.

Specifications

Product	DISRUPTOR GENIE	Digital DISRUPTOR GENIE
Control	Analog	Digital
Display	--	Red LED
Accessories Included	Universal Microtube Holder, Pop-off Cup	
Capacity	12x 1.5/2.0 mL microtubes	
Speed Range	3000 rpm 2850 rpm (for 50Hz models)	1000 rpm - 3000 rpm 1000 rpm - 2850 rpm (for 50Hz models)
Speed Alarms	--	On/Off, Adjustable
Timer	0-15 minutes	1-99 minutes
Power	120V, 0.65A, 60Hz 230V, 0.5A, 50Hz 100V, 1.0A, 50Hz	120V, 0.65A, 60Hz 230V, 0.5A, 50Hz 100V, 1.0A, 50/60Hz
Dimensions (L x W x H)	165 x 122 x 190 mm / 6.5 x 4.8 x 7.5 inch	
Net Weight	4.3 kg / 9.5 lb.	

Control Panel

DISRUPTOR GENIE

A	Timer Dial
B	Timer Control Knob

Digital DISRUPTOR GENIE

A	Time Display
B	Speed Display
C	Speed Alarm On / Off Indicator
D	Speed Alarm Button
E	Timer Up / Down Buttons
F	Start / Stop Button
G	Speed Up / Down Buttons

Applications

The OHAUS DISRUPTOR GENIE is a vigorous mechanical cell disruption device. Applications include cell disruption of yeast, bacteria, plant and animal tissue and for cell re-suspension for DNA mini-preps. The DISRUPTOR GENIE may be used in a cold room.

The DISRUPTOR GENIE consists of a VORTEX-GENIE mixer base with a unique and patented accessory attached directly to the shaft eccentric of the mixer. The DISRUPTOR GENIE dramatically increases the disruption efficiency of the disruption beads because of the unique and patented design of the sample holder.

Twelve sample tubes are placed in the Upper Ring of the attachment and allowed to project through corresponding holes in the "Floating" Ring below.

The holes in the "Floating" Ring are smaller than the orbit path of the shaft. Therefore, as each of the twelve tubes travel in an orbit corresponding in size with the orbit path of the shaft, the bottoms of the tubes repeatedly and vigorously collide within the confines of the holes in the "Floating" Ring. This action, in combination with the action of the disruption beads used inside the sample tubes, results in a significant increase in the impact frequency and force. The disruption efficiency is increased dramatically.

Selection of bead size and material is important. The following guide is for reference only. It is strongly suggested that you use this guide only as a starting point to determine the optimum conditions for the particular cell disruption operation that you wish to perform. The reference to **Typical Run Time** should result in a yield of 85 - 90% release of the cellular material for each of the conditions noted. Normally, 1 part of glass or zirconia (ZrO₂) beads to 2 parts of total liquid biomass is the minimum condition recommended for effective cell disruption. You can use more beads i.e. 4 parts of beads to 3 parts of biomass providing that adequate agitation of the slurry is possible. Generally, the more beads that are used, the faster the disruption of cells.

Type of Cell	Bead Size	Bead Material	Typical Run Time
Bacteria/ Spores	0.1mm	Borosilicate glass	2 to 3 min / full speed
		Zirconia (ZrO ₂)	1 to 2 min / full speed
Yeast, Mycelia, Microalgae Trypsinized cultures, Cyanobacteria	0.5mm	Borosilicate glass	3 to 4 min / full speed
		Zirconia (ZrO ₂)	1.5 to 2 min / full speed

During cell disruption, friction inside the sample tubes will generate heat. The DISRUPTOR GENIE is capable of being operated in a cold room or an environmental chamber. In most short

runs, 3 minutes or less, it is feasible that if you refrigerate the sample tubes, beads and sample before the disruption run, you will be able to perform the disruption at room temperature.

The small size of the microtubes may allow sufficient heat to dissipate so that sample integrity is not impaired. The operator must make the decision based on the type of beads, the sample content and size and the length of time of the run. In any event, it is advisable to keep the tubes, beads and sample chilled before each run.

Operation Instructions

DISRUPTOR GENIE

- 1.0 - Place the Disruptor Genie on a sturdy, level work surface.
- 2.0 - Plug the line cord to a properly grounded three-pronged electrical outlet
- 3.0 - Load up to twelve microtubes in the tube holder. Space them evenly if using less than twelve tubes.
- 4.0 - Close the cover, making sure that the stem of the tube holder goes into the cap in the cover.
- 5.0 - The DISRUPTOR GENIE is now ready for use. Turning the timer knob clockwise starts the mixer. It will run for up to 15 minutes, depending on how far you turn the dial. Turn the dial back counterclockwise to stop the unit before the time has expired, otherwise the DISRUPTOR GENIE will turn itself off after the amount of time set. Alternatively, turning the timer knob firmly counterclockwise from the "off" position until a click is heard, will turn on the DISRUPTOR GENIE in untimed mode. Firmly turn the timer knob clockwise back to the "off" position to stop the unit. At room temperature, do not operate the mixer for more than three (3) hours continuously.
- 6.0 - Accessory Tube Holders may be ordered separately and used on the DISRUPTOR GENIE. (Part No. 83041467 - Universal Microtube Holder and Part No. 83041730 - 1.5ml Snap-top Tube Holder) To remove and replace the Tube Holder, **first...make sure the DISRUPTOR GENIE is turned off**. Pop off the Tube Holder by inserting your fingers under the top disk and lift the Tube Holder with a prying motion. Align the flat surface of the DISRUPTOR GENIE mixing shaft with the flat surface of the shaft receptacle on the underside of the new Tube Holder. Press downward until the new Tube Holder snaps into place. The DISRUPTOR GENIE is now ready for use with the new Tube Holder.

Digital DISRUPTOR GENIE

- 1.0 - Place the Digital DISRUPTOR GENIE on a sturdy, level work surface.
- 1.1 - Connect the power cord to a properly grounded three-pronged receptacle.
- 1.2 - Load up to twelve microtubes in the tube holder. Space them evenly if using less than twelve tubes.

1.3 - Close the cover, making sure that the stem of the tube holder goes into the cap in the cover.

1.4 - The Digital DISRUPTOR GENIE is now ready for use.

MODES OF OPERATION

2.0 - The Digital DISRUPTOR GENIE has the following modes of operation:

2.1 - "TIMED"

Press the UP/DOWN arrow button beneath the **TIME** window to set the desired cycle time from 1 minute to 99 minutes. Next, set the speed control by pressing the UP/ DOWN arrow beneath the **SPEED** window to any value between 1000 and 3000 (2850 for 50hz models). A setting of 1000 is the least vigorous mixing and 3000 the most vigorous. Press the **START/STOP** button to begin the mixing action. The mixer will run until the set time has elapsed. To stop the cycle before the set time has elapsed, press the **START/STOP** button. To initiate a new timed cycle, press the **START/STOP** button again.

2.2 - "TIMED with SPEED ALARM"

Set the **TIME** and **SPEED** as described in 2.1. Press the **SPEED ALARM** button. The mixer will toggle between **SPEED ALARM OFF** and **ON** and the **SPEED** display will flash. Set the speed control by pressing the UP/DOWN arrow beneath the **SPEED** window. It is recommended to set the **SPEED ALARM** a minimum of 100 RPM slower than the mixer speed. Should the mixer not be able to maintain the speed within the ALARM SPEED setting, the mixer will beep.

2.3 - "CONTINUOUS"

The "Continuous" mode allows for unattended operation that will continue indefinitely until you stop the cycle by pressing the **START/STOP** button. To initiate the "Continuous" mode, press the DOWN arrow button beneath the **TIME** window until the illuminate display reaches "- -". Upon pressing the **START/STOP** button, the mixer will start and continue to run, unattended, until you once again press the **START/STOP** button. To restart the cycle, press the **START/STOP** button again.

2.4 - "CONTINUOUS with SPEED ALARM"

Set the "Continuous" **TIME** as described in 2.3. Set the **SPEED** as described in 2.1. Set the **SPEED ALARM** as described in 2.2.

2.5 - "SLEEP"

The mixer will automatically go into a "SLEEP" mode if the mixer is plugged in but not used for 10 minutes. The **TIME** and **SPEED** illuminated displays will go blank and the ON or OFF illuminated display will blink. Press the **START/STOP** button to "WAKE UP" the mixer.

Maintenance

The OHAUS DISRUPTOR GENIE mixer is built for long, trouble-free, dependable service. No lubrication or other technical user maintenance is required.

The unit should be given the care normally required for any electrical appliance. Avoid wetting or unnecessary exposure to fumes. Spills should be removed promptly. DO NOT use a cleaning agent or solvent on the front panel which is abrasive or harmful to plastics, nor one which is flammable. Always ensure the power is disconnected from the unit prior to any cleaning. If the unit ever requires service, contact your Ohaus representative.

Technical Support Information

For service assistance or technical support, please visit our website at www.ohaus.com to locate and contact the Ohaus office nearest you.

Serial Number: _____

Date of Purchase: _____

Supplier: _____

Technical Data

Equipment Rating

- Indoor use only
- Altitude: 0 to 2000 m
- Environmental: 0°C to 38°C (32°F to 100°F), max 95% relative humidity, non-condensing
- Overvoltage category (Installation category): II

Compliance

Compliance to the following standards is indicated by the corresponding mark on the product.

Mark	Standard
	This product complies with the applicable harmonized standards of EU Directives 2011/65/EU (RoHS), (EU)2015/863, 2014/30/EU (EMC) and 2014/35/EU (LVD). The EU Declaration of Conformity is available online at www.ohaus.com/ce .
	This product complies with the EU Directive 2012/19/EU (WEEE). Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. For disposal instructions in Europe, refer to www.ohaus.com/weee .
	EN 61326-1
	CAN/CSA-C22.2 No. 61010-1-12, CAN/CSA-C22.2 No. 61010-2-051 UL 61010-1, UL 61010-2-051

Global Notice

Warning: This is a Class A product. In a domestic environment this product may cause radio interference in which case the user may be required to take adequate measures.

ISED Canada Compliance Statement:

CAN ICES-003(A) / NMB-003(A) FCC Notice

FCC Compliance Statement:

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

ISO 9001 Registration

The management system governing the production of this product is ISO 9001 certified.

FCC Supplier Declaration of Conformity

Unintentional Radiator per 47CFR Part B

Trade Name: OHAUS CORPORATION

Model: DISRUPTOR GENIE e-D2..., Digital DISRUPTOR GENIE e-D2D...

Party issuing Supplier's Declaration of Conformity

Troemner LLC

201 Wolf Drive

Thorofare, NJ, 08086

USA

Phone: +1 856 686 1600

Responsible Party – U.S. Contact Information

Ohaus Corporation

8 Campus Drive, Ste. 105

Parsippany, NJ 07054

United States

Phone: +1 973 377 9000

Web: www.ohaus.com

LIMITED WARRANTY

OHAUS products are warranted against defects in materials and workmanship from the date of delivery through the duration of the warranty period. During the warranty period OHAUS will repair, or, at its option, replace any component(s) that proves to be defective at no charge, provided that the product is returned, freight prepaid, to OHAUS.

This warranty does not apply if the product has been damaged by accident or misuse, exposed to radioactive or corrosive materials, has foreign material penetrating to the inside of the product, or as a result of service or modification by other than OHAUS. In lieu of a properly returned warranty registration card, the warranty period shall begin on the date of shipment to the authorized dealer. No other express or implied warranty is given by OHAUS Corporation. OHAUS Corporation shall not be liable for any consequential damages.

As warranty legislation differs from state to state and country to country, please contact OHAUS or your local OHAUS dealer for further details.

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P/N: 31018507 A © 2026 Ohaus Corporation, all rights reserved