

# Instructions For Use

## Biomek 4000

RNAAdvance Blood Application

Total RNA Isolation from PAXgene  
Preserved Blood

B39494AA  
October 2013



Beckman Coulter, Inc.  
250 S. Kraemer Blvd.  
Brea, CA 92821 U.S.A.



**Biomek 4000 RNAdvance Blood Application**

**Instructions for Use**

PN B39494AA (October 2013)

© 2013 Beckman Coulter, Inc.

All rights reserved

The Beckman Coulter stylized logo, Biomek, SPRI are trademarks of Beckman Coulter, Inc. and are registered in the USPTO.

*All other trademarks, service marks, products, or services are trademarks or registered trademarks of their respective holders.*

Find us on the World Wide Web at:

[www.beckmancoulter.com](http://www.beckmancoulter.com)

Made in USA

# Revision History

*This document applies to the latest software listed and higher versions. When a subsequent software version changes the information in this document, a new issue will be released.*

**AA/Initial Issue, October 2013**

Biomek Software version 4.1 or higher

**Revision History**

# Safety Notice

Read all product manuals and consult with Beckman Coulter-trained personnel before attempting to operate instrument. Do not attempt to perform any procedure before carefully reading all instructions. Always follow product labeling and manufacturer's recommendations. If in doubt as to how to proceed in any situation, contact your Beckman Coulter Representative.

## Alerts for Warning, Caution, Important, and Note

---

### **WARNING**

WARNING indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

### **CAUTION**

CAUTION indicates a potentially hazardous situation, which, if not avoided, may result in minor or moderate injury. It may also be used to alert against unsafe practices.

**IMPORTANT** IMPORTANT is used for comments that add value to the step or procedure being performed. Following the advice in the Important adds benefit to the performance of a piece of equipment or to a process.

**NOTE** NOTE is used to call attention to notable information that should be followed during installation, use, or servicing of this equipment.

**Safety Notice**

Alerts for Warning, Caution, Important, and Note

# Contents

[Revision History](#), iii

[Safety Notice](#), v

[Introduction](#), xi

## **CHAPTER 1:** [Instrument Requirements](#), 1-1

[Instrument Requirements](#), 1-1

[Deck Configuration](#), 1-1

[System Configuration Settings](#), 1-3

[Customer Default Settings](#), 1-3

[To adjust user defaults and settings](#), 1-3

## **CHAPTER 2:** [Method Operation](#), 2-1

[Installation](#), 2-1

[Operating the RNAdvance Blood Method](#), 2-1

[Starting the RNAdvance Blood Method and](#)

[Establishing Default Protocol Settings](#), 2-1

[To open and start the method](#), 2-1

[Configuring the RNAdvance Blood Method](#), 2-5

[Quick Start Tab](#), 2-6

[To complete the Quick Start Tab Run Settings](#), 2-7

[Columns Tab](#), 2-9

[To complete the Columns Tab](#), 2-10

[Procedure Tab](#), 2-11

[To complete the Procedures Tab](#), 2-14

[Running the RNAdvance Blood Method](#), 2-15

[To run the RNAdvance Blood Method](#), 2-15

[Recovering from an Interrupted Method](#), 2-17

[Recovery Tab](#), 2-18

[To Recover the RNAdvance Blood Method](#), 2-19

[Abbreviations](#)

[Glossary](#)

[Index](#)

[Beckman Coulter, Inc.](#)

[Customer End User License Agreement](#)


# Illustrations

- |      |   |
|------|---|
| 1.1  | <a href="#">Recommended Deck Layout</a> , 1-1                       |
| 1.2  | <a href="#">Configure Method Outline</a> , 1-3                      |
| 1.3  | <a href="#">RNAdvance Blood Customer Default Values</a> , 1-4       |
| 2.1  | <a href="#">Method Outline</a> , 2-2                                |
| 2.2  | <a href="#">User Defaults and Settings, Top portion shown</a> , 2-3 |
| 2.3  | <a href="#">Default GUI Values, Bottom portion shown</a> , 2-4      |
| 2.4  | <a href="#">Quick Start Tab</a> , 2-6                               |
| 2.5  | <a href="#">Labware Selector</a> , 2-8                              |
| 2.6  | <a href="#">Columns Tab showing Sample Plate Selection</a> , 2-9    |
| 2.7  | <a href="#">Procedure Tab, Top portion shown</a> , 2-11             |
| 2.8  | <a href="#">Lock/Unlock Password</a> , 2-14                         |
| 2.9  | <a href="#">Instrument Setup</a> , 2-15                             |
| 2.10 | <a href="#">Editing or Expiring the Timer</a> , 2-17                |
| 2.11 | <a href="#">Recovery Tab</a> , 2-18                                 |

# Tables

- |     |  |
|-----|--|
| 1.1 | Requirements for RNAdvance Blood Method, 1-1             |
| 1.2 | RNAdvance Blood Customer Default Value Descriptions, 1-5 |
| 2.1 | Describing the RNAdvance Blood Primary Screen, 2-7       |
| 2.2 | Describing the Columns Screen, 2-10                      |
| 2.3 | Describing the Procedures Screen, 2-12                   |
| 2.4 | Procedure Default, Minimum and Maximum Values, 2-12      |
| 2.5 | Reagent and Deck Setup Guide, 2-16                       |
| 2.6 | Describing Recovery Tab, 2-19                            |

# Introduction



## Intended Use

---

The product is not intended or validated for use in the diagnosis of disease or other conditions.

## Overview

---

This introductory section contains the following topics:

- [How to Use This Manual](#)
- [About This Manual](#)

## How to Use This Manual

---

Use this manual to configure, test and operate the RNAdvance Blood method. It contains information on:

- Instrument Requirements
- System Configuration Settings
- Method Operation

## About This Manual

---

**NOTE** Screens in this manual may differ slightly from the screens in your RNAdvance Blood method.

The information in this *Instructions for Use* manual is organized as follows:

### **Chapter 1, Instrument Requirements**

Contains deck configuration and general settings for operating the RNAdvance Blood method.

### **Chapter 2, Method Operation**

Contains the settings and procedures for operating the RNAdvance Blood method.

The manual also includes a List of Abbreviations, a Glossary, and an Index.

## Conventions

- Selections that appear on a screen are in **bold face**.
- Information that is to be typed (entered) is in ***bold face and italics*** font.
- The software path to a specific function or screen appears with the greater than ( > ) symbol between succeeding screen options, like this: Select **Control Panel > Network and Sharing**.
- Links to information in another part of the document for additional information are in blue. To access the linked information, select the blue text.

# Instrument Requirements

## Instrument Requirements

Biomek 4000 with software version 4.1 or higher. Contact your Beckman Coulter Representative for the Method CD for your platform, call 1-800-369-0333, or email [reagentsupport@beckman.com](mailto:reagentsupport@beckman.com) if your sales contact information is unknown.

### Deck Configuration

See [Figure 1.1](#) for the recommended deck layout.

**Figure 1.1** Recommended Deck Layout



**Table 1.1** Requirements for RNAAdvance Blood Method

Type	Qty	Description	Part number
Tools	1	Gripper Tool	987371
	1	MP200 Eight-Tip Pipette Tool	986146
	1	MP1000 Eight-Tip Pipette Tool	A91112

**Table 1.1** Requirements for RNAdvance Blood Method

Type	Qty	Description	Part number
<b>ALPs</b>	1	Liquid Waste ALP	B21398
	1	Disposal ALP	609751
	5	Labware Holder	609120
	4	Tip Rack Holder	391910
	1	Static Peltier Device	A93938
	1	Orbital Shaker	379448
	1	Peltier Adaptor Plate (96 well round bottom)	A49568
	1	Off Deck Tool Rack Kit (left)	B21395
<b>Magnet Plate</b>	1	Super Magnet Plate	A32782
<b>Reservoirs</b>	2	Reservoir Frame	372795
	2	Half Reservoir	372786
	3	Quarter Reservoir	372790
	1	IMReservoir96 (Fisher Scientific)	50-995-860
	1	Quarter Reservoir (divided by length)	372788
<b>Consumables</b>	3	Biomek Span-8 P1000 Tips, Pre-Sterile with Barrier	B01124
	1	Biomek AP96 P250 Tips, Pre-Sterile with Barrier	717253
	1	96 Well PCR Plate <sup>a</sup>	AB-2800
	1	RK Riplate Deepwell Plate (Worldwide Medical)	99181000

a. ABGene 2800 or equivalent.

## System Configuration Settings

---

### Customer Default Settings


The system configuration settings section contains additional selections. Adjust the following fields as needed.

#### To adjust user defaults and settings

---


- 1 Open the Biomek Software.
  
- 2 Open the RNAdvance Blood method by clicking **File > Open**. The method opens and the method outline displays. See [Figure 1.2](#).

**Figure 1.2** Configure Method Outline



- 3 From the method outline click **User Defaults and Settings**. See [Figure 1.2](#). The screen in [Figure 1.3](#) displays.

Figure 1.3 RNAAdvance Blood Customer Default Values



**Table 1.2** RNAdvance Blood Customer Default Value Descriptions

Screen Element	Purpose	Do this
<b>A Default Protocol Settings</b>	These values can be adjusted for customer specified processes. These values will feed into the application interface as the default values.	Change values as needed based on your specific protocol.
<b>B Reset to Default Values</b>	Resets all values (except the Peltier Offset values), back to the factory default values.	Click the button to revert all values to the default values.
<b>C LIMS Folder</b>	File location to save Laboratory Information Management System (LIMS) data.	Enter the file path where you would like application data saved. For example: <b>c:\Program Files\LIMS</b>
<b>D Peltier Offsets</b>	During installation the Peltier Offset temperature values are determined based on the particular Peltier installed. Only change these values when recommended by Beckman Coulter.	Do not change.

**NOTE** It is recommended that you save the new version of the method under a different name to preserve the original. This allows you to go back to the original as a starting point if needed.

**IMPORTANT** Changing default protocol settings may cause invalid method results. Beckman Coulter's warranty applies only for unchanged Customer Default Protocol Settings.

To restore the Beckman Coulter Agencourt Default settings, click the **Reset to Default Values** Button.

**IMPORTANT** The Reset to Default Values button does not reset the Peltier temperature offsets. These offsets are determined during installation and are meant to ensure your samples reach the desired temperature during method execution. The actual Peltier temperature is equal to the sum of the set temperature, plus the offset temperature. The formula would be as follows:


$$\text{Sample Temp} + \text{Peltier Temp Offset} = \text{Actual Peltier Temp}$$

For example, to determine the actual temperature of the Peltier for Lysis:

If the Lysis Temperature is set to 55°C, and the Lysis Peltier Temperature Offset is set to 40°C, the actual temperature of the Peltier would be:

$$55^\circ\text{C} + 40^\circ\text{C} = 95^\circ\text{C}$$





## CHAPTER 2

# Method Operation

## Installation

---

Your Beckman Coulter Field Applications Scientist performs all Biomek and RNAdvance Blood Application installation and setup services. For service, see your Beckman Coulter Representative.

## Operating the RNAdvance Blood Method

---

**NOTE** For more information about Agencourt RNAdvance Blood and the Agencourt RNAdvance Blood protocol, see [www.beckmancoulter.com](http://www.beckmancoulter.com). The automated method has been optimized for automation and may differ from the manual protocol.

### Starting the RNAdvance Blood Method and Establishing Default Protocol Settings

The RNAdvance Blood Method is shipped with a set of standard protocol values. But these values may need to be modified for your laboratory. These values populate the RNAdvance Blood method application.

#### To open and start the method


---

- 1** Open the Biomek software.
  
- 2** Open the RNAdvance Blood Project by clicking **Project > Open Project**.

- 
- 3** Open the RNAdvance Blood method by clicking **File > Open**. Select **RNAdvanceMethod**.


**NOTE** Your method structure may look larger and may not have + and – icons to expand and collapse the structure. To change the appearance of the method structure, click **Options > Preferences** and then click **View** in the menu bar on the left. Use this option to change the appearance of the method structure.

**Figure 2.1** Method Outline



- 
- 4** From the method outline, click **User Defaults and Settings**. See [Figure 2.1](#). The User Defaults and Settings screen displays.


**NOTE** You can skip steps 4 and 5 if already configured previously.

**Figure 2.2** User Defaults and Settings, Top portion shown

## Method Operation

### Operating the RNAdvance Blood Method

**Figure 2.3** Default GUI Values, Bottom portion shown



- 5** Change or ensure that each variable is the correct default value. These values can be changed at run time. For details on this screen see [Table 1.2](#).

**IMPORTANT** Changing default protocol settings may cause invalid method results. Beckman Coulter's warranty applies only for unchanged Customer Default Protocol Settings.

To restore default settings, click the **Reset to Default Values** Button.

**IMPORTANT** The Reset to Default Values button does not reset the Peltier temperature offsets. These offsets are determined during installation and are meant to ensure your samples reach the desired temperature during method execution. The actual Peltier temperature is equal to the sum of the set temperature, plus the offset temperature. The formula would be as follows:

$$\text{Sample Temp} + \text{Peltier Temp Offset} = \text{Actual Peltier Temp}$$

For example, to determine the actual temperature of the Peltier for Lysis:

If the Lysis Temperature is set to 55°C, and the Lysis Peltier Temperature Offset is set to 40°C, the actual temperature of the Peltier would be:

$$55^{\circ}\text{C} + 40^{\circ}\text{C} = 95^{\circ}\text{C}$$

- 
- 6** From the Biomek software, click **Execution > Run** .
- 

## Configuring the RNAdvance Blood Method

The graphical user interface appears and allows you to make some selections. To achieve optimal performance for different sample types, you may need to adjust some selections.

The RNAdvance Blood Method contains a graphical user interface to ensure proper processing of samples. This interface contains four tabs, three of which step you through the setup process until the method is ready to run.

The four tabs include:

**Quick Start** — Use this Tab to configure the sample and finals plates, and to track reagent lot numbers. See the [Quick Start Tab](#) on page 2-6.

**Columns** — Use this Tab to define the number and location of active columns on the sample plates. See the [Columns Tab](#) on page 2-9.

**Procedure** — Use this Tab to configure the protocols. See the [Procedure Tab](#) on page 2-11.


**Recovery** — Use this Tab only if you need to restart the method and recover it at the point it was interrupted. See the [Recovery Tab](#) on page 2-18.

## Method Operation

### Operating the RNAdvance Blood Method

#### Quick Start Tab

Figure 2.4 Quick Start Tab



- A. Quick Start Tab
- B. Input Values
- C. Forward Button
- D. Labware Selector
- E. Output Values
- F. Tracking
- G. Deck Display
- H. Run Button
- I. Abort Button

**Table 2.1** Describing the RNAdvance Blood Primary Screen

Screen Element	Purpose	Do this	Notes
<b>A Quick Start Tab</b>	Allows you to define the sample inputs, outputs and whether to track the method.	Complete all fields.	N/A
<b>B Input Values</b>	Allows you to define the sample volume of the run.	Complete all fields.	N/A
<b>C Forward Button</b>	Allows you to advance to the next screen.	Click the button to advance.	Use this button or just click the next tab.
<b>D Labware Selector</b>	Allows you to determine the define the elution plate.	Select the desired elution plate type.	
<b>E Output Values</b>	Allows you to define the elution and transfer volumes.	Select whether to transfer to a final plate and then select the correct elution and (if transferring) transfer volumes.	Most laboratories choose to transfer the eluate to a new plate.
<b>F Tracking</b>	Allows you to track RNAdvance Blood reagents automatically by Lot numbers.	To track RNAdvance Blood reagents automatically, click Yes and record the User ID and the applicable lot numbers.	N/A
<b>G Deck Display</b>	Displays the correct deck setup based on your inputs.	Check the Deck Display against the Biomek 4000 Deck.	N/A
<b>H Run Button</b>	Allows you to start running the method.	Click the Run button.	This Button appears on all four tabs.
<b>I Abort Button</b>	Allows you to stop the method.	While the Application GUI is displayed, click this button to abort.	This Button appears on all four tabs.

**To complete the Quick Start Tab Run Settings**


**1** On the **Quick Start Tab** select the sample volume in  $\mu\text{L}$ .

**2** Select whether to transfer the eluate to the final plate.

**NOTE** Most laboratories choose to transfer the eluate to a new plate.

- 
- 3 Click the **Labware Selector**. The Labware Selector screen displays. See [Figure 2.5](#).

**Figure 2.5** Labware Selector




- 
- 4 Select the Elution Labware Plate Type used in your RNAdvance Method. Click **OK**.
- 
- 5 If transferring, select the Elution and Transfer Volume in  $\mu\text{L}$ . The default volume that Beckman Coulter recommends is  $40 \mu\text{L}$ .
- 
- 6 Select the Transfer Volume in  $\mu\text{L}$ . The default volume that Beckman Coulter recommends is  $30 \mu\text{L}$ .
- NOTE** Beckman Coulter recommends a transfer volume of at least  $10 \mu\text{L}$  less than the elution volume to avoid any bead carryover.
- 
- 7 Select whether the application should write a log file. If **Yes**, follow the steps below. If **No**, skip to step 8.
- Enter the **User ID**.
  - Enter the **RNAdvance Kit Lot #**.
  - Enter the **Elution Buffer Lot #**.
  - Enter the **Ethanol Lot #**.
  - Enter the **Isopropanol Lot #**.
  - Enter the **DNase Lot #**.
  - Enter the **DNase Buffer Lot #**.
  - Enter **Notes** if needed.
- NOTE** If this feature is used, the RNAdvance Blood application will record a log file to the computer folder defined on the Configure Method screen. See [Customer Default Settings](#) on page 1-3 for information on how to access this screen.

- 8 Click the **Columns Tab** or the button. The Columns Tab Setup Screen is displayed.

## Columns Tab

**Figure 2.6** Columns Tab showing Sample Plate Selection



- A. Columns Tab  
B. Back Button

- C. 96-Well Plate Column

**WARNING**

**Orbital shaker hazard.** Distribute samples evenly across the plate to prevent the plate from becoming unbalanced. An unbalanced plate may propel hazardous substances from the orbital shaker. Refer to the *Biomek 4000 ALPs and Accessories Manual* for complete safety information.

**Table 2.2** Describing the Columns Screen


Screen Element	Purpose	Do this	Notes
<b>A Columns Tab</b>	Allows you to define active columns.	Select columns in any order. Rows cannot be selected.	Any combination of individual columns can be made active, but cells cannot be made active as rows or as individual wells.
<b>B Back Button</b>	Allows you to return to the previous tab.	Select this button to return to the Quick Start Tab.	N/A
<b>C 96-Well Plate</b>	Allows you to visualize and select columns more easily.	Click on the heading buttons above the column to select and deselect columns to be run on the Biomek 4000.	<a href="#">Figure 2.6</a> shows all columns selected. Your application may be different.

**To complete the Columns Tab**

- 1 On the **Columns Tab**, select the active columns by clicking the appropriate 96-Well Plate Column Heading. To highlight all or none of the cells, click the top-left Well Plate Column Heading. This acts as a toggle.
  
- 2 Click the **Procedure Tab** or the ➤ button. The Procedure Tab Setup Screen displays.

## Procedure Tab

Figure 2.7 Procedure Tab, Top portion shown



- A. Procedure Tab
- B. Lock/Unlock Button
- C. Reset Default Values
- D. Procedure Values

**Method Operation**

## Operating the RNAdvance Blood Method

**Table 2.3** Describing the Procedures Screen

Screen Element	Purpose	Do this	Notes
<b>A Procedure Tab</b>	Allows you to define and confirm all values in the RNAdvance Blood method.	Select all correct values. Changed values are highlighted in yellow.	The starting values are determined by the User Default and Settings step. See <a href="#">Customer Default Settings</a> on page 1-3.
<b>B Lock/Unlock Indicator</b>	Indicates whether the method protocol inputs are locked with a password.	To lock the protocol values, click the lock tab and enter a password.	Only the values on the Procedures Tab are locked.
<b>C Reset Default Values</b>	Returns all values to the defaults established in the User Defaults and Settings step. See <a href="#">Customer Default Settings</a> on page 1-3.	Click this button to return all values to the default values.	All manual entries will be lost without recourse.
<b>D Procedure Values</b>	Allows you to change values in the RNAdvance Blood method.	Select the correct procedure values. Changed values are highlighted in yellow.	For a list of all procedure value defaults, minimums and maximums, see <a href="#">Table 2.4</a> .

The values for these steps originate from the Default Protocol Settings. You can access these settings from the User Defaults and Settings step. See [Customer Default Settings](#) on page 1-3.

Use this Tab to customize the automated protocol at runtime. Most RNAdvance Blood variables can be adjusted. All protocol settings are written to memory when the run starts and recalled for subsequent runs.

[Table 2.4](#) lists all automated method steps and the default, minimum and maximum values.

**Table 2.4** Procedure Default, Minimum and Maximum Values

Step	Step Description	Units	Default	Minimum	Maximum
<b>Lysis</b>					
1	Lysis Volume	µL	300	100	500
2	Proteinase K Solution Volume	µL	40	20	100
3	Lysis Shake speed on the Orbital Shaker	RPM	1200	1000	1200
	Lysis Shake time	minutes	3	1	60
4	Lysis Sample temperature <sup>a</sup>	°C	55	20	100
	Lysis Incubation time	minutes	25	1	60
<b>Bind 1</b>					
5	Bind 1 Volume	µL	410	200	500
6	Bind 1 Tippmixing	minutes	15	1	20
7	Bind 1 Incubation time	minutes	5	1	15
8	Bind 1 Settle time	minutes	10	5	20
<b>Wash</b>					

**Table 2.4** Procedure Default, Minimum and Maximum Values (Continued)

<b>Step</b>	<b>Step Description</b>	<b>Units</b>	<b>Default</b>	<b>Minimum</b>	<b>Maximum</b>
9	Wash Volume	µL	800	500	1000
10	Wash Shake 1 speed on the Orbital Shaker	RPM	1300	800	1300
11	Wash Shake 1 time	minutes	4	1	10
12	Wash Shake 2 speed on the Orbital Shaker	RPM	1000	800	1100
13	Wash Shake 2 time	minutes	2	1	10
14	Wash Settle time	minutes	8	5	15
<b>Ethanol 1</b>					
15	Ethanol 1 Volume	µL	750	500	1000
16	Ethanol 1 Shake speed	RPM	1100	800	1200
17	Ethanol 1 Shake time	minutes	2	1	5
18	Ethanol 1 Settle time	minutes	3	1	10
19	Ethanol 1 Dry time	minutes	5	1	15
<b>DNase</b>					
17	DNase Solution Volume	µL	100	20	150
19	DNase Shake speed on the Orbital Shaker	RPM	1000	800	1200
	DNase Shake time	minutes	2	1	5
20	DNase Sample Temperature <sup>b</sup>	°C	37	20	100
	DNase Incubation time	minutes	15	1	20
<b>Rebind with Bind 2</b>					
21	Bind 2 Volume	µL	200	100	500
22	Bind 2 Shake Speed on the Orbital Shaker	RPM	1000	800	1200
	Bind 2 Shake time	minutes	8	1	10
23	Bind 2 Incubate time	minutes	5	1	10
24	Bind 2 Settle time	minutes	5	1	10
<b>Ethanol 2</b>					
25	Ethanol 2 Cycles	time(s)	2	1	3
26	Ethanol 2 Volume	µL	750	500	1000
27	Ethanol 2 Settle time	minutes	4	1	5
28	Ethanol 2 Shake Speed on the Orbital Shaker	RPM	800	800	1200
29	Ethanol 2 Dry time	minutes	5	0	10
<b>Elution</b>					
30	Elution volume	µL	40	20	100

**Table 2.4** Procedure Default, Minimum and Maximum Values (*Continued*)

Step	Step Description	Units	Default	Minimum	Maximum
31	Elution Shake speed on the Orbital Shaker	RPM	1000	800	1200
32	Elution Shake time	minutes	2	1	15
33	Elution Incubation time	minutes	1	0	5
34	Elution Settle time	minutes	2	1	5
<b>Transfer</b>					
35	Eluent Transfer volume	µL	35	20	100

a. Set Lysis Sample Temperature to the desired Temperature in the well. See [D Peltier Offsets](#) in [RNAdvance Blood Customer Default Values](#).

b. Set DNase Sample Temperature to the desired Sample Temperature in the well. See [D Peltier Offsets](#) in [RNAdvance Blood Customer Default Values](#).

**IMPORTANT** Changing default protocol settings may cause undesirable results. Beckman Coulter's warranty applies only for unchanged Default Protocol Settings.

### To complete the Procedures Tab

- 1 On the **Procedures Tab** change or ensure that each variable is correct. Use the up and down arrow keys to select values. A value with a yellow background indicates that the value has changed from the default values.

**NOTE** To change the default values so that the variables here populate according to custom values, see [Customer Default Settings](#) on page 1-3.

- 2 **Optional:** To lock the values of the procedure, which prevents inadvertent changes to variables, click the **Lock/Unlock Button** . The Lock/Unlock screen displays.

**Figure 2.8** Lock/Unlock Password

**NOTE** For first-time use, leave the first field, **Enter Current Password**, blank.

- a. Enter a new password.
- b. Re-enter the password.

- c. Click **OK**. All variables become non-editable. To change locked values, click the Lock/Unlock Button and enter the password.


**NOTE** If a password is forgotten, contact Reagent Support at [reagentsupport@beckman.com](mailto:reagentsupport@beckman.com).

- 3 **Optional:** To reset all variables to the default values, click the

Reset Default Values Button  . All changed variables revert to default values.

Once the method begins running the software prompts you to set up reagents, tools and labware. See [Table 2.5](#) for the sequence of this setup. Follow the prompts on the screen for actual values.

**Figure 2.9** Instrument Setup



### Running the RNAAdvance Blood Method

**IMPORTANT** To ensure best possible outcomes, the RNAAdvance Blood bottle containing Bind 1 must be shaken and mixed so that the beads are homogenous in solution. Run the procedure immediately following the placement and filling of RNAAdvance Blood reagents. Failure to run the method immediately following the filling of reagents may cause the method to produce poor results.

After entering all protocol variables, positioning labware, and filling all reagent reservoirs, the RNAAdvance Blood Method is ready to run.

#### To run the RNAAdvance Blood Method

- 1 Confirm that all values in the user interface are correct for all RNAAdvance Blood application variables. See [Configuring the RNAAdvance Blood Method](#).

- 2 From any of the user interface tabs, click the green run arrow .


**Table 2.5** Reagent and Deck Setup Guide

Step	Reagent/Tool/ Labware type	Labware Name	Notes
Mix Reagents	DNase	N/A	Nuclease Free Water DNase I buffer DNase I Use freshly prepared DNase I solution
	Ethanol	N/A	Use freshly prepared 70% ethanol
Put Tools on the Deck	Gripper	N/A	N/A
	MP1000	N/A	N/A
Put Labware on the Deck	p1000 Barrier Tips	ReagentTips SampleTips_1 BloodTips	N/A
	p200 Barrier Tips	ElutionTips_1	N/A
	96R Super Magnet Plate	Mag_1	N/A
	p1000 Barrier Tips	SampleTips_1	N/A
	96-well plate	ElutionPlate_1	N/A
	Half Modular Reservoir	Res-Bind-WashBuffer	Bind 1 in Section 1 Washbuffer in Section 2
	AGCT_96RitterDeep Square_GSL	PrepPlate1_1	Sample in selected columns 1-8
	RNABlood_Mod_Res	Res-LY-PK-DN-RB-EL	Lysis in section 1 PK in section 2 DNase in section 3 Bind2 in section 4 Elute in section 5 Sections 6 and 7 are empty
	IMReservoir96	Res-Ethanol	Ethanol in section 1

The Biomek instrument may appear to be idle at times, during which a timer appears. When this occurs an incubation time is being observed before the next step can be performed. See [Figure 2.10](#).

#### Options for the timer include

- Editing the timer duration to the precise number of hours, minutes and seconds.
- Adding one minute to the elapsed time.
- Expiring the timer immediately.

**Figure 2.10** Editing or Expiring the Timer

When the method ends, the tips are unloaded and the Biomek Start arrow returns to green.

## Recovering from an Interrupted Method

There are several reasons a method may be interrupted:

- Power failure
- Tips did not load properly
- Insufficient reagent
- Instrument crash

The RNAAdvance Blood Method provides ways to recover or restart the method efficiently:


- Restart from the point where the method stopped.
- Reset all steps for the plate.
- Reset the deck: return all moveable labware to their starting positions.

## Method Operation

### Operating the RNAdvance Blood Method

## Recovery Tab

Figure 2.11 Recovery Tab



- A. Recovery Tab
- B. Plate Column
- C. Reset Current Selection
- D. Reset Current Plate
- E. Set All Sections of Current Plate to Match Current Selection
- F. Finished Tasks
- G. Unfinished Tasks
- H. Reset the Deck

**Table 2.6** Describing Recovery Tab

Screen Element	Purpose	Do this
<b>A Recovery Tab</b>	Allows you to see where the RNAdvance Blood Method was interrupted to more efficiently recover it.	Use this tab to reset the specific columns to be run again.
<b>B Plate Column</b>	Allows you to visualize the individual columns of the Plate.	Click the column tab to check the status of each column.
<b>C Reset Current Selection</b>	Allows you to restart the current column to the beginning of the method.	Click this button to restart the selected column to the beginning of the method.
<b>D Reset Current Plate</b>	.Allows you to restart the method from the beginning.	Click this button to restart the RNAdvance Blood method from the beginning.
<b>E Set All Sections of Current Plate to Match Current Selection</b>	Allows you to set all active columns to the state of the currently selected column.	Click this button to set all active columns to start from the state of the currently selected column.
<b>F Finished Tasks</b>	Cleared checkboxes show completed tasks by column.	To rerun a step for a column, check the appropriate checkbox.
<b>G Unfinished Tasks</b>	Checked checkboxes show unfinished tasks. The green bar indicates the task next in line.	To skip a step for a column, clear (uncheck) the appropriate checkbox.
<b>H Reset the Deck</b>	Moves all labware to the starting position.	Click this button to show how the deck should be reset.

**To Recover the RNAdvance Blood Method**

- 1 Address the source of the interruption.

---

- 2 Restart the RNAdvance Blood Method. The **Recovery** tab will be displayed.

---

- 3 Confirm that the state of the deck on the Biomek 4000 matches the Deck Display. If the deck on the Biomek 4000 does NOT match the Deck Display, correct the deck on the Biomek 4000 to match it by dragging labware in the Deck Display.

---

- 4 Click any Tasks that show finished, but need to be recovered.  
For example, if you can see that column eight does not have Bind 1 reagent, but the Recovery Tab shows that it was added, click that task in the Column 8 tab to make it an unfinished task.

---

- 5 Click the **Run** button .

## **Method Operation**

Operating the RNAdvance Blood Method

# Abbreviations

**µL** — microliter

**ALP** — Automated Labware Positioner

**DNA** — Deoxyribonucleic Acid

**GUI** — Graphical User Interface

**mL** — milliliter

**LIMS** — Laboratory Information Management System

**m** — meter

**RNA** — Ribonucleic Acid

**SPRI** — Solid Phase Reverse Immobilization

**Vol** — Volume

## Abbreviations

# Glossary

**assay** — procedure of repeat testing to determine the assigned value for a given lot and level of control.

**beads** — in SPRI technology, magnetic, uniform microparticles.

**dead volume** — in an automated system, the amount or volume of a sample or reagent that cannot be picked up by the pipette tip.

**elution buffer** — buffer which elutes DNA from magnetic particles.

**ethanol wash** — washes the magnetic beads with 85% ethanol to remove contaminants.

**gripper** — a tool that includes mechanical fingers to grip labware. Gripper fingers grasp labware along the long side and move the labware from one location on the Biomek deck to another. The gripper contains two fingers: a double gripper located to the front and a single gripper located to the back of the tool.

**PAXgene** — Preservative and stabilizer for nucleic acids in whole blood.

**reservoir** — one-well labware receptacle holding liquid to be used in a method.

**samples** — PAXgene preserved human blood.

**supernatant** — a liquid lying above a solid residue after crystallization, precipitation, centrifugation, or other process.

## Glossary

# Index

## A

- ABGene 2800 plate, [1-2](#)
- abort button
  - description of screen element, [2-7](#)
- Agencourt
  - configuration of, [2-5](#)
  - RNAadvancce protocol, [2-1](#)
  - RNAAdvance primary screen, [2-7, 2-12](#)
  - SuperMagnet Plate, [1-2](#)
- ALPs
  - defined, [Abbreviations-1](#)
- arrow
  - green, [2-15](#)
- assay
  - defined, [Glossary-1](#)

## B

- beads
  - defined, [Glossary-1](#)
- Biomek
  - as instrument requirement, [1-1](#)
  - occasionally idle, [2-16](#)
  - timer, [2-16](#)
- boxes
  - of tips as instrument requirement, [1-2](#)

## C

- columns tab
  - description of, [2-5](#)
  - screen elements, [2-9](#)
- configuration
  - RNAAdvance method, [2-5](#)
- configuring
  - system settings, [1-3](#)
- consumables, [1-2](#)
- customer default protocol settings, [2-2](#)

## D

- dead volume
  - defined, [Glossary-1](#)
- default protocol settings, [2-2](#)
- default values
  - description of screen element, [2-12](#)
- description
  - lock indicator, [2-7](#)

## E

- elution buffer
  - defined, [Glossary-1](#)
- ethanol wash
  - defined, [Glossary-1](#)

## F

- finished tasks
  - description of screen element, [2-19](#)
- forward button
  - description of screen element, [2-7](#)

## G

- general settings, [1-3](#)
- green arrow, [2-5, 2-15](#)
- gripper
  - defined, [Glossary-1](#)
- gripper tool, [1-1](#)
  - as instrument requirement, [1-1](#)

## I

- incubation
  - timing of, [2-16](#)
- input values
  - description of screen element, [2-7](#)
- instrument requirements, [1-1](#)

**L**

- lock indicator
  - description of, [2-7](#)
- lock/unlock indicator
  - description of screen element, [2-12](#)

**M**

- magnet plates
  - as instrument requirement, [1-2](#)
- manual
  - conventions, [xii](#)
  - navigation tools and tips, [xii](#)
- manual protocol, [2-1](#)
- method
  - changing appearance of structure, [2-2](#)
  - configuration of, [2-5](#)
  - running, [2-1](#)
- methods
  - causes of interruptions, [2-17](#)
  - recovering, [2-17](#)
  - ways to recover, [2-17](#)
- MP 200 tool, [1-1](#)
  - as instrument requirement, [1-1](#)

**O**

- opening
  - project, [2-1](#)
- output
  - description of screen element, [2-7](#)

**P**

- plates
  - thermocycling, [1-2](#)
- procedure tab
  - description of, [2-5](#)
  - screen elements, [2-11](#)
- procedure values
  - description of screen element, [2-12](#)
- project
  - opening, [2-1](#)
- protocol
  - default settings, [2-2](#)
  - manual, [2-1](#)

**Q**

- quarter reservoir

- as instrument requirement, [1-2](#)
- quick start tab
  - description of, [2-5](#)
  - screen elements, [2-6](#)

**R**

- recovering method, [2-17](#)
- recovery tab
  - description of, [2-5](#)
  - screen elements, [2-18](#)
- requirements
  - instrument, [1-1](#)
- reservoir
  - as instrument requirement, [1-2](#)
  - defined, [Glossary-1](#)
  - quarter, [1-2](#)
- reservoir frame, [1-2](#)
- reset all plates
  - description of screen element, [2-19](#)
- reset deck
  - description of screen element, [2-19](#)
- reset this plate
  - description of screen element, [2-19](#)
- RNAdvance
  - primary screen definitions, [2-7, 2-12](#)
- run button
  - description of screen element, [2-7](#)
- running method, [2-1](#)

**S**

- screen elements
  - columns tab, [2-9](#)
  - procedure tab, [2-11](#)
  - quick start tab, [2-6](#)
  - recovery tab, [2-18](#)
- settings
  - configuration, [1-3](#)
  - default protocol, [2-2](#)
  - general, [1-3](#)
- supernatant
  - defined, [Glossary-1](#)
- system configuration, [1-3](#)

**T**

- thermocycling plate
  - as instrument requirement, [1-2](#)
- timer

of Biomek, [2-16](#)  
tips  
    boxes of as instrument requirement, [1-2](#)  
    unloading, [2-17](#)  
tool  
    gripper, [1-1](#)  
    MP 200, [1-1](#)  
    wash 8, [1-1](#)  
tracking  
    description of screen element, [2-7](#)  
trademarks, [ii](#)

## **U**

unfinished tasks  
    description of screen element, [2-19](#)

## **W**

wash 8 tool, [1-1](#)  
    as instrument requirement, [1-1](#)



# Beckman Coulter, Inc.

## Customer End User License Agreement

This Product contains software that is owned by Beckman Coulter, Inc. or its suppliers and is protected by United States and international copyright laws and international trade provisions. You must treat the software contained in this Product like any other copyrighted material. This license and your right to use the Product terminate automatically if you violate any part of this agreement.

This is a license agreement and not an agreement for sale. Beckman Coulter hereby licenses this Software to you under the following terms and conditions:

**You May:**

1. Use this software in the computer supplied to you by Beckman Coulter;
2. Maintain one copy of this software for backup purposes (the backup copy shall be supplied by Beckman Coulter);
3. After written notification to Beckman Coulter, transfer the entire Product to another person or entity, provided you retain no copies of the Product software and the transferee agrees to the terms of this license agreement.

**You May Not:**

1. Use, copy or transfer copies of this Software except as provided in this license agreement;
2. Alter, merge, modify or adapt this Software in any way including disassembling or decompiling;
3. Loan, rent, lease, or sublicense this Software or any copy.

**Limited Warranty**

Beckman Coulter warrants that the software will substantially conform to the published specifications for the Product in which it is contained, provided that it is used on the computer hardware and in the operating system environment for which it was designed. Should the media on which your software arrives prove defective, Beckman Coulter will replace said media free of charge within 90 days of delivery of the Product. This is your sole remedy for any breech of warranty for this software.

Except as specifically noted above, Beckman Coulter makes no warranty or representation, either expressed or implied, with respect to this software or its documentation including quality, performance, merchantability, or fitness for a particular purpose.

**No Liability for Consequential Damages**

In no event shall Beckman Coulter or its suppliers be liable for any damages whatsoever (including, without limitation, damages for loss of profits, business interruption, loss of information, or other pecuniary loss) arising out of the use of or inability to use the Beckman Coulter Product software. Because some states do not allow the exclusion or limitation of liability for consequential damages, the above limitation might not apply to you.

**General**

This agreement constitutes the entire agreement between you and Beckman Coulter and supersedes any prior agreement concerning this Product software. It shall not be modified except by written agreement dated subsequent to the date of this agreement signed by an authorized Beckman Coulter representative. Beckman Coulter is not bound by any provision of any purchase order, receipt, acceptance, confirmation, correspondence, or otherwise, unless Beckman Coulter specifically agrees to the provision in writing. This agreement is governed by the laws of the State of Florida.

[www.beckmancoulter.com](http://www.beckmancoulter.com)



© 2013 Beckman Coulter, Inc.  
All Rights Reserved