Vi-CELL MetaFLEX

Safety Notices and Basic Troubleshooting

Read all product manuals and consult with Beckman Coulter-trained personnel before unpacking, setting up or operating this equipment. Pay attention to all warning and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

Do not attempt to perform any procedure before carefully reading all instructions. Always follow product labeling and manufacturer's recommendations. If in doubt as to how to proceed in any situation, contact your Beckman Coulter Representative.

To make sure that the protection provided by this equipment is not impaired, do not use or install this equipment in any manner other than that specified in this manual.

Beckman Coulter, Inc. urges its customers and employees to comply with all national health and safety standards such as the use of barrier protection. This may include, but is not limited to, protective eyewear, gloves, and suitable laboratory attire when operating or maintaining this or any other automated laboratory analyzer.

For labeling updates, go to www.beckmancoulter.com and download the most recent manual or system help for your instrument.

Safety Notices

Alerts for Warning and Caution

WARNING indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

In this document the signal word WARNING is only used to indicate the possibility of personal injury. It is not used to indicate the possibility of erroneous data.

 CAUTION

CAUTION indicates a potentially hazardous situation, which, if not avoided, may result in minor or moderate injury. It may also be used to alert against unsafe practices.

In this document the signal word **CAUTION** is only used to indicate the possibility of damage to the instrument or of erroneous data.

 WARNING

Risk of operator injury if:

- All doors, covers, and panels are not closed and secured in place prior to and during instrument operation.
- Instrument alarms and error messages are not acknowledged and acted upon.
- You mishandle broken parts.
- Doors, covers, and panels are not opened, closed, removed, and/or replaced with care.
- Improper tools are used for troubleshooting.

To avoid injury:

- Keep doors, covers, and panels closed and secured in place while the instrument is in use.
- Take full advantage of the safety features of the instrument.
- Acknowledge and act upon instrument alarms and error messages.
- Report any broken parts to your Beckman Coulter Representative.
- Open/remove and close/replace doors, covers, and panels with care.
- Use the proper tools when troubleshooting.

 CAUTION

System integrity could be compromised and operational failures could occur if:

- This equipment is used in a manner other than specified. Operate the instrument as instructed in the Product Manuals.
- You introduce software that is not authorized by Beckman Coulter into your computer. Only operate your system's computer with software authorized by Beckman Coulter.
- You install software that is not an original copyrighted version. Only use software that is an original copyrighted version to prevent virus contamination.
- You connect external devices such as thumb drives, external hard drives, and CD/DVDs. Ensure that all external devices are free from viruses before connecting.

 CAUTION

If you purchased this product from anyone other than Beckman Coulter or an authorized Beckman Coulter distributor, and, it is not presently under a Beckman Coulter service maintenance agreement, Beckman Coulter cannot guarantee that the product is fitted with the most current mandatory engineering revisions or that you will receive the most current information bulletins concerning the product. If you purchased this product from a third party and would like further information concerning this topic, contact your Beckman Coulter Representative.

 WARNING

California Proposition 65:

This product may contain chemicals known to the State of California to cause cancer and birth defects or other reproductive harm.

General Warning and Cautions

 WARNING

Risk of infection. Dispose and handle all used sampling devices, quality control (QC) ampoules, Solution Packs, Sensor Cassettes, Inlet Probes, Inlet Gasket Holders, Inlet Connector Gaskets and Inlet Modules as biohazardous waste. Follow your local regulations.

Use universal precautions when working with pathogenic materials. Means must be available to decontaminate the instrument and to dispose of biohazardous waste.

 WARNING

Risk of infection. Make sure you do not prick or scratch yourself on the Inlet Probe.

Use universal precautions when working with pathogenic materials. Means must be available to decontaminate the instrument and to dispose of biohazardous waste.

 WARNING

Risk of infection. Only let authorized personnel collect and work with biologic samples. Make sure to wear gloves.

Use universal precautions when working with pathogenic materials. Means must be available to decontaminate the instrument and to dispose of biohazardous waste.

 WARNING

Risk of infection. Make sure that you wear gloves during replacement and maintenance procedures.

Use universal precautions when working with pathogenic materials. Means must be available to decontaminate the instrument and to dispose of biohazardous waste.

 CAUTION

Risk of operator injury. The sensor cassette door and sample inlet probe are automated moving parts. Use caution around these areas to avoid injury.

Precautionary labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed. A symbol, if noted on the instrument, will be included with a warning or caution statement in the manual.

Symbol/icon	Explanation
	Caution / Warning
	Biohazard
	Possible laser radiation inside.
	Sample mixer
	Keep dry
	Do not use if package is damaged

Symbol/icon	Explanation
	Do not re-use
	Use by
	Contains sufficient for <n> tests
	Temperature limit
	Lot no.
	Catalog no. (product code)
	Consult instructions for use
	Manufacturer
	Keyboard
	CE = Conformité Européenne A “CE” mark indicates that a product has been assessed before being placed on the market, and has been found to meet European Union safety, health, and/or environmental protection requirements. It is a mandatory marking for various product categories that are placed into service within the European Economic area.
	COM gate (scanner/barcode reader)
	VGA (monitor)
	Mouse
	Network
	Off
	On
	cNRTLus certification
	The RCM mark is for products that comply with Australian Communications Media Authority (ACMA) EMC requirements.

Symbol/icon	Explanation
	Fuse symbol
	USB
	This symbol indicates that Beckman Coulter Medical ApS and its distributors within the European Union (EU) and associated states have taken the necessary steps to comply with the "DIRECTIVE 2012/19/EU OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 4 July 2012 on waste of electrical and electronic equipment (WEEE)". Equipment marked with this symbol must not be disposed of as household waste but as electronic waste in accordance with local legislation. Please note that equipment contaminated with potentially infectious substances, such as body fluids, must be decontaminated before recycling. If this is not possible, the equipment must be disposed of as biohazardous material. Contact your local Beckman Coulter representative for instructions.
	Marks compliance with SJ/T 11363-2006 (China RoHS). The number in the symbol shows the environmentally friendly use period in years.
	Marks compliance with SJ/T 11363-2006 (China RoHS). The product contains no restricted substances above the prescribed thresholds.
	UKCA=UK Conformity Assessed A "UKCA" mark indicates that a product has been assessed before being placed in the UK market, and has been found to meet UK safety, health, and/or environmental protection requirements.

Product Description

Risk of erroneous results. Low pO_2 levels can have an effect on the linearity of glucose measurements. This can lead to incorrect low glucose results. Note that cGlu linearity is not specified when the pO_2 level is less than 10 mmHg (1.3 kPa).

Power-Supply Cords

The power-supply cord and plug of the analyzer must comply with national regulations. If the regulations are not complied with, the equipment may be damaged.

External devices connected to the analyzer must be in compliance with the standard UL 60950 for US and IEC 60950 for Europe. If you do not do this, the equipment may be damaged.

Turning On the Instrument

CAUTION

1. If the system fails to start properly, check first to see whether the power cable and connection cables are properly connected.
2. Never shut off the power or disconnect a data cable while the system is performing a task. Doing so can result in data loss or damage to the system.

Ampoule-Based QC Measurements

WARNING

Risk of biohazardous contamination if you have skin contact with the ampoule and its contents. The ampoule could contain residual biological material and must be handled with care. Clean up spills immediately. Dispose of the ampoule in accordance with your local regulations and acceptable laboratory procedures.

Anticoagulants

CAUTION

Risk of erroneous results. Do not use EDTA as an anticoagulant when working with blood samples, as it will cause incorrect pH, $p\text{CO}_2$, cNa^+ , cK^+ and cCa^{2+} results and have an effect on subsequent cCa^{2+} measurements.

CAUTION

Risk of instrument damage. Do not use EDTA as anticoagulant as it will decrease the lifetime of the calcium sensor.

Long-Term Shutdown of the Analyzer

WARNING

Risk of biohazardous contamination if you have skin contact with the Sensor Cassette. The Sensor Cassette could contain residual biological material and must be handled with care. Clean up spills immediately. Dispose of the Sensor Cassette in accordance with your local regulations and acceptable laboratory procedures.

 WARNING

Risk of biohazardous contamination if you have skin contact with the Solution Pack and its contents. The Solution Pack could contain residual biological material and must be handled with care. Clean up spills immediately. Dispose of the Solution Pack in accordance with your local regulations and acceptable laboratory procedures.

Parameter Settings

 CAUTION

Risk of incorrect measurement results. Changes made to the offset and/or slope of parameters will have an effect on results and change some performance characteristics. If you do not want the changes to have an effect on QC results too, make sure the Apply parameter corrections to QC check button is deselected in the Miscellaneous setup screen.

Laser Related Hazards

 WARNING

Risk of personal injury. The laser beam can cause eye damage if viewed either directly or indirectly from reflective surfaces (such as a mirror or shiny metallic surfaces). To prevent eye damage, avoid direct exposure to the laser beam. Do not view it directly or with optical instruments. This product uses a Class I laser to read the bar codes of the samples. This product complies with 21CFR subchapter J, section 1040.12 and IEC 60825-1.

Cleaning

WARNING

Risk of infection. Make sure that you wear gloves during replacement and maintenance procedures.

Use universal precautions when working with pathogenic materials. Means must be available to decontaminate the instrument and to dispose of biohazardous waste.

Replacements

WARNING

Risk of infection. The used Inlet Gasket Holder could have been in contact with biologic samples and must be handled as potentially infectious.

Use universal precautions when working with pathogenic materials. Means must be available to decontaminate the instrument and to dispose of biohazardous waste.

WARNING

Risk of infection. The used Inlet Probe could have come in contact with biologic samples and must be handled as potentially infectious.

Use universal precautions when working with pathogenic materials. Means must be available to decontaminate the instrument and to dispose of biohazardous waste.

WARNING

Risk of infection. The used Inlet Connector Gasket has been in contact with biologic samples and must be handled as potentially infectious.

Use universal precautions when working with pathogenic materials. Means must be available to decontaminate the instrument and to dispose of biohazardous waste.

WARNING

Risk of infection. The used Inlet Module has been in contact with biologic samples and must be handled as potentially infectious.

Use universal precautions when working with pathogenic materials. Means must be available to decontaminate the instrument and to dispose of biohazardous waste.

 WARNING

Risk of personal injury. A shock hazard exists if the power cable is connected. Turn off the Cytometer and disconnect the primary power cable before performing these procedures.

 CAUTION

Risk of instrument damage. Ensure that you use the correct fuse type and rating (2 x F2.5A, 250 V) to prevent the risk of fire and/or improper instrument operation.

Installation Environment Validation

 CAUTION

Risk of instrument damage. Place the instrument on a level surface. Failing to do so places the system in danger of toppling and can result in damage. Take all necessary precautions throughout the process of storing or transporting the instrument.

 CAUTION

Risk of instrument damage. Use an outlet with power fail protection to prevent damage to the Sensor Cassette during extended power outages.

 CAUTION

Risk of instrument damage. Ensure that the power source is properly grounded. Improper grounding can cause electric shock and damage the system.

 CAUTION

Risk of instrument damage and/or erroneous results. To ensure reliability, the system must be operated in the specified environment, within the required temperature and humidity ranges. If the ambient temperature or humidity level falls outside the ranges mentioned above, use appropriate air conditioning.

Basic Troubleshooting


IMPORTANT This section covers basic troubleshooting for the Vi-CELL MetaFLEX instrument. For additional troubleshooting, refer to the Vi-CELL MetaFLEX Bioanalyzer Instruction for Use manual.

Fluid Transport Error

Probable Cause

- Air introduced with the sample
- Inhomogeneous solution
- Unable to pump solution
- Leak detected
- Cal backlog error


1 Clean the inlet.

Start the video guidance and follow the instructions on the screen.

OR

- Select **Analyzer status > Other activities > Inlet check**. The Inlet cleaning screen appears.
- Start the guided action.
- Check if the probe is bent. If bent, replace it.

- Dampen a lint-free tissue with tap water.


- e. Remove any deposits from the inlet gasket.
 - f. Check the position of the inlet gasket in the Inlet Gasket Holder.
 - g. Select **Action Completed**.
 - h. Select **Test again**.
 - i. Perform a calibration and/or QC, if necessary.
- If the problem persists, proceed to Step 2.

2 Replace Inlet Gasket Holder and Connector Gasket.


Start the video guidance and follow the instructions on the screen.

OR

- a. Select **Analyzer status > Other activities > Inlet check**. The Inlet replacement screen appears.
- b. Lift the inlet handle to the capillary position.


- c. Pull out the Inlet Gasket Holder.

- d. Put a new Inlet Gasket Holder in place.
 - e. Pull off the inlet module.
 - f. Pull out the Connector Gasket.
 - g. Dampen a new Connector Gasket with water.
 - h. Wiggle the Connector Gasket fully into place.
 - i. Put on the inlet module and select **Action Completed**.
 - j. Select **Done**.
 - k. Perform a calibration and/or QC, if necessary.
- If the problem persists, proceed to Step 3.

IMPORTANT The lid on the Solution Pack must be parallel to the top of the compartment that it is placed in.

3 Check the Solution Pack.

Start the video guidance and follow the instructions on the screen.


OR

- a. Select **Analyzer status > Consumables > Replace > Replace solution pack**. The Replacing Solution Pack screen appears.

Remove the Solution Pack when ejected.

- b. Verify that the two blue tabs are in the tab holes.

If not, press the lid down with both hands until the tabs click into the holes.

- c. Reinsert the Solution Pack and follow the instructions on the screen. Select **Test again**, if available.

- d. Perform a calibration and/or QC, if necessary.

If the problem persists, proceed to Step 4.

4 Flush the analyzer.

Start the video guidance and follow the instructions on the screen.

OR


- a. Select **Analyzer status > Consumables > Replace > Replace solution pack**. The Replacing Solution Pack screen appears.

Remove the Solution Pack when ejected.

- b. Fill a Flush Device with 50% tap water and 50% air.

IMPORTANT Ensure that the Flush Device is securely connected.

- c. Connect the Flush Device to the waste connector.

- d. Hold the Flush Device with the plunger upwards and inject a little tap water.

NOTE It may be necessary to pull and push the syringe plunger back and forth a couple of times.

- e. Hold the Flush Device with the plunger downwards and inject a little air.

NOTE It may be necessary to pull and push the syringe plunger back and forth a couple of times.

- f. Repeat Steps d and e several times, finishing with tap water.
- g. Push the Inlet Gasket Holder back to expose the probe. Fluidic path is cleaned when an unbroken stream of water comes out of the inlet.
- h. Reinsert the Solution Pack and follow the instructions on the screen. Select **Test again**, if available.
- i. Perform a calibration and/or QC, if necessary.

Single Parameter Error

Probable Cause

- Sensitivity error
- Status error
- Response error
- Impedance error

IMPORTANT Only use the table below if a Fluid transport error has not been reported. If a Fluid Transport error has been reported, refer to [Fluid Transport Error](#).

Use the following table when you troubleshoot single parameter errors after Sensor Cassette start-up.

Parameter	Error reading (Same error on consecutive calibrations or QCs)	Ranges	Corrective Action
pH, K ⁺ , Na ⁺ , Ca ²⁺ , Cl ⁻	Sensitivity error	[85-105%] Cl: [75-105%]	Replace the Sensor Cassette.
	Status error	[-50-250 mV] Cl: [-50-100mV]	Replace the Sensor Cassette.
	Impedance error	N/A	Replace the Sensor Cassette.
	Response error	N/A	Replace the Sensor Cassette.
pCO ₂	Sensitivity error	[60-105%]	Replace the Sensor Cassette.
	Status error	[-50-250mV]	Replace the Sensor Cassette.
	Impedance error	N/A	Replace the Sensor Cassette.
	Response error	N/A	Replace the Sensor Cassette.
pO ₂	Sensitivity error	[85-110%]	Replace the Sensor Cassette.
	Status error	[-20-20]	1. Replace the Solution Pack. 2. Replace the Sensor Cassette.
	Substrate thickness	N/A	Replace the Sensor Cassette.
	pO ₂ QC error	N/A	Replace the Solution Pack.
Glu, Lac	Sensitivity error	[100-2000 pA]	Replace the Sensor Cassette.
	Status error	[0-3000 pA]	1. Replace the Solution Pack. 2. Replace the Sensor Cassette.
	Response error	N/A	Replace the Sensor Cassette.

If none of these actions solved your problem, [contact us](#).

Contact Us

If you have any questions, contact our Customer Support Center.

- Worldwide, find us via our website at
www.beckmancoulter.com/customersupport/support.
- In the USA and Canada, call us at 1-800-369-0333.
- Outside of the USA and Canada, contact your local Beckman Coulter Representative.

Trademarks

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

All other trademarks, service marks, products, or services are trademarks or registered trademarks of their respective holders.

Beckman Coulter, Inc.
250 S. Kraemer Blvd.
Brea, CA 92821 U.S.A.

Beckman Coulter Eurocenter S.A.
22, rue Juste-Olivier
Case Postale 1044
CH - 1260 Nyon 1, Switzerland
Tel: +41 (0) 22 365 36 11

© 2022 Beckman Coulter, Inc.
All Rights Reserved

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