Stem-Kit Reagents

REF IM3630 50 tests

TABLE OF CONTENTS

English	2
Français (FR)	14
Deutsch (DE)	27
Italiano (IT)	40
Español (ES)	53
Português Portugal (PT-PT)	66
Dansk (DA)	79
Svenska (SV)	91
Norsk (NO)	103
Ελληνικά (EL)	115
Magyar (HU)	129
Polski (PL)	142
Čeština (CZ)	155
Türkçe (TR)	168
Русский (RU)	180
eesti keel (ET)	195
Română (RO)	207
Latviešu (LV)	220
Українська (UK)	232
Português Brasil (PT-BR)	246
APPENDIX	259
REFERENCES	271

Stem-Kit Reagents

REF IM3630 50 tests

For *In Vitro* Diagnostic Use

INTENDED USE

Stem-Kit Reagents consist of a two-color fluorescent (FITC, PE) murine monoclonal antibody reagent, a two-color murine fluorescent (FITC, PE) isoclonic control, an absolute count reagent, a cell viability reagent, and a lysing reagent. It is intended "For In Vitro Diagnostic Use" for the simultaneous identification and enumeration of CD45+ and dual-positive CD45+ CD34+ cell population percentages and absolute counts in biological specimens by flow cytometry. Biological specimens include fresh normal or mobilized peripheral whole blood, and fresh or thawed apheresis products, cord blood and bone marrow. Cell population measurements may be obtained using either an automated method or a manual method for gating and analysis.

Refer to this Stem-Kit Reagents package insert for complete instructions if using the manual method. Refer to the stemCXP System Guide provided with the stemCXP 2.0 software kit and stemCXP version 2.0, WIN 7 for complete instructions if using the automated method.

SUMMARY

Hematopoietic Progenitor Cells (HPC)

The Stem-Kit Reagents are designed to identify the human hematopoietic progenitor cells (HPC) using the following criteria: "True CD34+ cells (a) express CD34 antigen, (b) express CD45 antigen with staining intensity characteristic of blast cells (i.e., readily detectable but at lower levels than lymphocytes and monocytes), and (c) exhibit low side-angle and low to intermediate forward angle light scatter characteristics of blast cells" (1). Exclusion of dead CD34+ HPC from viable CD34+ HPC enumeration is achieved using the 7-AAD Viability Dye (2).

CD34

The CD34 antigen is a family of differentially glycosylated type I transmembrane single chain glycoproteins (3, 4, 5). It is expressed on virtually all hematopoietic precursor cells including multipotent stem cells. CD34 epitopes have been separated into three classes based upon sensitivity to neuraminidase and to *Pasteurella haemolytica* glycoprotease (6). The CD34 monoclonal antibody (mAb) used in Stem-Kit Reagents is named 581. It specifically recognizes a class III epitope of the CD34 molecule (7) as defined by enzymatic cleavage patterns (6, 8). This 581 monoclonal antibody was assigned to CD34 during the 5th HLDA Workshop on Human Leukocyte Differentiation Antigens in Boston, Massachusetts, USA, in 1993 (7).

CD45

The CD45 antibody recognizes members of the CD45 family of pan leukocyte antigens with molecular weights of 180, 190, 210, and 220 kD (9, 10, 11). It is also known as the leukocyte common antigen (LCA). CD45 antigen is expressed on every type of hematopoietic cell except mature erythrocytes and their immediate progenitors (11, 12). It has not been detected in differentiated nonhematopoietic tissue (12, 13, 14, 15). The CD45 mAb used in Stem-Kit Reagents is named J33. It was assigned to CD45 during the 3rd HLDA Workshop on Human Leukocyte Differentiation Antigens in Oxford, England, in 1986 (16).

7-AAD

7-Amino-Actinomycin D (7-AAD) Viability Dye is an analog of actinomycin D that contains a substituted amino group at the 7 position of the chromophore. 7-AAD Viability Dye intercalates between cytosine and guanine bases of the DNA (15). Actinomycins are biologically active compounds containing a 2-amino-4, 6 dimethylphenoxazine-3 chromophore and two cyclic pentapeptides (17). Actinomycins form stable complexes with double-stranded deoxyribonucleic acid (DNA) but neither with double-stranded ribonucleic acid (RNA), nor with RNA-DNA hybrids or with single-stranded DNA or RNA (15, 18).

The spectral properties of 7-AAD Viability Dye make this molecule particularly suitable for flow cytometry analysis (17). The maximum absorption of the 7-AAD / DNA complex is situated in the green spectral region and thus is suitable for an argon laser equipped-cytometer (excitation wavelength of 488 nm) (15). The deep red fluorescence emission of 7-AAD Viability Dye (635 to 675 nm) makes the use of this probe easier in combination with fluorescein isothiocyanate (FITC) and phycoerythrin (PE) conjugated antibodies because, in contrast to propidium iodide (PI), the 7-AAD / DNA complex shows minimal spectral overlap with FITC and PE (17).

Apoptotic, necrotic cells, and / or damaged cells are sources of interference in the analysis of viable cells by flow cytometry. Nonviable cells shall be evaluated and discriminated by focusing on 7-AAD-positive staining when viable cells remain unstained (negative) (19).

Flow Cytometer Systems

Flow cytometer systems (FC500) apply the principles of flow cytometry to analyze biological specimens in order to identify various cell populations determined by the specific monoclonal antibodies and fluorochromes used. The present test uses the FC500 flow cytometer systems equipped with CXP System Software to analyze samples.

PRINCIPLE

This test depends on the ability of a monoclonal antibody to bind to the surface of cells expressing discrete antigenic determinants. The CD45-FITC / CD34-PE reagent in Stem-Kit Reagents is a combination of two murine monoclonal antibodies; each conjugated to a specific fluorochrome and specific for a different cell surface antigen. Specific cell surface staining is accomplished by incubating duplicate samples of a biological specimen with the two-color CD45-FITC / CD34-PE reagent.

An additional test of the same sample is stained with the CD45-FITC / IsoClonic Control-PE reagent to check the non-specific binding of the CD34-PE monoclonal antibody. 7-AAD Viability Dye, a nucleic acid dye that binds to accessible base pairs (cellular DNA), is used to distinguish between viable and nonviable cells. The red blood cells in each sample are lysed with NH₄Cl Lysing Solution prepared at 1X. Stem-Count Fluorospheres or Flow-Count Fluorospheres (the absolute count reagent) are added, and the remaining cells analyzed by flow cytometry. Process control and quality control, for the assay, are provided by the instrument and sample quality control reagents: Flow-Check Fluorospheres, Flow-Set Fluorospheres, appropriate compensation reagents, and Stem-Trol Control Cells.

EXAMPLES OF CLINICAL APPLICATIONS

A chronology of hematopoietic progenitor cells (HPC) scientific reports shows early studies relied upon nucleated cell counts for quantification. Subsequent studies involving colony forming unit (CFU) assays provided rough approximations of progenitor cell populations. The time intensive CFU assays were followed by the discovery of the CD34 cell surface molecule. CD34 was found to be present on the surface of early hematopoietic cells and hematopoietic colony-forming cells in bone marrow and peripheral blood (20). Monoclonal antibodies directed against the CD34 molecule were then developed and utilized in flow cytometry. CD34 determinations by flow cytometry are now widely accepted as a means to assess hematopoietic progenitor and stem cells.

In bone marrow of normal individuals, 1 to 3% of the cells expressing the CD34 antigen are capable of reconstituting long-term multi-lineage hematopoiesis (3, 20). In peripheral blood of normal individuals, less than 0.1% of circulating cells express the CD34 antigen (1, 21). The use of Stem-Kit reagents provides a simple, rapid, and sensitive flow cytometric method to enumerate hematopoietic progenitor and stem cells in peripheral blood, leukopheresis products, bone marrow, and cord blood.

REAGENT CONTENTS

Stem-Kit Reagents contain the following:

- CD45-FITC / CD34-PE Reagent (45/34) – 1 vial.
- CD45-FITC / IsoClonic Control-PE Reagent (45/CTRL) – 1 vial.
- Stem-Count Fluorospheres – 1 vial.
- 7-AAD Viability Dye – 1 vial.
- 10X NH₄Cl Lysing Solution – 2 vials.

1. These reagents must be used as a set as supplied in the kit.

Flow-Count Fluorospheres, which are equivalent to Stem-Count Fluorospheres, can be used for the absolute count determination if using the specific assayed concentration.

2. The CD45-FITC / CD34-PE is a ready-to-use two-color reagent provided in 2.2 mL of phosphate-buffered saline (PBS) containing 2 mg/mL bovine serum albumin (BSA) and 0.1% sodium azide as a preservative.
3. The CD45-FITC / IsoClonic Control-PE is a ready-to-use two-color reagent provided in 1.2 mL of PBS containing 2 mg/mL BSA and 0.1% sodium azide as a preservative. The CD45-FITC / IsoClonic Control-PE reagent is a mixture of unconjugated and PE-conjugated (CD34 identical) 581 monoclonal antibody. This mixture acts as an appropriate, isomorphologic control to check the nonspecific binding of the CD34-PE mAb on each and every sample.

NOTE: The CD45-FITC / CD34-PE and the CD45-FITC / IsoClonic Control-PE reagents are not to be separated. When the tests and the controls are performed, care should be taken that both reagents have the same lot number on the vial label.

Summary of antibodies characteristics

Specificity	CD45	CD34
Clone	J33	581
Hybridoma	NS1 x Balb/c	NSO x Balb/c
Immunogen	Human LAZ 221 ALL cell line	Human KG-1a cell line
Ig Chain	IgG1, κ light chains	IgG1, κ light chains
Species	Mouse	Mouse
Conjugation	Fluorescein isothiocyanate (FITC)	R Phycoerythrin (PE)
Molar ratio	5 – 8 moles FITC / mole Ig	0.5 – 1.5 moles PE / mole Ig

Specificity	CD45	CD34
Excitation wavelength	488 nM (blue-green)	488 nM (blue-green)
Fluorescence emission	504 to 541 nm (mainly green)	568 to 590 nm (mainly orange)

- The Stem-Count Fluorospheres or Flow-Count Fluorospheres consist of 10 µm (nominal diameter) polystyrene fluorescent microspheres in a 20 mL aqueous suspension medium, containing a surfactant and 1% formaldehyde. Each fluorosphere contains a dye that has a fluorescence emission range of 525 nm to 700 nm when excited at 488 nm. The Stem-Count Fluorospheres, or Flow-Count fluorospheres, assayed concentration (fluorospheres / µL) is derived from multiple replicate analyses using a COULTER particle size analyzer. The lot specific value is stated on the vial label.
- The 7-AAD Viability Dye is provided in 99% phosphate-buffered saline (PBS: 0.01 M sodium phosphate, 0.145 M sodium chloride, pH 7.2), with 1% dimethyl sulfoxide (DMSO).
- The NH₄Cl Lysing Solution is provided as a 10X concentrated, 20 mL solution. It contains a mixture of ammonium chloride (NH₄Cl), potassium bicarbonate (KHCO₂), Disodium Ethylene Diamine Tetraacetic Acid (EDTA).

WARNING AND PRECAUTIONS

- Monoclonal antibody preparations contain 0.1% sodium azide. Sodium azide, under acidic conditions, yields hydrazoic acid, an extremely toxic compound. Disposal into sewage systems should be accompanied with large quantities of running water to avoid dangerous buildup in metal piping and potential explosions.
- Specimens, samples and all material coming in contact with them should be handled as if capable of transmitting infection and disposed of with proper precautions.
- Never pipet by mouth and avoid contact of samples or reagents with skin and mucous membranes.
- Do not use reagents beyond the expiration dated on the vial label.
- Stem-Count Fluorospheres or Flow-Count Fluorospheres contain 1% formaldehyde. Avoid contact with skin and eyes as formaldehyde may cause irreversible effects to these tissues. Do not breathe vapors. Wear appropriate safety equipment such as gloves and eye protection.
- Do not expose Stem-Kit Reagents to strong light or to heat during storage or use.
- Avoid evaporation and leakage of Stem-Count Fluorospheres or Flow-Count Fluorospheres, or erroneous results may occur. Tightly cap the vials after use and store upright.
- Ensure the fluorospheres are completely resuspended before use. Avoid excessive mixing to minimize the formation of air bubbles. Do not pipet air bubbles or erroneous results may occur.
- Use a calibrated pipet (positive displacement or repeater pipets, preferred) for dispensing samples and Stem-Count Fluorospheres or Flow-Count Fluorospheres, or erroneous results may occur.
- Use pipetting techniques recommended by the pipet manufacturer to ensure accurate and precise pipetting of samples and Stem-Count Fluorospheres or Flow-Count Fluorospheres, or erroneous results may occur.
- Stained and lysed samples must be stored on ice and analyzed on the flow cytometer according to the information in section "Processed Sample stability" otherwise erroneous results can occur.
- Ensure that at least 1,000 Stem-Count Fluorospheres or Flow-Count Fluorospheres are counted, or the automatic conversion of results to absolute counts will not occur.
- Each lot of Stem-Count Fluorospheres or Flow-Count Fluorospheres has a specific concentration of fluorospheres. Ensure the correct assayed concentration is used when determining absolute count results.
- Preparation methods requiring a wash step should not be used due to an unknown loss of cells.
- The use of incubation or mixing times, or temperatures other than those specified, may give erroneous results.
- 7-AAD Viability Dye contains 0.005% 7-AAD. Pure 7-AAD is a potential carcinogen. Although this compound is highly diluted, we recommend to avoid contact with skin and eyes, to wear suitable protective clothing and gloves and appropriate eye / face protection.
- 7-AAD Viability Dye contains 1% Dimethyl Sulfoxide (DMSO). Pure DMSO is an irritant. Although this compound is diluted, we recommend avoiding contact with skin and eyes, to wear suitable protective clothing and gloves, and appropriate eye / face protection.
- As a chromophore, 7-AAD is sensitive to light. Store the 7-AAD Viability Dye within the kit's cardboard packaging container to avoid continuous exposure to light during storage.
- Minimize exposure of mixed samples to light during incubation.
- Use general good laboratory practices when handling reagents.

GHS HAZARD CLASSIFICATION

Stem-Count Fluorospheres

DANGER

H316

Causes mild skin irritation.

H317

May cause an allergic skin reaction.

H350

May cause cancer.

P201

Obtain special instructions before use.

P280

Wear protective gloves, protective clothing and eye/face protection.

P308+P313

IF exposed or concerned: Get medical advice/attention.

P333+P313

If skin irritation or rash occurs: Get medical advice/attention.

P362+P364

Take off contaminated clothing and wash it before use.

Formaldehyde 0.5 - 1.0%

reaction mass of:

5-chloro-2-methyl-4-isothiazolin

-3-one [EC# 247-500-7] and

2-methyl-4-isothiazolin-3-one [EC# 220-239-6](3:1) < 0.05%

7-AAD Viability Dye

WARNING

H316

Causes mild skin irritation.

P332+P313

If skin irritation occurs: Get medical advice/attention.

Dimethyl Sulfoxide 0.5 - 2%

Safety Data Sheet is available at beckman.com/techdocs

STORAGE CONDITIONS AND STABILITY

The unopened reagent is stable until the expiration date stated on the vial when stored at 2 – 8°C in the absence of light. Do not use after the expiration date. Do not freeze. The Stem-Kit Reagents are stable for 30 days after opening of Stem-Count Fluorospheres or Flow-Count Fluorospheres vial. Opened vials must be recapped immediately after use and stored at 2 – 8°C in the absence of light. Do NOT aliquot Stem-Count Fluorospheres or Flow-Count Fluorospheres from its original vial.

IMPORTANT: Risk of erroneous results if leakage of the Stem-Count Fluorospheres or Flow-Count Fluorospheres occurs. Once opened Stem-Count Fluorospheres or Flow-Count Fluorospheres must be stored in an upright position to prevent the possibility of leakage. Reagents should not be used if signs of leakage are observed.

Evidence of Deterioration

Evidence of deterioration may occur in the following situations:

Stem-Kit Reagents, except Stem-Count Fluorospheres

- A change in normal appearance as listed below

Reagent	Normal Appearance
CD45-FITC / CD34-PE	Clear, slightly pink liquid
CD45-FITC / IsoClonic Control-PE	Clear to slightly pink liquid
7-AAD Viability Dye	Clear, slightly pink to reddish liquid
NH ₄ CL Lysing Solution	Clear, colorless liquid

Stem-Count Fluorospheres

- The appearance, during acquisition, of secondary fluorescent populations containing more than 20% of the total population of fluorospheres.
- A major variation in values obtained for counts or samples.

IMPORTANT: Reagents should not be used if there is evidence of deterioration because incorrect results could be obtained. Call your local Beckman Coulter Representative if there is evidence of deterioration.

Reagent Preparation

1. The CD45-FITC / CD34-PE and CD45-FITC / IsoClonic Control-PE two-color reagents are used directly from the vials without any dilution.
2. Stem-Count Fluorospheres or Flow-Count Fluorospheres are used directly from the vial without any dilution. All of the aluminum from the heat seal must be removed from the top of the bottle when first opened. Proper mixings are required. When first handling the Stem-Count Fluorospheres or Flow-Count Fluorospheres vial at the beginning of the day, vortex for 3 – 5 seconds.

Before pipetting from the vial, gently agitate by 3 – 5 vial inversions.

After initial vortexing at the beginning of the day, return the fluorospheres vial to 2 – 8°C until use.

Avoid evaporation and leakage by tightly capping vial after use.

3. 7-AAD Viability Dye is used directly from the vial without any dilution.
4. A 1X NH₄Cl Lysing Solution must be prepared daily: add nine volumes of deionized water to one volume of 10X concentrated NH₄Cl Lysing Solution. Do not refrigerate. Discard the unused 1X NH₄Cl Lysing Solution at the end of the day.

Bring all reagents to room temperature (18 – 25°C) prior to use.

SPECIMEN HANDLING

Specimen Requirements

Specimens handling summary:

Specimen Type	Gating Method	Storage Temperature	Specimen Age
Apheresis	Algorithm	2-8°C	24 Hours
	Manual	2-8°C	24 Hours
Bone Marrow	Algorithm	2-8°C	24 Hours
	Manual	2-8°C	24 Hours
Cord Blood	Algorithm	18-25°C	24 Hours
	Manual	18-25°C	24 Hours
Normal Whole Blood	Algorithm	18-25°C	20 hours
	Manual	18-25°C	20 hours
Mobilized Whole Blood	Algorithm	18-25°C	20 hours
	Manual	18-25°C	20 hours

1. Biological specimens include fresh, normal or mobilized peripheral whole blood, and fresh or thawed apheresis products, cord blood and bone marrow.
2. Manipulated or purified hematopoietic progenitor cells (HPC) can only be analyzed with the manual gating method.
3. Collect apheresis products, peripheral blood specimens, cord blood specimens or bone marrow specimens aseptically into a sterile evacuated blood collection tube with ACD formula A or K3 EDTA anticoagulant.
4. Whole blood and mobilized whole blood specimens collected in anticoagulant tubes can be stored for up to 20 hours at room temperature (18-25 °C) before staining.
5. Cord blood specimens collected in anticoagulant tubes can be stored for up to 24 hours at room temperature (18-25 °C) before staining.
6. Collected Apheresis products and Bone Marrow specimens can be stored for up to 24 hours at 2-8 °C before staining.
7. Prior to sample preparation, ensure that the white blood cell (WBC) concentration is no greater than 30 x 10⁹ WBC cells / L. If necessary, dilute the sample with HBSS or equivalent to reach an average value of 15 x 10⁹ WBC cells / L. Record the dilution factor for the calculation of the final HPC absolute count.
8. Low protein concentrations may lead to Stem-Count Fluorospheres or Flow-Count Fluorospheres clumping and give erroneously high absolute cell counts. Dilution of specimens with a low protein concentration must be performed using a protein-enriched buffer.
9. All specimens should be processed with Stem-Kit Reagents including the 7-AAD Viability Dye.

Procedure for Specimen Processing

Perform tests in duplicate, following “ The ISHAGE Guidelines for CD34+ Cell Determination by Flow Cytometry ” (1).

NOTE: For Stem-Trol Control Cells processing, please refer to the Stem-Trol Control Cells package insert for complete instructions, if using the manual method. Refer to the stemONE System Guide provided with the stemONE System Software or stemCXP System Guide provided with stemCXP System Software for complete instructions, if using the automated analysis method.

1. For each specimen sample, label three tubes as follows:
 - 45/34/7-AAD #1
 - 45/34/7-AAD #2
 - 45/CTRL/7-AAD
 2. Pipet 20 µL of CD45-FITC / CD34-PE Reagent into tubes: 45/34/7-AAD #1 and 45/34/7-AAD #2.
 3. Pipet 20 µL of CD45-FITC / IsoClonic Control-PE Reagent into the tube: 45/CTRL/7-AAD.
 4. Pipet 20 µL of 7-AAD Viability Dye into all tubes.
- IMPORTANT:** Risk of incomplete lysis if blood specimen remains on the top or side of the test tube. Use care when pipetting to prevent the blood from touching the top or side of the test tube. Clean the tube with a cotton swab, if necessary, to remove all traces of blood specimen from the top or side of the test tube.
5. Using a positive displacement pipet or repeater pipet, pipet 100 µL of the patient sample into the tube.
 6. Immediately vortex the tube for 5 seconds.
 7. Incubate all tubes at room temperature (18 – 25°C) for 20 minutes, protected from light.
 8. After incubation, add 2 mL of prepared 1X NH₄Cl Lysing Solution to each tube and vortex immediately for 5 seconds.

NOTE: For details about preparing the Lysing Solution, see section 7.2 above, entitled REAGENT PREPARATION.

IMPORTANT: Risk of erroneous results if air bubbles are pipetted. Excessive mixing of Stem-Count Fluorospheres or Flow-Count Fluorospheres can create air bubbles. Do not excessively mix the Stem-Count Fluorospheres or Flow-Count Fluorospheres, and do not pipet air bubbles into the sample tubes.

9. Incubate at room temperature for 10 minutes, protected from light.
10. Gently mix the Stem-Count Fluorospheres or Flow-Count Fluorospheres by inverting the vial 3 to 5 times before using.
11. Prior to acquisition, and immediately after step 9 using a positive displacement or repeater pipet, pipet 100 µL of Stem-Count Fluorospheres or Flow-Count Fluorospheres into each of the three tubes.
12. Immediately vortex each tube for 5 seconds.

NOTE: For details about preparing the Lysing Solution, see section 7.2 above, entitled REAGENT PREPARATION.

Specimen Staining Summary

Reagents / Samples	45/34/7-AAD #1	Test Tube Label 45/34/7-AAD #2	45/CTRL/7-AAD
CD45-FITC / CD34-PE	20 µL	20 µL	
CD45-FITC / IsoClonic Control-PE			20 µL
7-AAD Viability Dye	20 µL	20 µL	20 µL
Specimen	100 µL	100 µL	100 µL
Vortex – Incubate at room temperature for 20 minutes. Protect from light.			
1X NH ₄ Cl Lysing Solution	2 mL	2 mL	2 mL
Vortex – Incubate at room temperature for 10 minutes. Protect from light.			
Stem-Count Fluorospheres	100 µL	100 µL	100 µL
Vortex			
For Collected Apheresis products, Cord Blood samples and Bone Marrow samples, do not wait more than 60 minutes on ice prior to acquisition.			
For whole blood and mobilized whole blood samples, do not wait more than 45 minutes on ice prior to acquisition.			

Reagents / Samples	45/34/7-AAD #1	Test Tube Label 45/34/7-AAD #2	45/CTRL/7-AAD
Analyze the thawed sample immediately.			
Protect from light			

Processed Sample stability:

- For Collected Apheresis products, Cord Blood samples and Bone Marrow samples, after the addition of Stem-Count or Flow-Count fluorospheres do not wait more than 60 minutes on ice prior to acquisition.
- For whole blood and mobilized whole blood samples, after the addition of Stem-Count or Flow-Count fluorospheres do not wait more than 45 minutes on ice prior to acquisition.
- For Frozen sample acquisition should be done immediately after the addition of Stem-Count or Flow-Count fluorospheres.
- Processed samples must be stored on ice protected from light, and limit the RT exposure as much as possible during acquisition.

CONTENTS

Sodium azide preservative may form explosive compounds in metal drain lines. See NIOSH Bulletin: Explosive Azide Hazard (8/16/76).

To avoid the possible build-up of azide compounds, flush wastepipes with water after the disposal of undiluted reagent. Sodium azide disposal must be in accordance with appropriate local regulations.

MATERIALS NEEDED BUT NOT SUPPLIED WITH KIT:

1. Deionized water.
2. Hank's Balanced Salt Solution (HBSS) or equivalent
3. Stem-Trol Control Cells (Ref. IM3632).
4. Flow-Check Fluorospheres (Ref. 6605359)
5. Flow-Set Fluorospheres (Ref. 6607007)
6. CYTO-COMP Cell Kit (Ref. 6607023)
7. QuickCOMP 2 Kit (Ref. 177018)
8. As an alternative to Stem-Count Fluorospheres: Flow-Count Fluorospheres (Ref. 7547053)
9. Ice water bath
10. Plastic test tubes (12 x 75 mm)
11. Calibrated repeater pipet (20 μ L, 100 μ L, 2 mL) and tips or calibrated positive displacement pipets (20 μ L, 100 μ L, 2 mL) and tips
12. Calibrated standard pipette (20 μ L, 100 μ L, 2 mL) and tips.
13. Vortex mixer.
14. Timer.
15. Flow cytometer.
16. stemCXP 2.0 software kit, stemCXP version 2.0, WIN 7, ONLY for automated analysis of samples on FC500 flow cytometer system equipped with CXP Software (Version 2.0).

RESULTS

Analysis Protocol Setup

The flow cytometer must be equipped to detect forward scatter, side scatter and four fluorescence channels. For the FL3 channel (for Stem-Count or Flow-Count Fluorospheres monitoring) use a 620 nm band pass filter. For the FL4 channel (for 7-AAD Viability Dye monitoring) use a 675 nm long pass filter.

Ensure that the flow cytometer is properly aligned and standardized for fluorescence intensity.

Ensure that the flow cytometer is properly adjusted for fluorescence compensation.

NOTE: For specimen analysis, all region positions are approximate until the end of acquisition of the 45/34/7-AAD #1 tube. Regions are optimized according to the following protocol and then NOT adjusted for the second and third tubes in the panel.

A series of 8 histograms and gating scheme must be followed for specimen analysis (see below). The same series of 8 histograms and gating scheme must be followed when using CD34+ control cells. Regions A, B, C, and D (see the heading "Region Creation" for region definition) shall be adapted to include the CD34+ control cells characteristic cluster.

Histogram Creation

Create histograms as follows:

1. Create Histogram 1 as FL1 CD45-FITC vs. Side Scatter.
2. Create Histogram 2 as FL2 CD34-PE vs. Side Scatter.
3. Create Histogram 3 as FL1 CD45-FITC vs. Side Scatter.
4. Create Histogram 4 as Forward Scatter vs. Side Scatter.
5. Create Histogram 5 as FL1 CD45-FITC vs. FL2 CD34-PE.
6. Create Histogram 6 as Forward Scatter vs. Side Scatter.
7. Create Histogram 7 as Time vs. FL3 Stem-Count or Flow-Count Fluorospheres.
8. Create Histogram 8 as FL4 7-AAD Viability Dye vs. Side Scatter.

Histograms 1 to 4 are intended to characterize CD34+ HPC, a process that may be delayed until the analysis step. These first four histograms are set up according to the ISHAGE Guidelines for CD34+ Cell Determination by Flow Cytometry (1, 2). Histograms 5 to 7 are intended to monitor parameters that are of importance during the acquisition step. These include the forward scatter discriminator, the number of CD45+ events to be collected, and the correct fluorosphere singlets accumulation.

Histogram 8 is intended to discriminate and analyze viable events from nonviable events.

Region Creation

Create regions as follows:

1. Histogram 1 – Create a rectilinear Region A on Histogram 1 to include all CD45+ leukocytes and eliminate platelets, red blood cell debris, and aggregates.
2. Histogram 1 – Create an amorphous Region E on Histogram 1 to include the lymphocytes (bright CD45, low side scatter).
3. Histogram 2 – Create a rectilinear Region B on Histogram 2 to include all CD34+ events. Set a stop count of 75,000 events (CD45+ events) for Histogram 2.
4. Histogram 3 – Create an amorphous Region C on Histogram 3 to include all clustered CD45dim events.
5. Histogram 4 – Create an amorphous Region D on Histogram 4 to include all clustered events with intermediate side scatter and intermediate to high forward scatter.
6. Histogram 5 – Create a Quadstat Region I on Histogram 5 to verify the lower limit of CD45 expression on CD34+ events.
7. Histogram 5 – Create an amorphous Region H on Histogram 5 to surround all Stem-Count or Flow-Count Fluorospheres, including doublets. The region H should be located at the top right corner of Histogram 5.

NOTE: Be sure that Region H is drawn as an AMORPHOUS region.

8. Histogram 5 – Create rectilinear Region K (named "listgate –") on Histogram 5. Set Region K boundaries at the first log decade of the two parameter histogram. This Region K allows the elimination of double negative events during the acquisition in order to analyze 75,000 relevant CD45+ events.
9. Histogram 6 – Copy Region D from Histogram 4 as Region F in Histogram 6.
10. Histogram 7 – Create a rectilinear Region G on Histogram 7 to include the Stem-Count Fluorospheres or Flow-Count Fluorospheres singlets only. Region G can be labeled as "CAL" to allow automatic calculation of absolute numbers of CD34+ HPC (refer to the instrument manual for further details).
11. Histogram 8 – Create a rectilinear Region J on Histogram 8 to include viable leukocytes (7-AAD Viability Dye negative events).

Gate Creation

Create gates as follows:

1. Histogram 1 – Assign "J" to Histogram 1 to display all viable events. If 7-AAD is not used assign "-H" region to Histogram 1 to display all events excluding all Stem-Count Fluorospheres or Flow-Count Fluorospheres including doublets.

2. Histogram 2 – Assign “ AJ ” (“ A ” if 7-AAD is not used) to Histogram 2 to display viable (total if 7-AAD is not used) CD34+ events.
3. Histogram 3 – Assign “ A ” and “ B ” and “ J ” (“ A ” and “ B ” if 7-AAD is not used) to Histogram 3 to display all viable (total if 7-AAD is not used) CD45+ / CD34+ events.
4. Histogram 4 – Assign “ A ” and “ B ” and “ C ” and “ J ” (“ A ” and “ B ” and “ C ” if 7-AAD is not used) to Histogram 4 to display viable (total if 7-AAD is not used) CD45+ / CD34+ clustered events with low side scatter while excluding platelets and debris. Events from Region ABCDJ are real viable CD34+HPC.
5. Histogram 5 – Ungated to display all events.
6. Histogram 6 – Assign “ EJ ” (“ E ” if 7-AAD is not used) to Histogram 6 to display viable (total if 7-AAD is not used) lymphocytes as a visual check on the discriminator.
7. Histogram 7 – Assign “ H ” to Histogram 7 to display all Stem-Count Fluorospheres or Flow-Count Fluorospheres including doublets.
8. Histogram 8 – Assign “ A ” to Histogram 8 to display CD45+ events.

Flow Cytometry Setting

1. Verify that the flow cytometer is properly aligned and optimized for light scatter and fluorescence intensities of the FL1, FL2, FL3, and FL4 detectors according to manufacturer's and laboratory guidelines. Verify that color compensation is set for standard operation. Refer to the instrument's manual for further instructions.
2. Vortex test tubes for 5 seconds.
3. Perform data acquisition on the flow cytometer. A minimum of 75,000 CD45+ events must be analyzed.
4. For specimen analysis: Adjust the discriminator and regions by analyzing sample tube 45/34/7-AAD #1. Save the protocol, and analyze the 45/34/7-AAD #2 and the 45/CTRL/7-AAD tubes with no further adjustments.

Analysis Example

The histograms shown in the APPENDIX are displayed in an ascending number order as displayed on the protocol.

Calculation of CD34+ HPC for Specimen

The number of CD34+ HPC is calculated as the average of the region D statistical results for the two replicate samples (i.e. 45/34/7-AAD #1 and 45/34/7-AAD #2).

The CD34+ HPC results are considered valid if the results obtained by the CD45-FITC / IsoClonic Control-PE reagent is less than 10% of the CD34+ HPC results.

1. Average the results obtained from the 45/34/7-AAD #1 and 45/34/7-AAD #2 tubes for each specimen (i.e. average of Region D statistics printout results). The number of CD34+ HPC may fall within 10% of the mean for the duplicate samples. The value obtained with the 45/CTRL/7-AAD control tube must represent less than 10% of the average value obtained from the tests tubes to validate the results

NOTES: If the mean for the absolute count obtained for the duplicate CD34+ HPC falls outside the accepted range (i.e. exceeds 10%), determine whether the pipetting has been correctly performed, especially for the Stem-Count Fluorospheres or Flow-Count Fluorospheres. Visually check the flow cytometer reports and the quality of the staining and lysing processes in your test tubes. Verify the way the region boundaries were set.

If the value obtained with the 45/CTRL/7-AAD falls outside of the range of accepted values (i.e. exceeds 10% of the absolute count mean obtained for the duplicate), determine whether the pipetting has been correctly performed, especially for the Stem-Count Fluorospheres or Flow-Count Fluorospheres. Visually check the flow cytometer reports and the quality of the staining and lysing processes in your test tubes. Verify the way the region boundaries were set.

If necessary, repeat the acquisition step with the addition of a BLANK tube after the first duplicate tube 45/34/7-AAD #1 and #2 and before the 45/CTRL/7-AAD tube. If necessary, repeat the whole preparation with new test tubes adding the BLANK tube. Acquire and analyze both the new and the previous series of tubes. Compare the two results. If there is no agreement, call your local Beckman Coulter Representative.

2. If the sample has been diluted, the result obtained above MUST be multiplied by the appropriate dilution factor (Dil). The final result obtained is the absolute CD34+ HPC count (cells / μ L).

$\text{CD34+ HPC Absolute Count (cells / } \mu\text{L)} = \text{CD34+ HPC Absolute Count Printout Result (cells / } \mu\text{L)} \times \text{Dil}$

3. For apheresis packs, the total number of CD34+ HPC per pack can be calculated from the following formula:

$$\text{CD34+ HPC Total Number per Pack} = \text{CD34+ HPC Absolute Count (cells / } \mu\text{L)} \times 10^6 \times \text{Pack Volume (L)}$$

Where:

4. For percentage of CD34+ HPC within the CD45+ leukocytes value calculation, refer to the following formula :

$$\% \text{ CD34+ HPC} = \frac{\text{CD34+ HPC Absolute Count (cells / } \mu\text{L)} \text{ (i.e. region D Printout Results)}}{\text{CD45+ Leukocytes Absolute Count (cells / } \mu\text{L)} \text{ (i.e. Region A Printout Results)}} \times 100$$

QUALITY CONTROL

Ensure that the flow cytometer is properly aligned according to the manufacturer's recommendations as listed in the instrument manuals. Target intensities have been established to optimize light scatter and fluorescence instrument settings for analyzing samples stained with the CD45-FITC / CD34-PE monoclonal antibody reagent and CD45-FITC / Isoclonic Control-PE reagent.

The fluorochromes fluorescein isothiocyanate (FITC), phycoerythrin (PE), and 7-Amino Actinomycin D dye (7-AAD) emit at different wavelengths, but they do have some spectral overlap that must be corrected by electronic compensation. Optimal compensation levels can be established by analyzing, in a dualparameter histogram, a whole blood sample spiked with Stem-Trol Control Cells stained with the 7-AAD Viability Dye and appropriate compensation reagents when using CXP Software. In this analysis, adjustments are made to ensure that no staining occurs in the dual-color quadrant with any individual fluorochrome.

7-AAD Viability Dye (7-Amino Actinomycin D) is used to identify nonviable cells by flow cytometric analysis to allow enumeration of the viable cells of interest. The dye is taken up by dead cells (bright events) that are then excluded from viable cells (negative events) by gating (see Histogram 8).

Before samples are analyzed, an assayed control (Stem-Trol Control Cells) should be stained to verify antibody reactivity.

PERFORMANCE

The performance characteristics provided are based on Stem-Kit Reagents used with an COULTER EPICS XL/XL-MCL flow cytometer with System II Software. For examples of performance on a FC 500 with CXP Software, please refer to the stemCXP System Guide.

Expected Values for Reference Ranges

Normal whole-blood specimens (n = 117) were collected from a population of males (n = 58) and females (n = 59) to establish a normal reference range. This study was performed at an external evaluation site. Replicate tests for each specimen were prepared according to the Stem-Kit Reagents package insert. Manual gating and analysis was carried out to obtain CD45+, CD34+, and CD34% values.

The statistical analysis suggests a non-Gaussian distribution of the collected data. Based on these findings, a nonparametric (robust) approach was used to analyze the data. CD45+, CD34+, and CD34% values were obtained. Values are expressed in terms of absolute counts (cells / μL) for CD45+ and CD34+ and percentage of the total leukocyte count (CD34%) as shown in the table below.

Normal Whole Blood: Reference Ranges

Manual Method (Robust Analysis)	n	CD45+ (cells / μL)				CD34+ (cells / μL)				CD34+ (%)			
		Min	Max	Mean	SD	Min	Max	Mean	SD	Min	Max	Mean	SD
Manual	117	2919.5	11240.0	6469.2	1635.2	0.50	6.50	2.36	1.14	0.015	0.600	0.052	0.056

These values are intended to be representative only. Each laboratory should establish its own expected values from the local population of normal donors.

Linearity

Replicate measurements were made at each of 10 serial dilutions for 2 types of samples: CD34+-spiked peripheral whole-blood and fresh mobilized peripheral whole-blood. A dynamic range of 0 to 2,000 cells / μL and a low-end sensitivity range of 0 to 150 cells / μL were achieved for CD34+ HPC concentrations. Samples were prepared according to the Stem-Kit Reagents package insert and analyzed on a COULTER EPICS XL/XL-MCL flow cytometer. Values are expressed as absolute counts (cells / μL). Results were analyzed in terms of linear regression analyses for expected versus recovered absolute counts for the mean of the replicates less endogenous CD34+ cells. See the following tables.

The slopes of the regression lines provided with the plots for total CD34+ HPC absolute counts show the mathematical measure of linearity. This is supported by minimal potential assay bias as evidenced by the y-intercepts for the analyses. Expected absolute counts demonstrate high correlation with the recovered absolute counts.

CD34+ Spiked Whole Blood Results

See Appendix for graph.

Peripheral (Mobilized) Blood Results

See Appendix for graph.

Accuracy of Method (Tables) (Expected vs. Recovered Results)

Accuracy of Method is defined as the degree of agreement between expected and recovered results for both manual and automated gating and analysis methods. Fresh normal peripheral whole blood was spiked with purified CD34+ cells to a level of 1875 cells / μL . Cell counts were confirmed by a Malassez slide hemacytometer and serially diluted to different concentrations. Expected results are compared to recovered results for the mean of duplicate measurements using automated and manual gating and analysis for each dilution level. Testing for slope and ordinates at the intercept demonstrate no statistical difference between the automated gating and analysis and the expected results, nor between the manual gating and analysis and expected results.

Accuracy of Method - CD34+ (%)

See Appendix for graph.

Accuracy of Method - CD34+ Absolute Count (cells / μL)

See Appendix for graph.

Accuracy of Method - Total Leukocytes - CD45+ (cells / μL)

See Appendix for graph.

Accuracy of Method (Automated vs. Manual Results)

The degree of accuracy between the automated and manual methods of the stemONE System was studied at four external evaluation sites using fresh peripheral whole blood and fresh or thawed apheresis, cord blood, and bone marrow specimens. The confidence intervals (CI) were calculated with 95% limits.

Fresh peripheral whole blood (n = 59 from three sites), fresh apheresis products (n = 51 from three sites), fresh cord blood (n = 29 from two sites), fresh bone marrow (n = 33 from two sites), thawed apheresis products (n = 23 from one site), thawed cord blood (n = 31 from two sites), and thawed bone marrow (n = 20 from one site) specimens were prepared according to the Stem-Kit Reagents package insert and analyzed on COULTER EPICS XL/XL-MCL flow cytometers. Automated and manual gating and analysis were carried out to obtain CD45+, CD34+, and CD34% values.

Accuracy of Methods (Figures)

Fresh Peripheral Whole Blood Specimens (n = 59)

See Appendix for graph.

Fresh Apheresis Specimens (n = 51)

See Appendix for graph.

Fresh Cord Blood Specimens (n = 29)

See Appendix for graph.

Fresh Bone Marrow Specimens (n = 33)

See Appendix for graph.

Thawed Apheresis Specimens (n = 23)

See Appendix for graph.

Thawed Cord Blood Specimens (n = 31)

See Appendix for graph.

Thawed Bone Marrow Specimens (n = 20)

See Appendix for graph.

Precision

Interlaboratory Precision

Studies were performed within one clinical facility by three operators, each using a different flow cytometer. Three specimens representing the CD34+ dynamic range (0 – 2000 cells / μL) were identified:

Low range:	fresh, mobilized peripheral blood
Mid range:	apheresis product
High range:	apheresis product

These specimens were split into three aliquots, one aliquot of each level (low, mid, high) was provided to each of the three operators. Samples were prepared (n = 10) from each aliquot according to the Stem-Kit Reagents package insert. Each operator performed analysis on an EPICS XL flow cytometer using automated and manual gating and analysis to obtain CD45+, CD34+, and CD34% values. Values are expressed in terms of percentage (%) of the total leukocyte count and as absolute counts (cells / μ L).

Results of the Interlaboratory Precision Study, by level, are given in the following tables. The values obtained for CD45+ and CD34+ absolute counts and CD34 percentages showed little variation as evidenced by the means, standard deviations, and coefficients of variation for the replicate measurements.

See Appendix for graph.

Intralaboratory Precision

Specimens covering the dynamic range of the assay were obtained by using fresh peripheral whole blood, fresh bone marrow, and fresh and thawed: apheresis and cord blood from four external evaluation sites. Ten samples from each were prepared according to the Stem-Kit reagents package insert and analyzed on COULTER EPICS XL/XL-MCL flow cytometers using manual gating and analysis to obtain CD45+, CD34+, and CD34% values.

The values for CD45+, CD34+, and CD34% given in the next tables demonstrate that the method of analysis is reproducible as evidenced by the means, standard deviations, and coefficients of variation for these measurements. The data tables shown in APPENDIX demonstrate Intralaboratory Precision for the assay.

See Appendix for graph.

LIMITATIONS

Refer to the stemCXP System Guide provided with the stemCXP System Software for limitations when using the automated method.

- Stained and lysed samples must be stored on ice and analyzed on the flow cytometer according to the information in section "Processed Sample stability" otherwise erroneous results can occur.
- Do not dilute, aliquot, or freeze the reagents. Use only as packaged.
- The reagents are for flow cytometry use only.
- Thaw frozen specimens properly to ensure optimal staining.
- All red blood cells may not lyse under the following conditions: presence of nucleated red blood cells, abnormal protein concentration, or hemoglobinopathies. This may cause falsely decreased CD34% results due to unlysed red blood cells being counted as leukocytes. CD34+ and CD45+ absolute counts should not be affected.
- Certain patients may present special problems due to altered or very low numbers of certain cellular populations.
- Results obtained with flow cytometry may be erroneous if the gates are improperly set.
- Stem-Kit Reagents is designed to use 7-AAD Viability Dye to determine HPC viable cells. If your laboratory is using an alternative viability determination method, follow general good laboratory practices and your laboratory's quality control procedures to ensure the quality of the reported results.

TRADEMARKS

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

ADDITIONAL INFORMATION

www.beckmancoulter.com

Symbols Key

Glossary of Symbols is available at beckman.com/techdocs (document number B60062)

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



APPENDIX

- EXAMPLES

The following series of eight histograms correspond to the analysis of an apheresis sample stained with Stem-Kit Reagents.

NOTE: These Figures feature histograms obtained on a COULTER EPICS XL/ XL-MCL flow cytometer equipped with System II Software (version 3.0).

	Histogram 1 • Histogram 1 displays viable events from region J (see Histogram 8). If not using 7-AAD, Histogram 1 displays all events minus all Stem-Count Fluorospheres (i.e. « -H » see Histogram 5). NOTE : • Position Region A to include all CD45+ events (leukocytes) while excluding CD45- events. Position Region E to include only lymphocytes (bright CD45, low Side Scatter). • Region A is intended to serve as an appropriate denominator for viable WBC in the calculation of the percentage of CD34+HPC. Once the analysis is done, the absolute count of viable leukocytes is given on the statistic printout related to Region A.
	Histogram 2 • Histogram 2 displays events from region A and J. If not using 7-AAD, Histogram 2 displays events from region A. • Adjust Region B to surround CD34+ HPC, including CD34^{dim} events. The upper limit along the y-axis may be set to include CD34+ events with low to intermediate Side Scatter for CD34+ HPC.
	Histogram 3 • Histogram 3 displays events from A and B and J. If not using 7-AAD, Histogram 2 displays events from region A and B. • Adjust Region C to include cells forming a cluster with characteristic CD34+ HPC (i.e. low Side Scatter and low to intermediate CD45 staining). Note the characteristic shape of Region C. It is the cluster of events that determines where Region C is centered. Brightly FITC stained events (platelet aggregates and mature monomyeloid cells) must be excluded in the setting of this region.
	Histogram 4 • Histogram 4 displays events from A and B and C and J. If not using 7-AAD, Histogram 2 displays events from region A and B and C. • Lymph / Blast Region D identifies a cluster of events meeting all the fluorescence and light scatter criteria of ISHAGE Guidelines for CD34+ HPC. Note the characteristic shape of Region D. If present, platelet aggregates that stain weakly must be excluded to the left of Region D. • Once the analysis is complete, the absolute count of viable CD34+HPC is given on the statistic printout related to the CD34+ Stem-Trol Control Cells or CD34+ HPC (Region D).

5: 	Histogram 5 displays all events. NOTE: This histogram is useful to visualize the lower limit of CD45 expression within the CD34+ events. • Set Quadstat Region I2 to enclose all CD45+CD34+ events. Set the left boundary of Quadstat Region I2 to include all ICD34+CD45^{dim} events. Verify on Histogram 1 that the lower limit of Region A includes all CD34+ events (use Histogram 5 as a guide). • Amorphous Region H is drawn to include all Stem-Count Fluorospheres. Region H must include the last channel of the FL1 Log (on the right) and FL2 Log scales (on the top). • Set Region K (named "listgate -") to exclude most of the double negative events from the acquisition.
6: E J 	Histogram 6 • Histogram 6 displays events from E and J regions. If not using 7-AAD, Histogram 2 displays events from region E. • Verify whether the Forward and Side Scatter gain parameters of the flow cytometer are optimally set for the processed sample. If necessary, adjust the Forward Scatter voltage / gain so that the smallest lymphocytes scatter in the middle-left part of the histogram. Adjust the Forward Scatter discriminator to ensure that even the smallest lymphocytes scatter above the discriminator.
7: H 	Histogram 7 • Histogram 7 displays events from region H. • Region G encloses the Stem-Count Fluorospheres singlet population. Check that the fluorosphere singlets accumulate homogeneously and constantly over time. Region G can be labeled as CAL to allow automatic calculation of absolute numbers of CD34+ HPC. Type the correct assayed concentration of the current lot of Stem-Count Fluorospheres (refer to the next section of this insert or to the instrument manual for further details).
8: A 	Histogram 8 • Histogram 8 displays events from region A. • Viable events are included within Region J. The 7-AAD Viability Dye positive events are excluded from Region J.

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Linearity


CD34+ Spiked Whole Blood Results

Linearity Results: Dynamic Range (Absolute Counts cells / μ L), Manual Gating and Analysis

Dilution	Test 1	Test 2	Control	Mean	Recovered Mean - Endogenous HPC	Expected Spiked CD34+ Cells
1	1829	1937	0	1883	1881	1875
1/2	882	905	0	893.5	892	938
1/4	461	456	0	458.5	457	469
1/8	237	235	0	236	234	234
1/16	122	115	0	118.5	117	117
1/32	58	66	0	62	60	59

Dilution	Test 1	Test 2	Control	Mean	Recovered Mean - Endogenous HPC	Expected Spiked CD34+ Cells
1/64	37	33	0	35	33	29
1/128	19	18	0	18.5	17	15
1/256	11	11	0	11	9	7
1/512	7	6	0	6.5	5	4
Endogenous	2	2	0	2		

Linearity, Dynamic Range, Manual Gating and Analysis

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
Linearity Results: Low End Sensitivity (Absolute Counts cells / μ L), Manual Gating and Analysis

Dilution	Test 1	Test 2	Control	Mean	Recovered Mean - Endogenous HPC	Expected HPC
1	172	140	0	156	154	154
1/2	87	83	0	85	83	77
1/4	42	38	0	40	38	39
1/8	23	21	0	22	20	19
1/16	10	10	0	10	8	10
1/32	9	9	0	9	7	5
1/64	6	5	0	5.5	4	2
1/128	3	4	0	3.5	2	1
Endogenous	2	2	0	2		

Linearity, Low End Sensitivity, Manual Gating and Analysis


Accuracy of Method

Accuracy of Method (Expected vs. Recovered Results)

Automated Gating and Analysis		Manual Gating and Analysis	
Expected Results	Recovered Results	Expected Results	Recovered Results
1875	1868.5	1875	1881
938	871	938	891.5
469	448.5	469	456.5
234	234	234	234
117	115	117	116.5

Automated Gating and Analysis		Manual Gating and Analysis	
Expected Results	Recovered Results	Expected Results	Recovered Results
59	58	59	60
29	31.5	29	33
15	16	15	16.5
7	6.6	7	9
4	4	4	4.5

		Automated Gating and Analysis	Manual Gating and Analysis
Slope	Estimate	0.985	0.993
	p-value*	0.216	0.481
Intercept	Estimate	-3.738	-0.009
	p-value*	0.641	0.999
*T-tests for Slope = 1 and Intercept = 0			

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Accuracy of Method - CD34+ Absolute Count (cells / μ L)

Lab No. / Sample Type (Method)	n	CD34+ (%)						
								CI
		Min	Max	Mean	±1 SD	CV	Lower Limit	Upper Limit
Lab 1 / Fresh Peripheral Blood (Automated)	13	0	191	66.77	62.45	93.53	22.43	111.11
Lab 1 / Fresh Peripheral Blood (Manual)	13	1	208	73.88	66.31	89.75	26.80	10.97
Lab 2 / Fresh Peripheral Blood (Automated)	23	0	313	40.09	67.82	169.18	6.07	74.10
Lab 2 / Fresh Peripheral Blood (Manual)	23	0	312	39.00	67.64	173.74	5.07	72.93
Lab 3 / Fresh Peripheral Blood (Automated)	23	2	215	65.24	65.54	100.46	32.37	98.11
Lab 3 / Fresh Peripheral Blood (Manual)	23	2	223	64.41	67.57	104.90	30.52	98.30
Lab 1 / Fresh Apheresis Products (Automated)	7	202	860	463.07	254.93	55.05	177.03	749.11
Lab 1 / Fresh Apheresis Products (Manual)	7	215	856	474.50	257.18	54.20	185.93	763.07
Lab 2 / Fresh Apheresis Products (Automated)	26	9	382	131.44	105.90	80.57	81.92	180.97
Lab 2 / Fresh Apheresis Products (Manual)	26	6	394	123.35	105.79	85.77	73.87	172.82
Lab 3 / Fresh Apheresis Products (Automated)	18	17	1211	249.67	311.90	124.93	68.96	430.37
Lab 3 / Fresh Apheresis Products (Manual)	18	15	1299	250.00	329.25	131.70	59.24	440.76
Lab 4 / Thawed Apheresis (Automated)	23	4	5831	860.02	1327.55	154.36	194.15	1525.89
Lab 4 / Thawed Apheresis (Manual)	23	6	6192	808.50	1352.88	167.33	129.03	1487.07
Lab 2 / Fresh Cord Blood (Automated)	8	13	89	37.25	26.39	70.84	10.74	63.76
Lab 2 / Fresh Cord Blood (Manual)	8	9	79	33.38	25.16	75.39	8.10	58.65
Lab 3 / Fresh Cord Blood (Automated)	21	12	113	40.33	24.75	61.37	27.24	53.42
Lab 3 / Fresh Cord Blood (Manual)	21	13	113	37.64	23.89	63.45	25.01	50.27
Lab 3 / Thawed Cord Blood (Automated)	21	3	618	98.74	174.62	176.85	6.41	191.07
Lab 3 / Thawed Cord Blood (Manual)	21	4	631	102.52	175.73	171.40	9.60	195.44
Lab 4 / Thawed Cord Blood (Automated)	10	6	48	17.85	14.48	81.13	5.55	30.15
Lab 4 / Thawed Cord Blood (Manual)	10	4	49	17.70	15.31	86.51	4.70	30.70
Lab 1 / Fresh Bone Marrow (Automated)	14	1	1740	192.29	451.05	234.57	27.24	53.42
Lab 1 / Fresh Bone Marrow (Manual)	14	1	2078	216.96	540.76	249.24	25.01	50.27
Lab 3 / Fresh Bone Marrow (Automated)	19	28	1115	274.45	271.99	99.11	10.74	63.76
Lab 3 / Fresh Bone Marrow (Manual)	19	19	1188	279.24	297.32	106.48	8.10	58.65
Lab 2 / Thawed Bone Marrow (Automated)	20	33	293	145.05	80.03	55.17	101.50	188.60
Lab 2 / Thawed Bone Marrow (Manual)	20	30	278	139.03	78.15	56.21	96.50	181.55

Accuracy of Method - Total Leukocytes - CD45+ (cells / μ L)


Lab No. / Sample Type (Method)	n	CD34+ (cells / μ L)						
								CI
		Min	Max	Mean	±1 SD	CV	Lower Limit	Upper Limit
Lab 1 / Fresh Peripheral Blood (Automated)	13	1835	59769	25149.12	18440.48	73.32	12055.93	38242.30
Lab 1 / Fresh Peripheral Blood (Manual)	13	1906	60968	26547.69	19947.94	75.14	12384.18	40711.21
Lab 2 / Fresh Peripheral Blood (Automated)	23	382	56177	27171.80	19640.81	72.28	17320.47	37023.14
Lab 2 / Fresh Peripheral Blood (Manual)	23	387	56382	27242.26	19722.48	72.40	17349.96	37134.56
Lab 3 / Fresh Peripheral Blood (Automated)	23	2236	32481	17605.65	7236.74	41.10	13975.89	21235.42
Lab 3 / Fresh Peripheral Blood (Manual)	23	2414	32917	17646.89	7397.23	41.92	13936.63	21357.16
Lab 1 / Fresh Apheresis Products (Automated)	7	22692	44772	28876.00	7721.27	26.74	20212.30	37539.70
Lab 1 / Fresh Apheresis Products (Manual)	7	22808	45083	29027.93	7781.27	26.81	20296.90	37758.95
Lab 2 / Fresh Apheresis Products (Automated)	26	9265	37410	17147.08	6411.77	37.39	14148.54	20145.61
Lab 2 / Fresh Apheresis Products (Manual)	26	9244	37582	17162.75	6457.88	37.63	14142.65	20182.85
Lab 3 / Fresh Apheresis Products (Automated)	18	11378	33452	20213.42	5444.91	26.94	17058.80	23368.03
Lab 3 / Fresh Apheresis Products (Manual)	18	11734	33070	20514.58	5378.22	26.22	17398.61	23630.56
Lab 4 / Thawed Apheresis (Automated)	23	3939	33245	16326.83	7960.22	48.76	12334.18	20319.47
Lab 4 / Thawed Apheresis (Manual)	23	3921	32479	15798.93	7804.07	49.40	11884.61	19713.26
Lab 2 / Fresh Cord Blood (Automated)	8	4550	17188	8958.38	4692.19	52.38	4442.92	13671.83
Lab 2 / Fresh Cord Blood (Manual)	8	4566	17200	8966.81	4709.09	52.52	4236.39	13697.24
Lab 3 / Fresh Cord Blood (Automated)	21	8764	18703	14032.79	2710.84	19.32	12599.39	15466.18
Lab 3 / Fresh Cord Blood (Manual)	21	8882	19046	14094.02	2751.92	19.53	12638.90	15549.15
Lab 3 / Thawed Cord Blood (Automated)	21	484	30866	7835.90	7082.77	90.39	4090.77	11581.04
Lab 3 / Thawed Cord Blood (Manual)	21	542	31789	8719.83	7192.23	82.48	4916.82	12522.85
Lab 4 / Thawed Cord Blood (Automated)	10	2134	4931	3113.25	1118.39	35.92	2163.65	4062.85

Lab No. / Sample Type (Method)	n	CD34+ (cells / μ L)						
								CI
		Min	Max	Mean	± 1 SD	CV	Lower Limit	Upper Limit
Lab 4 / Thawed Cord Blood (Manual)	10	1518	4322	2724.40	1061.63	38.97	1823.00	3625.80
Lab 1 / Fresh Bone Marrow (Automated)	14	938	37342	8968.36	8829.40	98.45	12599.39	15466.18
Lab 1 / Fresh Bone Marrow (Manual)	14	961	38666	9075.61	9151.80	100.84	12638.90	15549.15
Lab 3 / Fresh Bone Marrow (Automated)	19	5070	26396	13369.29	6192.49	46.32	4244.92	13671.83
Lab 3 / Fresh Bone Marrow (Manual)	19	4820	27284	10646.42	5912.00	55.53	4236.39	13697.24
Lab 2 / Thawed Bone Marrow (Automated)	20	2025	7977	3977.68	1358.94	34.16	3238.23	4717.12
Lab 2 / Thawed Bone Marrow (Manual)	20	1884	7742	3787.73	1322.52	34.92	3068.10	4507.35


Accuracy of Method (Figures)

Fresh Peripheral Whole Blood Specimens (n = 59)


Accuracy, Fresh Peripheral Whole Blood: Average CD34%

Accuracy, Fresh Peripheral Whole Blood: Average CD34 Absolute Counts


Accuracy, Fresh Peripheral Whole Blood: Average Total Leukocytes


Fresh Apheresis Specimens (n = 51)


Accuracy, Fresh Apheresis Specimens: Average CD34%

Accuracy, Fresh Apheresis Specimens: Average CD34 Absolute Counts


Accuracy, Fresh Apheresis Specimens: Average Total Leukocytes


Fresh Cord Blood Specimens (n = 29)


Accuracy, Fresh Cord Blood: Average CD34%

Accuracy, Fresh Cord Blood: Average CD34% Absolute Counts


Accuracy, Fresh Cord Blood: Average Total Leukocytes


Fresh Bone Marrow Specimens (n = 33)


Accuracy, Fresh Bone Marrow Specimens: Average CD34%

Accuracy, Fresh Bone Marrow Specimens: Average CD34 Absolute Counts


Accuracy, Fresh Bone Marrow Specimens: Average Total Leukocytes


Thawed Apheresis Specimens (n = 23)


Accuracy, Thawed Apheresis Specimens: Average CD34%

Accuracy, Thawed Apheresis Specimens: Average CD34 Absolute Counts


Accuracy, Thawed Apheresis Specimens: Average Total Leukocytes


Thawed Cord Blood Specimens (n = 31)


Accuracy, Thawed Cord Blood Specimens: Average CD34%

Accuracy, Thawed Cord Blood Specimens: Average CD34 Absolute Counts


Accuracy, Thawed Cord Blood Specimens: Average Total Leukocytes


Thawed BoneSpecimens (n = 20)


Accuracy, Thawed Apheresis Specimens: Average CD34 Absolute Counts

Accuracy, Thawed Bone Marrow Specimens: Average CD34 Absolute Counts

Accuracy, Thawed Bone Marrow Specimens: Average CD34 Total Leukocytes

Precision

Interlaboratory Precision

Interlaboratory Precision: Low Level

Instrument No. / Method	n	CD45+ (cells / μ L)			CD34+ (cells / μ L)			CD34+ (%)		
		Mean	± 1 SD	CV	Mean	± 1 SD	CV	Mean	± 1 SD	CV
1 / Manual	10	10637	774	7.28	12	1.25	10.70	0.12	0.01	10.90
2 / Manual	10	10318	424	4.11	11	1.29	11.80	0.11	0.01	10.51
3 / Manual	10	10187	260	2.55	11	1.43	13.49	0.11	0.01	12.90

Interlaboratory Precision: Mid-Range Level

Instrument No. / Method	n	CD45+ (cells / μ L)			CD34+ (cells / μ L)			CD34+ (%)		
		Mean	± 1 SD	CV	Mean	± 1 SD	CV	Mean	± 1 SD	CV
1 / Manual	10	13217	410	3.10	165	4.77	2.89	1.26	0.03	2.13
2 / Manual	10	12786	468	3.66	169	7.43	4.40	1.33	0.05	3.96
3 / Manual	10	12759	357	2.79	164	6.80	4.14	1.30	0.05	4.09

Interlaboratory Precision: High Level

Instrument No. / Method	n	CD45+ (cells / μ L)			CD34+ (cells / μ L)			CD34+ (%)		
		Mean	± 1 SD	CV	Mean	± 1 SD	CV	Mean	± 1 SD	CV
1 / Manual	10	443240	676	3.05	5080	12.42	4.88	1.17	0.05	4.26
2 / Manual	10	449540	647	2.88	5220	11.27	4.33	1.18	0.04	3.38
3 / Manual	10	476560	1373	5.76	5540	15.56	5.62	1.18	0.04	3.19

Interlaboratory Precision: High Level

Manual Method Lab No. / Specimen Type	n	CD45+ (cells / μ L)			CD34+ (cells / μ L)			CD34 (%)		
		Mean	± 1 SD	CV	Mean	± 1 SD	CV	Mean	± 1 SD	CV
Lab 1 / Fresh Peripheral Blood	10	12636	204	1.62	24	2.5	10.08	0.193	0.019	9.76
Lab 2 / Fresh Peripheral Blood	10	6080	675	11.11	5	1.0	21.60	0.074	0.015	20.20
Lab 2 / Fresh Peripheral Blood	10	19861	437	2.20	37	3.0	8.08	0.187	0.016	8.53
Lab 3 / Fresh Peripheral Blood	10	20519	500	2.44	114	4.5	3.94	0.558	0.021	3.77
Lab 3 / Fresh Peripheral Blood	10	10043	432	4.30	4	0.8	19.17	0.044	0.008	18.41
Lab 2 / Fresh Apheresis Products	10	36209	1045	2.89	361	13.8	3.82	0.997	0.035	3.55

Manual Method Lab No. / Specimen Type	n	CD45+ (cells / μ L)			CD34+ (cells / μ L)			CD34 (%)		
		Mean	± 1 SD	CV	Mean	± 1 SD	CV	Mean	± 1 SD	CV
Lab 2 / Fresh Apheresis Products	10	3723	100	2.69	19	1.5	7.69	0.516	0.035	6.74
Lab 3 / Fresh Apheresis Products	10	12081	299	2.47	239	7.3	3.06	1.978	0.037	1.88
Lab 3 / Fresh Apheresis Products	10	11537	747	6.48	5	0.9	17.90	0.046	0.008	17.95
Lab 4 / Thawed Apheresis Products	10	9730	496	5.10	2573	205.8	8.00	26.410	0.836	3.17
Lab 3 / Fresh Cord Blood	10	17378	807	4.64	71	7.5	10.56	0.408	0.036	8.90
Lab 4 / Thawed Cord Blood	10	1639	435	26.54	12	3.5	29.89	0.718	0.101	14.08
Lab 4 / Thawed Cord Blood	10	3896	733	18.80	37	6.9	18.55	0.954	0.097	10.12
Lab 3 / Fresh Bone Marrow	10	10488	507	4.83	197	13.7	6.94	1.878	0.052	2.75

REFERENCES

1. Sutherland, D.R., Anderson, L., Keeney, M., Nayar, R., Chin-Yee, I. The ISHAGE guidelines for CD34+ cell determination by flow cytometry. 1996, *J. Hematotherapy*, 5, 213-226.
2. Gratama, J.W., Keeney, M., and Sutherland, D.R. Enumeration of Hematopoietic Stem and Progenitor Cells. *Current Protocols in Cytometry*, 1999, 6.4.1-6.4.22.
3. Sutherland, D.R., Keating, A. The CD34 antigen: structure, biology and potential clinical applications. 1992, *J. Hematotherapy*, 1, 115-129.
4. Krause, D.S., Fackler, M.J., Civin, C.I., Stratford, M. CD34: structure, biology and clinical utility. 1996, *Blood*, 87, 1-13.
5. To, L.B., Haylock, D.N., Simmons, P.J., Juttner, C.A. The biology and clinical uses of blood stem cells. 1997, *Blood*, 89, 2233-2258.
6. Sutherland D.R., Marsh, J.C.W., Davidson, J., Baker, M.A., Keating, A., Mellors, A. Differential sensitivity of CD34 epitopes to cleavage by *Pasteurella haemolytica* glycoprotease: Implications for purification of CD34-positive progenitor cells. 1992, *Exp. Hematol.*, 20, 590-599.
7. Greaves, M.F., Titley, I., Colman, S.M., Bühring, H.-J., Campos, L., Castoldi, G.L., Garrido, F., Gaudernack, G., Girard, J.-P., Inglès-Esteve, J., Invernizzi, R., Knapp, W., Lansdorp, P.M., Lanza, F., Merle-Béral, H., Parravicini, C., Razak, K., Ruiz-Cabello, F., Springer, T.A., van der Schoot, C.E., Sutherland, D.R. CD34 cluster workshop report. 1995, *Leucocyte Typing V*, eds. Schlossman, S.F., et al., Oxford University Press, 840-846.
8. Greaves, M.F., Brown, J., Molgaard, H.V., Spurr, N.K., Robertson, D., Delia, D., Sutherland, D.R. Molecular features of CD34: a hemopoietic progenitor cell-associated molecule. 1992, *Leukemia*, 6, 31-36.
9. Newman, W., Targan, S.R., and Fast, L.D. Immunobiological and immunochemical aspects of the T-200 family of glycoproteins. 1984, *Mol Imm*, 21, 1113-1121.
10. Fabre J.W. and Williams A.F. Quantitative serological analysis of a rabbit anti-rat lymphocyte serum and preliminary biochemical characterisation of the major antigen recognised. 1977, *Transplantation*, 23, 349-359.
11. Coffman, R. L., and Weissman, I.L. A B cell-specific member of the T200 glycoprotein family. 1981, *Nature*, B220 ; 289, 681-683.
12. Dalchau, R, and Fabre, J. W. Identification with a monoclonal antibody of a predominantly B lymphocyte- specific determinant of the human leukocyte common antigen Evidence for structural and possible functional diversity of the human leukocyte common molecule. 1981, *J Exp Med*, 153, 753-765.
13. Omary, M.B., Trowbridge, I.S., and Battifora, H.A. Human homologue of murine T200 glycoprotein. 1980, *J Exp Med*, 152, 842-852.
14. Dalchau, R, Kirkley, J. and Fabre, J.W. Monoclonal antibody to a human leukocyte-specific membrane glycoprotein probably homologous to the leukocyte-common (L-C) antigen of the rat. 1981, *Eur J Immunol*, 10, 737-744.
15. Zelenin, A.V., Poletaev, A.I., Stepanova, N.G., Barsky, V.E., Kolesnikov, V.A., Nikitin, S.M., Zhuze, A.L., Gnatchev, N.V. 7-Amino-actinomycin D as a specific fluorophore for DNA content analysis by laser flow cytometry. 1984, *Cytometry*, 5, 348-354.
16. Cobbold, S., Hale, G., Waldmann, H. Non-lineage, LFA-1 family, and leukocyte common antigens: new and previously defined clusters. 1987, *Leucocyte Typing III, White Cell Differentiation Antigens*, eds. McMichael, A.J., et al, Oxford University Press, 788-803.
17. Schmid, I., Krall, W.J., Uittenbogaart, C.H., Braun, J., Giorgi, J.V. Dead cell discrimination with 7-aminoactinomycin D in combination with dual color immunofluorescence in single laser flow cytometry. 1992, *Cytometry*, 13, 204-208.
18. Gill, J.E., Jotz, M.M., Young, S.G., Modest, E.J., Sengupta, S.K. 7-Amino-actinomycin D as a cytochemical probe. I. Spectral properties. 1975, *J. Histochem. Cytochem.*, 23, 793-799.
19. Lecoq, H., Ledru, E., Prévost, M.-C., Gougeon, M.-L. Strategies for phenotyping apoptotic peripheral human lymphocytes comparing ISNT, annexin-V, and 7-AAD cytofluorometric staining methods. 1997, *J. Immunol. Methods*, 209, 111-123.
20. Civin, C.I., Trischmann, T.M., Fackler, M.J., Bernstein, I.D., Bühring, H.J., Campos, L., Greaves, M.F., Kamoun, M., Katz, D.R., Landsdorp, P.M., Look, A.T., Seed, B., Sutherland, D.R., Tindle, R.W., Uchanska-Ziegler, B. "Report on the CD34 cluster workshop" in *Leucocyte Typing IV, White Cell Differentiation Antigens*, Knapp, W., et al., Eds., Oxford Univ. Press, p. 818-825.
21. Keeney, M., Chin-Yee, I., Weir, K., Popma, J., Nayar, R., Sutherland, D.R. Single Platform Flow Cytometric Absolute CD34+ Cell Counts Based on the ISHAGE Guidelines. 1998, *Cytometry (Communications in Clinical Cytometry)*, 34, 61-70.

Beckman Coulter do Brasil Com. e Imp. de Prod. de Lab. Ltda
Alameda Rio Negro, 500, 15º andar, Torre B Alphaville Industrial
CEP 06.454-00 Barueri, São Paulo, Brasil
CNPJ: 42.160.812/0001-44 Telefone: 0800-771-8818

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IMMUNOTECH SAS A Beckman Coulter Company, 130, Avenue de Lattre de Tassigny,

BP 177, 13276 Marseille Cedex 9, France, +(33) 4 91 17 27 27

www.beckman.com