

# Specification Sheet




## Product Description: Multiwell Plates

Cat #: 25-101, 25-106

Dimensions=mm

## Component Materials:

- 100% USP VI Virgin Crystal Class Polystyrene.
- Tray: Medical style recyclable PET
- Gamma resistant.



## Features:

- Available in individual recyclable PET
- Temperature range: -20°C to +50°C
- Autoclavability: No
- Non-pyrogenic
- Shelf life: 3 years after month of production

## STERILITY:

- Sterilized by gamma irradiation SAL10<sup>-6</sup>

## PRODUCT ISO 9000 CERTIFICATION:

- Manufactured in a class 100,000 room environment.
- Manufactured under Management System as per ISO9001:2008 and ISO13485

Toll Free: 800.789.5550 Fax: 888.789.0444 Web: [www.geneseesci.com](http://www.geneseesci.com) Email: [support@geneseesci.com](mailto:support@geneseesci.com)