

## Certificate of Analysis (Typical)

<b>Item Number</b>	<b>Product Description</b>	<b>Lot Number</b>
25-529, 25-529C	Albumin, Bovine Fraction V, 100 G	Typical
<b>Formula</b>	<b>Molecular Weight</b>	<b>CAS Number</b>
		9048-46-8

QC TEST/RELEASE DATE: MM/DD/YYYY

SUGGESTED RETEST DATE: MM/DD/YYYY

S.No	Test	Unit	Specifications	Test Value
1	Appearance		White to yellow with tan to green cast powder	conforms
2	Solubility		Clear to slightly hazy faint yellow to yellow green solution	conforms
3	Protein	%	≥ 98.0	conforms
4	Purity	%	≥ 98	conforms
5	Moisture	%	≤ 5.0	conforms
6	pH		6.5 - 7.5	conforms
7	Ash	%	≤ 2	conforms
8	Endotoxin	EU/mg	≤ 3	conforms
9	APC Average	cfu/g	≤ 1000	conforms
10	IgG Free		Conforms	conforms
11	Bioburden		Low Bioburden	conforms
12	BSE/TSE Free		BSE/TSE Free	conforms
13	Country of Origin		Report	conforms

This is to certify that samples of the lot number above were tested and found to comply with the stated specifications. Some or all test data on this COA have been transcribed from supplier's/vendor's/third parties' test report or COAs. RPI may or may not perform all testing on the finished product unless specified by customers. Genesee expressly disclaims all warranties, expressed or implied, including the implied warranties of merchantability and fitness for a particular use or purpose. Certain products (USP/FCC/NF/EP/BP/JP grades) are sold for use in Research and further manufacturing only. **NOT FOR DIRECT USE IN FOOD, DRUG, HUMANS, ANIMALS OR MEDICAL DEVICES.** Persons using this information should make their own determination regarding its suitability for their particular use.

**This is an electronically generated document and does not require signatures**

**Certified By**