


Zymo-Seq RiboFree® Universal cDNA Kit

One for All. Universal rRNA Depletion for Any Organism.


- **Universal Depletion:** Novel probe-free technology depletes rRNA from any organism.
- **Simplest Kit:** Prepare NGS-grade RiboFree cDNA from total RNA in as little as 1.75 hours.
- **Automation Friendly:** Streamlined protocol for increased scalability.

The Only Universal rRNA Depletion Use One cDNA Kit For Any Organism



RiboFree Universal Depletion enzymatically removes rRNA efficiently across different organisms validated by next-generation sequencing. 100 ng of total RNA from each species were used as input to prepare rRNA depleted cDNA and continued with library preparation, respectively. The standard 1-hour incubation for depletion reaction was implemented. Classification of the STAR-aligned reads was based on Ensembl annotations and RepeatMasker rRNA tracks from UCSC genome browser when applicable.

Simplest NGS-Grade cDNA Prep



cDNA ready for downstream analyses in as little as **1.75 hours!**

| Product | Cat. No. | Size |
|---------------------------------------|----------|----------|
| Zymo-Seq RiboFree® Universal cDNA Kit | R3001 | 12 preps |