



ACCESSORIES

Empty Acrylic Columns

PRODUCT	Empty Acrylic Columns
CAT.	AC8-3
Contains	3 Empty Acrylic Columns; 3 Column bodies, 6 End plugs and 6 Stop plugs.
Column Material	Acrylic resin
Frit Material	Two layers of mesh (Coarse and fine) ⁽¹⁾
End Plugs	Polypropylene (10 – 32 UNF female thread)
Max. Pressure (bar/psi)	3 / 42
Chemical Stability	Stable in all commonly used reagents
Capacity (ml of packed resin)	~8ml

⁽¹⁾This column can be used for any type of chromatography media with a particle size larger than 20 µm.
Larger amounts available on request.