

ACCESSORIES

Empty FPLC Columns

PRODUCT	Empty FPLC 8 ml	Empty FPLC 30 ml	Empty FPLC 45 ml	Empty FPLC 80 ml
Cat. No.	FPLC8-3	FPLC30-2	FPLC45-2	FPLC80-1
Contains	3 Columns, 9 frits & 6 caps	2 Columns, 6 frits & 4 caps	2 Columns, 6 frits & 4 caps	1 Column, 3 frits & 3 caps
Capacity (ml Of Packed Resin) (1)	8	30	45	80
Dimensions Diameter X Height(Mm)	~12 x 70	~21 x 87	~21 x 137	~26 x 144
Column Material	Polypropylene			
Frit Material	Polyethylene			
Frit Pore Material	10 µm			
Max Pressure (Psi/Bar/Mpa)	200 / 14 / 1.38			
Chemical Stability	Stable in all commonly used reagents			
Connections (2)	Standard connection usually supplied with the common FPLC instruments (Such as AKTATM)			
Accessories	FPLC LUER Lock-male/ CAT. No. FPLCM-1 FPLC LUER Lock-female/ CAT. No. FPLCF-1			

(1)Packed volume depends on the type of resin – only used as a guide

(2)For fitting the column to the FPLC equipment. If the user does not have standard connections. ABT supplies them CAT. No. FPLCM-1/FPLCF-1

Larger amounts available on request.

ABT TS Rev. 2022/A