

ACCESORIES

Empty Columns

PRODUCT	Plastic Small Columns (20 Units)	Plastic Columns (50 Units)	Plastic XL Columns (50 Units)	Plastic XXL Columns (50 Units)
Cat. No.	CS-20	C-50	CXL-50	CXXL-50
Frite Pore Size	20 µm			
Column Material	Polypropylene			
Frit Material	Polyethylene			
Chemical Stability	Stable in all commonly used reagents			
Total Capacity (Ml)	1	12	35	60
Inner Diameter (Mm)	5.65	15.55	23.50	26.50
Outer Diameter (Mm)	7.65	17.85	26.10	29.60
Column Height (Mm)	65.55	90.15	122.90	155.0

⁽¹⁾Plastic columns include one frit and two caps (top and end). ABT supplies additional frits for all sizes of Plastic Columns.
For single-use application.

ABT TS Rev. 2022/A