

# pipetman®

THE PIPETTING STANDARD

**RECOGNIZED RELIABILITY AND PROVEN DURABILITY COMBINED  
WITH EASY PIPETTING AND HANDLING**



## Legendary Performance

Accurate and precise, PIPETMAN® air-displacement pipettes are recognized as one of the world standard. They come in a large range of volumes, including single and multichannel models, allowing you to handle various sample types with confidence.

## Easy Pipetting

PIPETMAN feels like PIPETMAN® Classic but with lower forces for enhanced pipetting habits with a standard pipette.

## Built-to-Last and Convenient

Made of PVDF and stainless steel, PIPETMAN offers premium durability. Built-to-last, it is the economical solution for years of pipetting.

## Simple to Maintain

Routine cleaning and maintenance are all that is required to keep PIPETMAN in top condition for years. The PIPETMAN design allows for easy access to user-serviceable parts, along with a fully autoclavable tip holder and ejector.



## Choose Your Pipette Stand

Different pipette stands and holders allow you to store pipettes vertically, protecting them from contamination.

CARROUSEL™



TRIO™



Universal  
Multichannel  
Stand



SINGLE™



### Pipette Stands

Description	Part Number
CARROUSEL™ Pipette Stand (7 pipettes)	F161401
TRIO™ Pipette Stand (3 pipettes)	F161405
Universal Multichannel Stand	F161417
SINGLE™ Pipette Holder	F161406

## Choose Your Pipette Kits

PIPETMAN single channel models are available in kits, offering the essential pipetting tools to start an experiment right away.



### PIPETMAN® Kits

Model	Part Number	Models
PIPETMAN® 4-Pipette Kit	F167360	P2 P20 P200 P1000
PIPETMAN® Starter Kit	F167900	P20 P200 P1000
PIPETMAN® Micro-volume Kit	F167800	P2 P10 P100
PIPETMAN® Large-volume Kit	F167920	P5000 P10mL

### PIPETMAN® Single Channel

Model	PIPETMAN® DIAMOND Tips	Part Number	Volume	
P2	D10 DL10	DF10ST DFL10ST	F144054M	0.2-2 µL
P10	D10 DL10	DF10ST DFL10ST	F144055M	1-10 µL
P20	D200	DF30ST	F144056M	2-20 µL
P100	D200	DF100ST	F144057M	10-100 µL
P200	D200 D300	DF200ST DF300ST	F144058M	20-200 µL
P1000	D1000 D1200	DF1000ST DF1200ST	F144059M	100-1000 µL
P5000*	D5000		F144066	500-5000 µL
P10mL*	D10mL		F144067	1-10 mL


### PIPETMAN® Multichannel

P8x10 P12x10	D10 DL10	DF10ST DFL10ST	F144068 F144069	1-10 µL
P8x20 P12x20	DL10 D200	DFL10ST DF30ST	F144070 F144071	2-20 µL
P8x200 P12x200	D200 D300	DF200ST DF300ST	F144072 F144073	20-200 µL
P8x300 P12x300	D200 D300	DF200ST DF300ST	F144074 F144075	30-300 µL

Models equipped with a tip ejector are only sold with a stainless steel tip ejector. Plastic tip ejector is available as an accessory. Refer to *PIPETMAN® User's Guide LT801122*.

\*P5000 and P10mL do not have tip ejectors.

Gilson maximum permissible errors are guaranteed only when PIPETMAN® pipettes are used with the recommended PIPETMAN DIAMOND Tips.



PIPETMAN DIAMOND Tips are designed to fit your PIPETMAN perfectly and guarantee maximum accuracy and precision.

For ordering information, visit [www.gilson.com](http://www.gilson.com)

[gilson.com/contactus](http://gilson.com/contactus)

LT800591/F | ©2023 Gilson S.A.S.  
All rights reserved. All trademarks  
are the property of Gilson, Inc.  
and its subsidiaries.

