

FLYSTUFF® FLY VIAL RELOADS

- Huge Time Saver
- Environmentally Responsible
- Excellent Value

Introducing a truly revolutionary breakthrough in Drosophila research products: the new fly vial reload system. The system saves valuable time, is environmentally responsible and economically priced. This patented product has forever changed how plastic vials are purchased and used by Drosophila research labs.

The reload system tray is precision plastic injection molded and features a self-contained divider/partition system. 100 plastic vials can be loaded in seconds into these unique trays. The vial stacking plates are made of RPET which have a recycling resin code of "1" and can therefore be readily incorporated into any recycling program.



Reload Pack Vials (Available in narrow only) Narrow (25x95mm) vials, packaged in 5 reload trays of 100 vials/unit


Cat #	Type	Material	O.D. x H	Master Case Qty	Units/Pallet
32-109RL	Narrow	Polystyrene	25 x 95mm	4 Units	56
32-113RL	Narrow	Polypropylene	25 x 95mm	4 Units	56

Fly Vial Reload Trays

For narrow plastic vials

Cat #	Color	List Price
59-207BLK	● Black	140.00
59-207BLU	● Blue	140.00
59-207GRE	● Green	140.00
59-207RED	● Red	140.00
59-207WH	○ White	140.00

For narrow vials (25x95mm). Packaged 12 trays/unit. Not autoclavable. If machine washed, avoid dry cycle. Plastic is constructed of a rugged, durable polymer blend. Fully submersible and bleachable. Dimensions: 10.375" x 10.375" x 2" (26.35cm x 26.35cm x 5cm).




Blue Tray Pictured



USA Phone: 800.789.5550
Int'l Phone: +1 858.536.8044
Fax: +1 858.536.8087
Order Online: www.geneseesci.com
Int'l Online: www.flystuff.com/distributors


①

cut open box and orient
 ouvrir le carton et le mettre debout
 開梱
 abre el cartón y orientalo
 den Karton oeffnen und senkrecht stellen
 把纸盒的方位摆正




③

load tray
 positionner le rack
 トレイを被せる
 coloca la bandeja
 das Tray auf die Reloads anpassen
 装上盘子



②

pull stack of vials
 sortir les recharges du carton
 積重ねバイアルを引き出す
 jala hacia afuera los tubos
 die Reloads heraus ziehen
 拉出整叠试管



④

flip tray
 retourner le rack et sa recharge
 トレイをひっくり返して完成
 voltea la bandeja
 das Tray mit seinem Reload umdrehen
 翻转过来

