

Certificate of Analysis

Product Name: Flystuff Cotton Balls, Degreased/scoured (fat-free)

Cat #: 51-101, 51-102B, 51-102W

Lot #: 202601

Quality Assurance:

Inspection Standard: YY/T 0330-2015 Medical Absorbent Cotton

Item	Technical Requirement	Result	Status
Appearance	Products should be spread out white in color, soft and elastic fiber, pollution-free, no damage, no impurities, odorless, tasteless.	Meets Specifications	PASS
Specifications	Error control at $\pm 10\%$	Meets Specifications	PASS
Alkalinity	No solution should be shown in pink	Meets Specifications	PASS
Suction Time	Less than 10 second	3.15s	PASS
Absorbency	Not less than 23.00g/g	23.09g	PASS
Loss on Drying	Humidity shall not exceed 8%	5.98%	PASS
Surface Active Substances	Liquid surfactant foam shall not cover the entire liquid surface.	Meets Specifications	PASS
Fluorescence	Not show strong blue fluorescence	Meets Specifications	PASS