

Cell Culture Dish

(Cat #: 25-200, 25-260, 25-202, 25-203)

1. Description

Genesee Scientific provides a wide range of Cell Culture Dishes with highest quality. All Cell Culture Dishes are produced with optically clear high quality polymers for microscopy.

2. Features

- ❖ Unique lid inner design for effective gas exchange
- ❖ External grip for better handling (Cat. No. 25-200, 25-260, 25-202)
- ❖ Excellent stackability
- ❖ USP guideline compliant (USP class VI tested)

3. Functional Test

Sterility	:	Passed
Cytotoxicity	:	Passed
Non-Pyrogenic	:	Passed
RNase/DNase free	:	Passed
DNA free	:	Passed

4. Product Range

Cat. No.	Dish Style d x h (mm)	Internal Dimension d x h (mm)	Growth Area (cm ²)	Working Volume (ml)	External Grip	Surface Treatment	Sterile	Packaging
25-200	35.00 x 10.00	35.00 x 9.60	9.40	3.00	+	+	+	10/500
25-260	60.00 x 15.00	52.80 x 12.80	21.50	5.00	+	+	+	10/600
25-202	100.00 x 20.00	87.48 x 16.40	57.50	12.50	+	+	+	10/300
25-203	150.00 x 20.00	138.57 x 15.40	148.00	35.00	-	+	+	5/100

5. Certification

Genesee Scientific certifies that the products identified above are inspected to be in compliance with product quality specification and requirements as documented in the manufacturer's ISO 9001:2015 Quality Management System.



6. Specifications & Technical Drawing

(1) Cell Culture Dish Ø35 (# 25-200)

- Specifications

Specifications				Material (Resin)	Packaging
Lid		Dish			
Diameter:	40.30 mm	Inner Diameter:	35.00 mm	Crystal-grade polystyrene (PS)	▪ Pieces/Sleeve: 10 ▪ Pieces/Case: 500
Height:	6.00 mm	Inner Height:	9.60 mm		
		Volume/Dish:	3.00 ml		
		Growth Area:	9.40 cm ²		

- Technical Drawing



(2) Cell Culture Dish Ø60 (# 25-260)

- Specifications

Dimensions		Material (Resin)	Packaging
Lid	Dish		
Diameter: 58.00 mm Height: 8.50 mm	Inner Diameter: 52.80 mm Inner Height: 12.80 mm Volume/Dish: 5.00 ml Growth Area: 21.50 cm ²	Crystal-grade polystyrene (PS)	- Pieces/Sleeve: 10 - Pieces/Case: 600

- Technical Drawing




(5) Cell Culture Dish Ø100 External Grip (# 25-202)

- Specifications

Dimensions		Material (Resin)	Packaging
Lid	Dish		
Diameter: 95.67 mm Height: 11.30 mm	Inner Diameter: 87.48 mm Inner Height: 16.40 mm Volume/Dish: 12.50 ml Growth Area: 57.50 cm ²	Crystal-grade polystyrene (PS)	- Pieces/Sleeve: 10 - Pieces/Case: 300

- Technical Drawing



(6) Cell Culture Dish Ø150 (# 25-203)

- Specifications

Dimensions		Material (Resin)	Packaging
Lid	Dish		
Diameter: 147.66 mm Height: 10.00 mm	Inner Diameter: 138.57 mm Inner Height: 15.40 mm Volume/Dish: 35.00 ml Growth Area: 148.00 cm ²	Crystal-grade polystyrene (PS)	- Pieces/Sleeve: 5 - Pieces/Case: 100

- Technical Drawing

