

Product Specification Sheet

Product Name: 6 Well Tissue Culture Plates

Cat #: 25-100 (Non-Treated) / 25-105 (Surface Treated)

Description

6 well, Flat Bottom, Surface Treated or Non-Treated, Sterile


Purpose

For use in applications where cell attachment is desired or not desired.

Materials

- Plate: GPPS (General Polystyrene)
- Lid: GPPS (General Polystyrene)

Dimensions (Unit: mm)



Features

- Flat bottom round well
- Lid with notched corners for good fit and easy orientation
- Wells are labeled with alphanumeric code for easy identification
- Assembled plate and lid height: 22.2 mm
- Max. volume (single well): 17.0 mL
- Working volume (single well): 1.9 ~ 2.9 mL
- Growth area (single well): 9.6 cm²
- Autoclavability: No
- DNase/RNase free and non-pyrogenic
- Sterilized by irradiation SAL10⁻⁶ (ISO 11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 room environment
- Manufactured under ISO 13485 and ISO 9001 quality management system

Product Specification Sheet

Product Name: 12 Well Tissue Culture Plates

Cat #: 25-101 (Non-Treated) / 25-106 (Surface Treated)

Description

12 well, Flat Bottom, Surface Treated or Non-Treated, Sterile


Purpose

For use in applications where cell attachment is desired or not desired.

Materials

- Plate: GPPS (General Polystyrene)
- Lid: GPPS (General Polystyrene)

Dimensions (Unit: mm)



Features

- Flat bottom round well
- Lid with notched corners for good fit and easy orientation
- Wells are labeled with alphanumeric code for easy identification
- Assembled plate and lid height: 22.7 mm
- Max. volume (single well): 6.8 mL
- Working volume (single well): 0.76 ~ 1.14 mL
- Growth area (single well): 3.85 cm²
- Autoclavability: No
- DNase/RNase free and non-pyrogenic
- Sterilized by irradiation SAL10⁻⁶ (ISO 11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 room environment
- Manufactured under ISO 13485 and ISO 9001 quality management system

Product Specification Sheet

Product Name: 48 Well Tissue Culture Plates

Cat #: 25-103 (Non-Treated) / 25-108 (Surface Treated)

Description

48 well, Flat Bottom, Surface Treated or Non-Treated, Sterile


Purpose

For use in applications where cell attachment is desired or not desired.

Materials

- Plate: GPPS (General Polystyrene)
- Lid: GPPS (General Polystyrene)

Dimensions (Unit: mm)



Features

- Flat bottom round well
- Lid with notched corners for good fit and easy orientation
- Wells are labeled with alphanumeric code for easy identification
- Assembled plate and lid height: 22.7 mm
- Max. volume (single well): 1.55 mL
- Working volume (single well): 0.19 ~ 0.29 mL
- Growth area (single well): 0.83 cm²
- Autoclavability: No
- DNase/RNase free and non-pyrogenic
- Sterilized by irradiation SAL10⁻⁶ (ISO 11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 room environment
- Manufactured under ISO 13485 and ISO 9001 quality management system

