


Set Up Your Pipettes Remotely

The Eppendorf Pipette Manager puts connected electronic pipettes to use for faster operation, increased reproducibility and digital documentation of every step.



Who doesn't enjoy greater freedom and convenience when it comes to pipetting? Be ahead of the curve! Switch to connected electronic pipettes and boost your individual pipetting skills while bringing teamwork up to a new level.



Connected electronic pipettes

Convert your Eppendorf Xplorer®, Xplorer plus or Xplorer plus Move It® pipettes into **connected electronic pipettes**.

Pipette Manager

Establish communication with connected electronic pipettes and tablets via WiFi technology.

Tablets

Connect tablets (Android & iOS) to work in parallel with other lab users.

Try our demo tool

What Connectivity Can Do for You

Be faster

Quicker input of volumes via the Pipette Manager. All settings will be transferred to connected electronic pipettes for an immediate start of pipetting.

Be more accurate

Get guidance by the software-embedded "Liquid Manager" when working with challenging liquid classes (e.g. viscous, volatile, or foaming) to increase the reproducibility of your pipetting steps.

Choose from predefined Eppendorf liquid types, e.g. acetone, or create your own.

Be confident

Automatically document each pipetting step and easily export your recordings as a PDF file via USB. The documentation contains details such as timestamp, serial number of the pipette, pipetting mode, volume as well as the selected liquid type.

Work more closely

Manage your pipette fleet more effectively and enhance collaboration in teams. Integrated features like "Pipette Set Creation" or "Liquid Type Sharing" support smooth and efficient operations.

Evolve Your Electronic Pipettes

Stay flexible with your electronic pipettes!

Easily convert your Eppendorf Xplorer/Xplorer plus or Xplorer plus Move It electronic pipettes into a connected device with the WiFi module.

Plug in, get connected and start working in less than a minute!



Ordering information

Description	Order no.
Pipette Manager , an external touch server enabling communication with connected electronic pipettes	1004 000 001
Eppendorf Xplorer® connect , WiFi module incl. battery for Eppendorf Xplorer electronic pipette family	4861 000 970

Note: The Pipette Manager is not available worldwide. Please contact your Eppendorf Sales Representative for more information.

Your local distributor: www.eppendorf.com/contact
 Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com/pipette-manager

Eppendorf®, the Eppendorf Brand Design, Eppendorf Xplorer® and Move It® are registered trademarks of Eppendorf SE, Germany. All rights reserved including graphics and images. Copyright © 2023 by Eppendorf SE. Order no. AA01 028 620/EN2/WEB/PDF