

Macro Volume Tip Compatibility Chart

5 mL Gilson / Universal Fit

Item #	Description	Color	Packaging	Eppendorf	Socorex	Gilson	Rainin Classic	Sartorius (Biohit)	Finnpipette Brand / Brantech
27-345B	Bulk	Clear	250/unit, 10 units/case						
27-345BF	Bulk, Filtered	Clear	100/unit, 10 units/case						
27-345R	Racked	Clear	50 tips/unit, 10 units/case			✓	✓		
27-345	Racked, Sterile	Clear	50 tips/unit, 10 units/case						
27-345F	Racked, Sterile, Filtered	Clear	50 tips/unit, 10 units/case						

5 mL Eppendorf Fit

Item #	Description	Color	Packaging	Eppendorf	Socorex	Gilson	Rainin Classic	Sartorius (Biohit)	Finnpipette Brand / Brantech
27-526	Bulk	Clear	250/unit, 10 units/case						
27-526BF	Bulk, Filtered	Clear	100/unit, 10 units/case						
27-526R	Racked	Clear	50 tips/unit, 10 units/case	✓	✓			✓*	
27-526RS	Racked, Sterile	Clear	50 tips/unit, 10 units/case						
27-526F	Racked, Sterile, Filtered	Clear	50 tips/unit, 10 units/case						

5 mL Thermo Finn / Brand Fit

Item #	Description	Color	Packaging	Eppendorf	Socorex	Gilson	Rainin Classic	Sartorius (Biohit)	Finnpipette Brand / Brantech
27-122	Bulk	Clear	5 bags of 100/unit, 2 units/cs						
27-122R	Racked	Clear	50 tips/unit, 10 units/case					✓	✓
27-122RS	Racked, Sterile	Clear	50 tips/unit, 10 units/case						
27-122F	Racked, Sterile, Filtered	Clear	50 tips/unit, 10 units/case						

10 mL Thermo Finn / Socorex Fit

Item #	Description	Color	Packaging	Eppendorf	Socorex	Gilson	Rainin Classic	Sartorius (Biohit)	Finnpipette Brand / Brantech
27-124	Bulk	Blue	200/unit, 5 units/case						
27-124R	Racked	Blue	25 tips/rack, 10 racks/case		✓				✓
27-124RS	Racked, Sterile	Blue	25 tips/rack, 10 racks/case						

10 mL Gilson / Eppendorf Fit

Item #	Description	Color	Packaging	Eppendorf	Socorex	Gilson	Rainin Classic	Sartorius (Biohit)	Finnpipette Brand / Brantech
27-346	Bulk	Clear	200/unit, 5 units/case						
27-346BF	Bulk, Filtered	Clear	50/unit, 10 units/case						
27-111R	Racked	Clear	25 tips/rack, 10 racks/case	✓	✓	✓	✓*		✓
27-111RS	Racked, Sterile	Clear	25 tips/rack, 10 racks/case						
27-112	Racked, Sterile, Filtered	Clear	25 tips/rack, 10 racks/case						

*Does not fit Proline pipettors, but does fit all other Biohit pipettors, including Proline Plus.

Tips shown in actual size.

