



Never Been Easier

Simply migrate your PCR programs with one step on Eppendorf Mastercycler® X50

Do you have unsatisfactory results after transferring protocols from other thermal cyclers and no time for optimizing your PCR program? The Program Migration Feature is an easy and convenient function to simulate a slower thermal cycler.

Just set your desired runtime of the PCR on the Mastercycler X50:

- > Automatic adjustment of temperature ramps through runtime-based algorithm
- > Optimized PCR protocols are ready to use once transferred to Mastercycler X50



www.eppendorf.com/next-stage