

Kit SDS Cover Sheet

Document ID: B93029-75: Version AF
Revision Date (year/month/day) 2023/12/22
Last Revision Date (year/month/day) 2023/05/17

Product information

Product name	Solution Pack
Part number	B93029

Components

Description	Bag 1 Bags 2 - 7
--------------------	---------------------

Transport information

Transportation of this product is not regulated under ICAO, IATA DGR, IMDG, US DOT, European ADR and RID or Canadian TDG.



SAFETY DATA SHEET

Document ID: B93029-75 Version AF
 Revision Date (year/month/day) 2023/12/22
 Last Revision Date (year/month/day) 2023/05/17

Section 1 Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

Product name Bag 1
Part number Component of P/N B93029

1.2 Relevant identified uses of the substance or mixture and uses advised against

Product use For Professional/Industrial Use. See product literature for details.

1.3 Details of the supplier of the safety data sheet

Manufacturer

Manufactured for
 Beckman Coulter, Inc.
 250 S. Kraemer Blvd
 Brea, CA 92821, U.S.A.
 Tel: 800-854-3633

Supplier

CANADA
 Beckman Coulter Canada LP
 7075 Financial Drive
 Mississauga, ON L5N 6V8
 Canada
 1-800-463-7828

UNITED KINGDOM
 Beckman Coulter (UK) Ltd.
 Oakley Court
 Kingsmead Business Park, London
 Road
 High Wycombe
 United Kingdom HP11 1JU
 01494 441181

AUSTRALIA
 Beckman Coulter Australia Pty Ltd
 23-27 Chaplin Drive
 Lane Cove NSW 2066
 Australia
 ABN 81 002 011 672
 24 Hour emergency contact phone
 number:
 1800 060 881

SWITZERLAND
 Beckman Coulter Eurocenter SA
 22, rue Juste-Olivier, Case Postale
 1044,
 CH-1260 Nyon 1, Switzerland.
 Telephone: +41 (0)22 365 36 11
 Monday through Friday, 9:00 am to
 7:00pm)

Section 1 Identification of the substance/mixture and of the company/undertaking (Continued)

NEW ZEALAND
Beckman Coulter NZ
Unit J, 33 Walmsley Road, Otahuhu,
Auckland 1062, New Zealand
Hours available: 08:30 - 17:00

ICELAND / ÍSLAND
Beckman Coulter AB
Ekbacksvägen 28
168 69 Bromma
Sweden
Phone No.: +46 80564 85 900
Hours available: 08.00-16.30

MALTA
DX Distributor:
Cherubino Ltd
DELFI Building, Sliema Road, Gzira,
GZR 1637
Telephone: +356 21343270
Hours available: 08:30 – 17:00
SDSNT@beckman.com

e-mail address

1.4 Emergency telephone number

Telephone number (24H)

Chemtrec Emergency Tel No. U.S.A. 800-424-9300, International (001) 703-527-3887

Distributor and emergency phone no.

Refer to attached list, Document ID: [472050](#), for local distributor and emergency phone numbers.

UNITED STATES - Emergency Phone (24h): Chemtrec (800) 424-9300, International (001) 703-527-3887

CANADA - Poison Centre: 1-844-764-7669; Centre antipoison du Québec: 1-800-463-5060

UNITED KINGDOM - For UK and Scotland: Emergency Call 999

IRELAND - National Poisons Information Centre Phone No.: Members of Public: +353 (01) 809 2166 (8:00 am to 10:00 pm 7 days a week); Healthcare Professionals: +353 (01) 809 2566 (24 hour service)

AUSTRALIA - 24 Hour emergency contact phone number: 1800 060 881

NEW ZEALAND - 24 Hour emergency number: 0800 446 109

Section 2 Hazards identification

2.1 Classification of the substance or mixture

Product description Mixture
Reflective Mylar; Bag containing gases; Pungent

Classification according to EC 1272/2008 (CLP/GHS)
Oxidizing Gases, Category 1, H270

Classification according to US-OSHA (HCS 29 CFR 1910.1200) and UN GHS
Oxidizing Gases, Category 1

2.2 Label elements

According to EC 1272/2008 (CLP/GHS), US-OSHA and UN GHS

Hazardous ingredients

Oxygen

Pictogram



Signal word

DANGER

Hazard statements

H270 May cause or intensify fire; oxidiser.

Precautionary statements

Prevention

P220 Keep away from clothing and other combustible materials.

Response

P370+P376 In case of fire: Stop leak if safe to do so.

Storage

P403 Store in a well-ventilated place.

Disposal

None

2.3 Other hazards

Results of PBT and vPvB assessment

PBT: Not applicable.

vPvB: Not applicable.

See Section 11 Toxicological Information for more detailed health information.

Section 3 Composition and information on ingredients

3.2 Mixtures

Hazardous ingredients:		Hazard classification of pure ingredients		
Chemical name	% by wt.	EU 1272/2008 CLP/GHS	GHS	Note

SAFETY DATA SHEET

Document ID: B93029-75 Version AF
Revision Date (year/month/day) 2023/12/22
Last Revision Date (year/month/day) 2023/05/17

Section 3 Composition and information on ingredients (Continued)

Oxygen CAS # 7782-44-7 EINECS # 231-956-9 Index # 008-001-00-8	40 - 45	Ox. Gas 1, H270 Press. Gas [CG], H280	Ox. Gas 1, H270 Press. Gas [CG], H280	
---	---------	--	--	--

See section 8 for available Occupational exposure limits
See Section 15 for additional regulatory information
See Section 16 for description of hazard class and hazard statements

Section 4 First aid measures

4.1 Description of first aid measures

Inhalation	If product is inhaled, move exposed individual to fresh air. If individual is not breathing, begin artificial respiration by trained personnel and obtain medical attention immediately.
Eye contact	If product enters eyes, rinse eyes gently with water as a precaution.
Skin contact	In case of skin contact, rinse with water as a precaution.
Ingestion	If product is ingested, rinse mouth with water. If irritation or discomfort occurs, obtain medical attention immediately.

4.2 Most important symptoms and effects, both acute and delayed

No adverse symptoms or effects have been identified.

4.3 Indication of any immediate medical attention and special treatment needed

No specific medical attention or treatment required.

Section 5 Firefighting measures

5.1 Extinguishing media In case of fire use carbon dioxide (CO₂), dry chemical, water spray or foam. For large fires use extinguishing media suitable for surrounding fire.

5.2 Special hazards arising from the substance or mixture Special fire and explosion hazards

May cause or intensify fire; oxidiser.

Hazardous combustion products

No combustion products posing significant hazards are expected from this product.

5.3 Advice for firefighters

Protective equipment Self-contained breathing apparatus is recommended for firefighters in all chemical fire situations.

Additional information No further relevant information available.

Section 6 Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

Personal precautions Observe general safety guidelines for protection; avoid eye and skin contact.

6.2 Environmental precautions Contain spill to prevent migration or evaporation.

6.3 Methods and material for containment and cleaning up

Spill and leak procedures Collect in containers and seal securely.

6.4 Reference to other sections Refer sections 8 and 13.

Section 7 Handling and storage

7.1 Precautions for safe handling Use good laboratory procedures; avoid eye and skin contact.

7.2 Conditions for safe storage, including any incompatibilities

To maintain product quality, store according to the instructions in the product labeling.

7.3 Specific end uses No further relevant information available.

Section 8 Exposure controls and personal protection

8.1 Control parameters

Exposure limits

US OSHA None established

ACGIH None established

ACGIH Biological Exposure Indices (BEI)

None established

DFG MAK None established

Ireland None established

IOELVs

None established

NIOSH None established

China None established

Croatia None established

Japan None established

Sweden (AFS 2015:7 and amendments)

None established

Turkey None established

SAFETY DATA SHEET

Document ID: B93029-75 Version AF
Revision Date (year/month/day) 2023/12/22
Last Revision Date (year/month/day) 2023/05/17

Section 8 Exposure controls and personal protection (Continued)

8.2 Exposure controls

Engineering controls	No special engineering controls are required. Use with good general ventilation.
Eye protection	Safety glasses or chemical goggles should be worn to prevent eye contact. Refer U.S. OSHA 29 CFR 1910.133, European Standard EN166 or appropriate government standards.
Skin protection	Wear impervious gloves such as Nitrile or equivalent and protective clothing. Refer to U.S. OSHA 29 CFR 1910.138, European Standard EN 374, EN 14605:2005+A1:2009 or appropriate government standards.
Respiratory protection	Under normal conditions, the use of this product should not require respiratory protection. If overexposure should occur and ventilation is not adequate to maintain airborne concentrations at acceptable levels, the use of respiratory protection should be evaluated by a qualified professional.

Section 9 Physical and chemical properties

9.1 Information on basic physical and chemical properties

Physical state	Bag containing gases	Density and/or relative density	Not determined
Color	Reflective Mylar	Solubility	
Odor	Pungent	Water	Not applicable
pH	Not applicable	Organic	Not applicable
Melting Point	Not determined	Partition coefficient n-octanol/water (log value)	Not applicable
Boiling point or initial boiling point and boiling range	Not determined	Auto-ignition temp.	Not determined
Flash point	Not determined	Decomposition temperature	Not determined
Flammability	Not determined	Vapor pressure	Not determined
		Kinematic viscosity	Not determined
Lower and upper explosion limit	Not determined		
Relative vapor density	Not determined		
Particle characteristics	Not determined		

9.2 Other information

Information with regard to physical hazard classes

No further relevant information available.

Section 9 Physical and chemical properties (Continued)

Other safety characteristics

No further relevant information available.

Section 10 Stability and reactivity

- 10.1 Reactivity** No further relevant information available.
- 10.2 Chemical stability** The product is stable in accordance with recommended storage conditions.
- 10.3 Possibility of hazardous reactions**
May cause or intensify fire; oxidiser.
- 10.4 Conditions to avoid** Keep away from heat, sparks and open flame.
- 10.5 Incompatible materials** No further relevant information available.
- 10.6 Hazardous decomposition products**
When stored as labeled, no known hazardous decomposition products are formed during the shelf-life of this product.

Section 11 Toxicological information

11.1 Information on hazard classes

Toxicity data for hazardous ingredients

Not applicable

Primary routes of exposure Eye contact, ingestion, inhalation, and skin contact.

Acute toxicity Not classified based on available data.

Skin corrosion/irritation Not classified based on available data.

Serious eye damage/irritation Not classified based on available data.

Respiratory or skin sensitisation Not classified based on available data.

Germ cell mutagenicity Not classified based on available data.

Carcinogenicity No ingredients in this product are listed as carcinogens by ACGIH, IARC, NTP, OSHA or 1272/2008 EC regulation.

Reproductive toxicity Not classified based on available data.

Specific target organ toxicity (STOT) – single exposure

Not classified based on available data.

Specific target organ toxicity (STOT) – repeated exposure

Not classified based on available data.

SAFETY DATA SHEET

Document ID: B93029-75 Version AF
Revision Date (year/month/day) 2023/12/22
Last Revision Date (year/month/day) 2023/05/17

Section 11 Toxicological information (Continued)

Aspiration hazard Not classified based on available data.

11.2 Information on other hazards

Endocrine disrupting properties

This product does not have substance(s) of endocrine disrupting properties for health according to REACH Article 57(f).

Other information No further relevant information available.

Section 12 Ecological information

12.1 Toxicity

Fresh water species No information available.

Microtox/organisms No information available.

Water flea No information available.

Fresh water algae No information available.

12.2 Persistence and degradability Not determined for the product.

12.3 Bioaccumulative potential Not determined for the product.

12.4 Mobility in soil Not determined for the product.

12.5 Results of PBT and vPvB assessment

Not determined for the product. PBT: Not applicable, vPvB: Not applicable.

12.6 Endocrine disrupting properties

This product does not have substance(s) of endocrine disrupting properties for environment according to REACH Article 57(f).

12.7 Other adverse effects No further relevant information available.

Section 13 Disposal considerations

13.1 Waste treatment methods

Product waste disposal

Chemical residues and remains should be routinely handled as special waste. This must be disposed of in compliance with anti-pollution and other laws of the country concerned. To ensure compliance we recommend that you contact the relevant (local) authorities and/or an approved waste-disposal company for information.

Additional information

Suggested European waste catalogue 18 01 07 - chemicals other than those mentioned in 18 01 06. Dispose in accordance with national, state and local waste regulations.

Section 14 Transport information

Transportation of this product is not regulated under ICAO, IATA DGR, IMDG, US DOT, European ADR and RID or Canadian TDG.

- 14.1 **UN/ID number:** Not regulated for transportation
- 14.2 **UN proper shipping name:** Not regulated for transportation
- 14.3 **Transport hazard class(es):** Not regulated for transportation
- 14.4 **Packing group:** Not regulated for transportation
- 14.5 **Environmental hazards:** Not regulated for transportation
- 14.6 **Special precautions for user:** None
- 14.7 **Maritime transport in bulk according to IMO instruments:** Not applicable

Section 15 Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

US Federal and State Regulations

SARA 313 (Section 313, Title III reporting requirements)

No ingredients listed.

CERCLA (The Comprehensive Environmental Response, Compensation, and Liability Act) 40 CFR 302.4

No ingredients listed.

California Proposition 65

Chemical which is known to the State of California to cause cancer

No ingredients listed.

Chemical which is known to the State of California to cause development toxicity

No ingredients listed.

Chemical which is known to the State of California to cause male reproductive toxicity

No ingredients listed.

Chemical which is known to the State of California to cause female reproductive toxicity

No ingredients listed.

Massachusetts Right To Know (RTK) List

CAS # 7727-37-9	Nitrogen
CAS # 124-38-9	Carbon Dioxide
CAS # 7782-44-7	Oxygen

Section 15 Regulatory information (Continued)

New Jersey Dept. of Health Right To Know (RTK) List

CAS # 7727-37-9	Nitrogen
CAS # 124-38-9	Carbon Dioxide
CAS # 7782-44-7	Oxygen

Pennsylvania Right To Know (RTK) List

CAS # 7727-37-9	Nitrogen
CAS # 124-38-9	Carbon Dioxide
CAS # 7782-44-7	Oxygen

EU Regulations

This SDS complies with EC Regulations 1907/2006 (REACH) and amendments.

Regulation (EU) 2019/1148 on the marketing and use of explosives precursors - Substances Subject to Suspicious Transactions Reporting

No ingredients listed.

Regulation (EU) 2019/1148 on the marketing and use of explosives precursors - Restricted Explosives Precursors

No ingredients listed.

REACH 1907/2006 EC - Candidate List of Substances of Very High Concern (SVHC)

No ingredients listed.

REACH 1907/2006 EC - Annex XVII – Restrictions on Certain Dangerous Substances

No ingredients listed.

REACH 1907/2006 EC - Annex XIV - list of substances subject to authorisation

No ingredients listed.

Refer to Section 3

UK Regulations

UK REACH Regulation (as Amended) - List of substances subject to authorisation

Refer to Section 3

Canada

This product meets WHMIS criteria for hazardous materials

China

Catalog of Hazardous Chemicals - Hazardous Chemicals

CAS # 7727-37-9	Nitrogen
CAS # 124-38-9	Carbon Dioxide
CAS # 7782-44-7	Oxygen

Inventory - China - Inventory of Existing Chemical Substances (IECSC)

All ingredients are listed or exempted.

SAFETY DATA SHEET

Document ID: B93029-75 Version AF
Revision Date (year/month/day) 2023/12/22
Last Revision Date (year/month/day) 2023/05/17

Section 15 Regulatory information (Continued)

Turkey

Turkey-REACH - KKDIK Regulation - Annex 17 – Restrictions

No ingredients listed.

International

UN/FAO/Rotterdam Convention - Chemicals Subject to Prior Informed Consent (PIC)

No ingredients listed.

15.2 Chemical safety assessment

A Chemical Safety Assessment has not been carried out.

Some hazardous ingredients listed in Section 15 are below the cutoff limits of 0.1% for carcinogen, mutagen and reproductive toxin and 1% for other health hazards required for reporting in Section 3.

Section 16 Other information

Beckman Coulter safety rating	Flammability: 3 Health: 0 Reactivity with water: 0 Physical contact: 0	Code 0=None 1=Slight 2=Caution 3=Severe
Revision changes	Updated sections 1, 2, 3, 4, 8 and 15	
Document version and issue/revision date	Revision Date (year/month/day) 2023/12/22 Last Revision Date (year/month/day) 2023/05/17 Document ID: B93029-75 Version: AF	
Hazard Classification Procedure	This mixture was classified using the calculation method for human health and environmental hazards. Physical hazards were determined based on the specification.	
Description of hazard class and hazard statements from Section 3	Press. Gas [CG] - Gases under pressure, Compressed Gas Ox. Gas 1 - Oxidizing Gases, Category 1 H270 - May cause or intensify fire; oxidiser. H280 - Contains gas under pressure; may explode if heated.	

Section 16 Other information (Continued)

Abbreviations and acronyms

ACGIH - American Conference of Governmental Industrial Hygienists (ACGIH)
ADR and RID - European Agreement Concerning the International Carriage of Dangerous Goods by Road and Rail
CLP - Classification, Labeling and Packaging
DFGMAK - Republic Germany's maximum exposure limit
EC50 - Concentration of a substance in an environmental medium expected to produce a certain effect in 50% of test organisms
GHS - Globally Harmonized System of Classification and Labelling of Chemicals (GHS)
HCS - Hazard Communication Standard
IARC - International Agency for Research on Cancer
IATA DGR - International Air Transport Association Dangerous Goods Regulation
ICAO - International Civil Aviation Organization
IDLH - Immediately Dangerous to Life or Health
IMDG - International Maritime Dangerous Goods
IMO - International Maritime Organization
IOELVs - European Unions' Indicative Occupational Exposure Limit Values
LC50 - Concentration of a substance in water causing death (50% of the tested population) to aquatic life
LD50 - Lethal Dose 50%
NIOSH - National Institute for Occupational Safety and Health
NTP - National Toxicology Program
OSHA - Occupational Safety and Health Administration
PBT - Persistent Bioaccumulative and Toxic substances
PEL - Permissible Exposure Limit
SARA - Superfund Amendments and Reauthorization Act
STEL – Short Term Exposure Limit
STLV - Short Term Limit Value
STV - Short Term Value
TDG - Canadian Transportation of Dangerous Goods Regulations
TLV - Threshold Limit Value
TWA – Time Weighted Average
UN GHS - United Nations Globally Harmonized System
US DOT - United States Department of Transportation
US OSHA - United States Occupational Safety and Health Administration
vPvB - very Persistent and very Bioaccumulative substances
WHMIS - Workplace Hazardous Material Information System

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

SAFETY DATA SHEET

Document ID: B93029-75 Version AF
Revision Date (year/month/day) 2023/12/22
Last Revision Date (year/month/day) 2023/05/17

Section 16 Other information (Continued)

For further information, please contact your local Beckman Coulter, Inc. representative.

WHILE BECKMAN COULTER, INC. BELIEVES THE INFORMATION CONTAINED HEREIN IS VALID AND ACCURATE, BECKMAN COULTER, INC. MAKES NO WARRANTY OR REPRESENTATION AS TO ITS VALIDITY, ACCURACY, OR CURRENCY. BECKMAN COULTER, INC. SHALL NOT BE LIABLE OR OTHERWISE RESPONSIBLE IN ANY WAY FOR USE OF EITHER THIS INFORMATION OR MATERIALS TO WHICH IT APPLIES. DISPOSAL OF HAZARDOUS MATERIALS MAY BE SUBJECT TO LOCAL LAWS OR REGULATIONS.



SAFETY DATA SHEET

Document ID: B93029-75 Version AF
 Revision Date (year/month/day) 2023/12/22
 Last Revision Date (year/month/day) 2023/05/17

Section 1 Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

Product name Bags 2 - 7
Part number Component of P/N B93029

1.2 Relevant identified uses of the substance or mixture and uses advised against

Product use For Professional/Industrial Use. See product literature for details.

1.3 Details of the supplier of the safety data sheet

Manufacturer

Manufactured for
 Beckman Coulter, Inc.
 250 S. Kraemer Blvd
 Brea, CA 92821, U.S.A.
 Tel: 800-854-3633

Supplier

CANADA
 Beckman Coulter Canada LP
 7075 Financial Drive
 Mississauga, ON L5N 6V8
 Canada
 1-800-463-7828

UNITED KINGDOM
 Beckman Coulter (UK) Ltd.
 Oakley Court
 Kingsmead Business Park, London
 Road
 High Wycombe
 United Kingdom HP11 1JU
 01494 441181

AUSTRALIA
 Beckman Coulter Australia Pty Ltd
 23-27 Chaplin Drive
 Lane Cove NSW 2066
 Australia
 ABN 81 002 011 672
 24 Hour emergency contact phone
 number:
 1800 060 881

SWITZERLAND
 Beckman Coulter Eurocenter SA
 22, rue Juste-Olivier, Case Postale
 1044,
 CH-1260 Nyon 1, Switzerland.
 Telephone: +41 (0)22 365 36 11
 Monday through Friday, 9:00 am to
 7:00pm)

Section 1 Identification of the substance/mixture and of the company/undertaking (Continued)

NEW ZEALAND
Beckman Coulter NZ
Unit J, 33 Walmsley Road, Otahuhu,
Auckland 1062, New Zealand
Hours available: 08:30 - 17:00

ICELAND / ÍSLAND
Beckman Coulter AB
Ekbacksvägen 28
168 69 Bromma
Sweden
Phone No.: +46 80564 85 900
Hours available: 08.00-16.30

MALTA
DX Distributor:
Cherubino Ltd
DELFI Building, Sliema Road, Gzira,
GZR 1637
Telephone: +356 21343270
Hours available: 08:30 – 17:00
SDSNT@beckman.com

e-mail address

1.4 Emergency telephone number

Telephone number (24H)

Chemtrec Emergency Tel No. U.S.A. 800-424-9300, International (001) 703-527-3887

Distributor and emergency phone no.

Refer to attached list, Document ID: [472050](#), for local distributor and emergency phone numbers.

UNITED STATES - Emergency Phone (24h): Chemtrec (800) 424-9300, International (001) 703-527-3887

CANADA - Poison Centre: 1-844-764-7669; Centre antipoison du Québec: 1-800-463-5060

UNITED KINGDOM - For UK and Scotland: Emergency Call 999

IRELAND - National Poisons Information Centre Phone No.: Members of Public: +353 (01) 809 2166 (8:00 am to 10:00 pm 7 days a week); Healthcare Professionals: +353 (01) 809 2566 (24 hour service)

AUSTRALIA - 24 Hour emergency contact phone number: 1800 060 881

NEW ZEALAND - 24 Hour emergency number: 0800 446 109

Section 2 Hazards identification

2.1 Classification of the substance or mixture

Product description Mixture
Reflective Mylar; Bag containing liquid; Odorless

Classification according to EC 1272/2008 (CLP/GHS)
Skin Sensitization Category 1, H317

Classification according to US-OSHA (HCS 29 CFR 1910.1200) and UN GHS
Skin Irritation Category 3
Aquatic Hazard Acute, Category 3

2.2 Label elements

According to EC 1272/2008 (CLP/GHS), US-OSHA and UN GHS
Hazardous ingredients

2-Methyl-4-isothiazolin-3-one
4-Morpholinepropanesulfonic Acid

Pictogram



Signal word

WARNING

Hazard statements

H316 Causes mild skin irritation.
H317 May cause an allergic skin reaction.
H402 Harmful to aquatic life.

Precautionary statements

Prevention

P261 Avoid breathing vapours.
P272 Contaminated work clothing should not be allowed out of the workplace.
P273 Avoid release to the environment.
P280 Wear protective gloves, protective clothing and eye/face protection.

Response

P302+P352 IF ON SKIN: Wash with plenty of soap and water.
P332+P313 If skin irritation occurs: Get medical advice/attention.
P333+P313 If skin irritation or rash occurs: Get medical advice/attention.
P362+P364 Take off contaminated clothing and wash it before use.

Storage

None

Disposal

P501 Dispose of contents/container in accordance with local/national regulations

2.3 Other hazards

Results of PBT and vPvB assessment

PBT: Not applicable.

vPvB: Not applicable.

See Section 11 Toxicological Information for more detailed health information.

SAFETY DATA SHEET

Document ID: B93029-75 Version AF
Revision Date (year/month/day) 2023/12/22
Last Revision Date (year/month/day) 2023/05/17

Section 3 Composition and information on ingredients

3.2 Mixtures

Hazardous ingredients:		Hazard classification of pure ingredients		
Chemical name	% by wt.	EU 1272/2008 CLP/GHS	GHS	Note
4-Morpholinepropanesulfonic Acid CAS # 1132-61-2 EINECS # 214-478-5 Index # Not available	< 3	Eye Irrit. 2, H319 STOT SE 3, H335 Skin Irrit. 2, H315	Eye Irrit. 2, H319 STOT SE 3, H335 Skin Irrit. 2, H315	
2-Methyl-4-isothiazolin-3-one CAS # 2682-20-4 EINECS # 220-239-6 Index # 613-326-00-9	< 0.1	Acute Tox. Dermal 3, H311 Acute Tox. Inhal. 2, H330 Acute Tox. Oral 3, H301 Aquatic Acute 1, H400 Aquatic Longterm 1, H410 Eye Dam. 1, H318 M-factor Acute=10 M-factor Chronic=1 Skin Corr. 1B, H314 Skin Sens. 1A, H317 EUH071 Specific Concentration Limit (SCL) Skin Sens. 1A H317 >= 0.0015% Acute Toxicity Estimates (ATE) ATE Dermal = 300 mg/kg ATE Inhalation - Vapors = 0.5 mg/L ATE Oral = 232 mg/kg	Acute Tox. Dermal 3, H311 Acute Tox. Inhal. 2, H330 Acute Tox. Oral 3, H301 Aquatic Acute 1, H400 Aquatic Longterm 1, H410 Eye Dam. 1, H318 Skin Corr. 1B, H314 Skin Sens. 1A, H317	

See section 8 for available Occupational exposure limits
See Section 15 for additional regulatory information
See Section 16 for description of hazard class and hazard statements

Section 4 First aid measures

4.1 Description of first aid measures

Inhalation

If product is inhaled, move exposed individual to fresh air. If individual is not breathing, begin artificial respiration by trained personnel and obtain medical attention immediately.

Eye contact

If product enters eyes, rinse eyes gently with water as a precaution.

Skin contact

In case of skin contact, rinse with plenty of water. Remove contaminated clothing and shoes. If pain or irritation occurs, obtain medical advice/attention.

Ingestion

If product is ingested, rinse mouth with water. If irritation or discomfort occurs, obtain medical attention immediately.

SAFETY DATA SHEET

Document ID: B93029-75 Version AF
Revision Date (year/month/day) 2023/12/22
Last Revision Date (year/month/day) 2023/05/17

Section 4 First aid measures (Continued)

4.2 Most important symptoms and effects, both acute and delayed

Causes mild skin irritation.
May cause an allergic skin reaction.
See Section 11 Toxicological Information for more detailed health information.

4.3 Indication of any immediate medical attention and special treatment needed

No further relevant information available. Refer to Section 4.1.

Section 5 Firefighting measures

5.1 Extinguishing media In case of fire use carbon dioxide (CO₂), dry chemical, water spray or foam.
For large fires use extinguishing media suitable for surrounding fire.

5.2 Special hazards arising from the substance or mixture Special fire and explosion hazards

No special hazards determined.

Hazardous combustion products

No combustion products posing significant hazards are expected from this product.

5.3 Advice for firefighters

Protective equipment

Self-contained breathing apparatus is recommended for firefighters in all chemical fire situations.

Additional information

No further relevant information available.

Section 6 Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

Personal precautions

Observe general safety guidelines for protection; avoid eye and skin contact.
Wear protective gloves, protective clothing and eye/face protection.

6.2 Environmental precautions

Contain spill to prevent migration.
Do not allow the undiluted product to enter sewers/surface or ground water.

6.3 Methods and material for containment and cleaning up

Spill and leak procedures

Absorb spilled material with an appropriate inert, non-flammable absorbent and dispose according to local regulations.

6.4 Reference to other sections

Refer sections 8 and 13.

Section 7 Handling and storage

7.1 Precautions for safe handling Use good laboratory procedures; avoid eye and skin contact.

Section 7 Handling and storage (Continued)

7.2 Conditions for safe storage, including any incompatibilities

To maintain product quality, store according to the instructions in the product labeling.

Store away from strong acids, strong bases, strong oxidizers and incompatible materials (section 10).

7.3 Specific end uses

No further relevant information available.

Section 8 Exposure controls and personal protection

8.1 Control parameters

Exposure limits

US OSHA None established

ACGIH None established

ACGIH Biological Exposure Indices (BEI)

None established

DFG MAK

2-Methyl-4-isothiazolin-3-one
CAS # 2682-20-4

0.4 mg/m³ Peak (mixture in ratio 1:3 with CAS 26172-55-4) (inhalable fraction); 0.2 mg/m³ TWA MAK (mixture in ratio 1:3 with CAS 26172-55-4) (inhalable fraction)

Ireland None established

IOELVs

None established

NIOSH None established

China None established

Croatia None established

Japan None established

Sweden (AFS 2015:7 and amendments)

None established

Turkey None established

8.2 Exposure controls

Engineering controls

No special engineering controls are required. Use with good general ventilation.

Eye protection

Safety glasses or chemical goggles should be worn to prevent eye contact.

Refer U.S. OSHA 29 CFR 1910.133, European Standard EN166 or appropriate government standards.

Section 8 Exposure controls and personal protection (Continued)

Skin protection	Wear impervious gloves such as Nitrile or equivalent and protective clothing. Refer to U.S. OSHA 29 CFR 1910.138, European Standard EN 374, EN 14605:2005+A1:2009 or appropriate government standards.
Respiratory protection	Under normal conditions, the use of this product should not require respiratory protection. If overexposure should occur and ventilation is not adequate to maintain airborne concentrations at acceptable levels, the use of respiratory protection should be evaluated by a qualified professional.

Section 9 Physical and chemical properties

9.1 Information on basic physical and chemical properties

Physical state	Bag containing liquid	Density and/or relative density	≈ 1
Color	Reflective Mylar	Solubility	
Odor	Odorless	Water	Not applicable
pH	6.8 - 7.5	Organic	Not applicable
Freezing point	Not determined	Partition coefficient n-octanol/water (log value)	Not determined
Boiling point or initial boiling point and boiling range	Not determined	Auto-ignition temp.	Not applicable
Flash point	Not applicable	Decomposition temperature	Not determined
Flammability	Not applicable	Vapor pressure	Not determined
		Kinematic viscosity	Not determined
Lower and upper explosion limit	Not applicable		
Relative vapor density	Not determined		
Particle characteristics	Not applicable		

9.2 Other information

Information with regard to physical hazard classes

No further relevant information available.

Other safety characteristics

No further relevant information available.

Section 10 Stability and reactivity

10.1 Reactivity	No further relevant information available.
10.2 Chemical stability	The product is stable in accordance with recommended storage conditions.
10.3 Possibility of hazardous reactions	No further relevant information available.
10.4 Conditions to avoid	To maintain product performance keep away from strong acids, strong bases, strong oxidizers. Avoid exposure to heat and direct sunlight.
10.5 Incompatible materials	No further relevant information available.
10.6 Hazardous decomposition products	When stored as labeled, no known hazardous decomposition products are formed during the shelf-life of this product.

Section 11 Toxicological information

11.1 Information on hazard classes

Toxicity data for hazardous ingredients

2-Methyl-4-isothiazolin-3-one
CAS # 2682-20-4

Dermal LD50 Rabbit 200 mg/kg (NLM_HSDDB); Inhalation LC50 Rat 0.11 mg/L 4 h (aerosol)(EU_CLH); Oral LD50 Rat 232 - 249 mg/kg (50% MIT in water)(EU_CLH); Oral LD50 Rat 120 mg/kg (50% MIT in water)(EU_CLH)

Primary routes of exposure Eye contact, ingestion, inhalation, and skin contact.

Acute toxicity Not classified based on available data.

Skin corrosion/irritation Causes mild skin irritation.

Serious eye damage/irritation Not classified based on available data.

Respiratory or skin sensitisation May cause an allergic skin reaction.

Germ cell mutagenicity Not classified based on available data.

Carcinogenicity No ingredients in this product are listed as carcinogens by ACGIH, IARC, NTP, OSHA or 1272/2008 EC regulation.

Reproductive toxicity Not classified based on available data.

Specific target organ toxicity (STOT) – single exposure
Not classified based on available data.

Specific target organ toxicity (STOT) – repeated exposure
Not classified based on available data.

Aspiration hazard Not classified based on available data.

Section 11 Toxicological information (Continued)

11.2 Information on other hazards

Endocrine disrupting properties

This product does not have substance(s) of endocrine disrupting properties for health according to REACH Article 57(f).

Other information

No further relevant information available.

Section 12 Ecological information

12.1 Toxicity

Fresh water species

No information available.

Microtox/organisms

No information available.

Water flea

No information available.

Fresh water algae

No information available.

12.2 Persistence and degradability Not determined for the product.

12.3 Bioaccumulative potential Not determined for the product.

12.4 Mobility in soil Not determined for the product.

12.5 Results of PBT and vPvB assessment

Not determined for the product. PBT: Not applicable, vPvB: Not applicable.

12.6 Endocrine disrupting properties

This product does not have substance(s) of endocrine disrupting properties for environment according to REACH Article 57(f).

12.7 Other adverse effects

No further relevant information available.

Section 13 Disposal considerations

13.1 Waste treatment methods

Product waste disposal

Chemical residues and remains should be routinely handled as special waste. This must be disposed of in compliance with anti-pollution and other laws of the country concerned. To ensure compliance we recommend that you contact the relevant (local) authorities and/or an approved waste-disposal company for information.

Additional information

Suggested European waste catalogue 18 01 07 - chemicals other than those mentioned in 18 01 06. Dispose in accordance with national, state and local waste regulations.

Section 14 Transport information

Transportation of this product is not regulated under ICAO, IATA DGR, IMDG, US DOT, European ADR and RID or Canadian TDG.

- 14.1 **UN/ID number:** Not regulated for transportation
- 14.2 **UN proper shipping name:** Not regulated for transportation
- 14.3 **Transport hazard class(es):** Not regulated for transportation
- 14.4 **Packing group:** Not regulated for transportation
- 14.5 **Environmental hazards:** Not regulated for transportation
- 14.6 **Special precautions for user:** None
- 14.7 **Maritime transport in bulk according to IMO instruments:** Not applicable

Section 15 Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

US Federal and State Regulations

SARA 313 (Section 313, Title III reporting requirements)

No ingredients listed.

CERCLA (The Comprehensive Environmental Response, Compensation, and Liability Act) 40 CFR 302.4

No ingredients listed.

California Proposition 65

Chemical which is known to the State of California to cause cancer

No ingredients listed.

Chemical which is known to the State of California to cause development toxicity

No ingredients listed.

Chemical which is known to the State of California to cause male reproductive toxicity

No ingredients listed.

Chemical which is known to the State of California to cause female reproductive toxicity

No ingredients listed.

Massachusetts Right To Know (RTK) List

No ingredients listed.

New Jersey Dept. of Health Right To Know (RTK) List

Section 15 Regulatory information (Continued)

No ingredients listed.

Pennsylvania Right To Know (RTK) List

No ingredients listed.

EU Regulations

This SDS complies with EC Regulations 1907/2006 (REACH) and amendments.

Water Hazard Class (Germany)

WGK 1, low water endangering

Regulation (EU) 2019/1148 on the marketing and use of explosives precursors - Substances Subject to Suspicious Transactions Reporting

No ingredients listed.

Regulation (EU) 2019/1148 on the marketing and use of explosives precursors - Restricted Explosives Precursors

No ingredients listed.

REACH 1907/2006 EC - Candidate List of Substances of Very High Concern (SVHC)

No ingredients listed.

REACH 1907/2006 EC - Annex XVII – Restrictions on Certain Dangerous Substances

Not applicable.

REACH 1907/2006 EC - Annex XIV - list of substances subject to authorisation

No ingredients listed.

Refer to Section 3

UK Regulations

UK REACH Regulation (as Amended) - List of substances subject to authorisation

Refer to Section 3

Canada

This product meets WHMIS criteria for hazardous materials

China

Catalog of Hazardous Chemicals - Hazardous Chemicals

No ingredients listed.

Inventory - China - Inventory of Existing Chemical Substances (IECSC)

All ingredients are listed or exempted.

Turkey

Turkey-REACH - KKDIK Regulation - Annex 17 – Restrictions

Section 15 Regulatory information (Continued)

No ingredients listed.

International

UN/FAO/Rotterdam Convention - Chemicals Subject to Prior Informed Consent (PIC)

No ingredients listed.

15.2 Chemical safety assessment

A Chemical Safety Assessment has not been carried out.

Some hazardous ingredients listed in Section 15 are below the cutoff limits of 0.1% for carcinogen, mutagen and reproductive toxin and 1% for other health hazards required for reporting in Section 3.

Section 16 Other information

Beckman Coulter safety rating	Flammability: 0 Health: 2 Reactivity with water: 0 Physical contact: 2	Code 0=None 1=Slight 2=Caution 3=Severe
--------------------------------------	---	---

Revision changes Updated sections 1, 2, 3, 4, 8 and 15

Document version and issue/revision date

Revision Date (year/month/day) 2023/12/22
Last Revision Date (year/month/day) 2023/05/17
Document ID: B93029-75
Version: AF

Hazard Classification Procedure

This mixture was classified using the calculation method for human health and environmental hazards. Physical hazards were determined based on the specification.

Description of hazard class and hazard statements from Section 3

Aquatic Acute 1 - Aquatic Hazard Acute, Category 1
Acute Tox. Dermal 3 - Acute Toxicity Dermal, Category 3
Acute Tox. Inhal. 2 - Acute Toxicity Inhalation, Category 2
Acute Tox. Oral 3 - Acute Toxicity Oral, Category 3
Eye Dam. 1 - Eye Damage Category 1
Eye Irrit. 2 - Eye Irritation Category 2
Aquatic Longterm 1 - Aquatic Hazard Long term, Category 1
Skin Corr. 1B - Skin Corrosion Category 1B
Skin Irrit. 2 - Skin Irritation Category 2
Skin Sens. 1A - Skin Sensitization Category 1A
STOT SE 3 - Specific Target Organ Toxicity Single Exposure Category 3
EUH071 - Corrosive to the respiratory tract.
H301 - Toxic if swallowed.
H311 - Toxic in contact with skin.
H314 - Causes severe skin burns and eye damage.
H315 - Causes skin irritation.

Section 16 Other information (Continued)

Abbreviations and acronyms

H317 - May cause an allergic skin reaction.
H318 - Causes serious eye damage.
H319 - Causes serious eye irritation.
H330 - Fatal if inhaled.
H335 - May cause respiratory irritation.
H400 - Very toxic to aquatic life.
H410 - Very toxic to aquatic life with long lasting effects.

ACGIH - American Conference of Governmental Industrial Hygienists (ACGIH)
ADR and RID - European Agreement Concerning the International Carriage of Dangerous Goods by Road and Rail
CLP - Classification, Labeling and Packaging
DFGMAK - Republic Germany's maximum exposure limit
EC50 - Concentration of a substance in an environmental medium expected to produce a certain effect in 50% of test organisms
GHS - Globally Harmonized System of Classification and Labelling of Chemicals (GHS)
HCS - Hazard Communication Standard
IARC - International Agency for Research on Cancer
IATA DGR - International Air Transport Association Dangerous Goods Regulation
ICAO - International Civil Aviation Organization
IDLH - Immediately Dangerous to Life or Health
IMDG - International Maritime Dangerous Goods
IMO - International Maritime Organization
IOELVs - European Unions' Indicative Occupational Exposure Limit Values
LC50 - Concentration of a substance in water causing death (50% of the tested population) to aquatic life
LD50 - Lethal Dose 50%
NIOSH - National Institute for Occupational Safety and Health
NTP - National Toxicology Program
OSHA - Occupational Safety and Health Administration
PBT - Persistent Bioaccumulative and Toxic substances
PEL - Permissible Exposure Limit
SARA - Superfund Amendments and Reauthorization Act
STEL – Short Term Exposure Limit
STLV - Short Term Limit Value
STV - Short Term Value
TDG - Canadian Transportation of Dangerous Goods Regulations
TLV - Threshold Limit Value
TWA – Time Weighted Average
UN GHS - United Nations Globally Harmonized System

SAFETY DATA SHEET

Document ID: B93029-75 Version AF
Revision Date (year/month/day) 2023/12/22
Last Revision Date (year/month/day) 2023/05/17

Section 16 Other information (Continued)

US DOT - United States Department of Transportation

US OSHA - United States Occupational Safety and Health Administration

vPvB - very Persistent and very Bioaccumulative substances

WHMIS - Workplace Hazardous Material Information System

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

For further information, please contact your local Beckman Coulter, Inc. representative.

WHILE BECKMAN COULTER, INC. BELIEVES THE INFORMATION CONTAINED HEREIN IS VALID AND ACCURATE, BECKMAN COULTER, INC. MAKES NO WARRANTY OR REPRESENTATION AS TO ITS VALIDITY, ACCURACY, OR CURRENCY. BECKMAN COULTER, INC. SHALL NOT BE LIABLE OR OTHERWISE RESPONSIBLE IN ANY WAY FOR USE OF EITHER THIS INFORMATION OR MATERIALS TO WHICH IT APPLIES. DISPOSAL OF HAZARDOUS MATERIALS MAY BE SUBJECT TO LOCAL LAWS OR REGULATIONS.