

## SAFETY DATA SHEET

### Water Bath Algae Inhibitor

ISSUE DATE: 12/15/14

REV.: 4.0

REV Date : 02/19/2025

ITEM NO: CP-G4800, CPB-190700

DOCUMENT #:  
SDS 4800-4.0

#### **SECTION 1 - CHEMICAL IDENTIFICATION**

PRODUCT NAME:	Algae Clear
PURPOSE:	Water Bath Algae Inhibitor
CP CATALOG #:	CP-G4800, CPB-190700
SUPPLIER'S NAME:	ChemPure Brand Chemicals
SUPPLIERS ADDRESS:	39103 Warren Road, Westland, MI 48185
EMERGENCY CONTACT:	CHEMTREC 1-800-424-9300 (within the U.S.) +1 703-741-5500 (Worldwide)
OTHER INFORMATION:	734-884-4773, INFO@CHEMPUREBRAND.COM

#### **SECTION 2 - HAZARDS IDENTIFICATION**

##### **2.1 Classification of the substance or mixture**

##### **GHS Classification in accordance with 29 CFR 1910 (OSHA HCS)**

H226	Flammable liquids	Category 3
H302	Acute toxicity, Oral	Category 4
H314	Skin corrosion/irritation	Category 1
H318	serious eye damage	Category 1

##### **2.2 GHS Label elements, including Hazard and Precautionary Statement(s)**

##### **Pictogram**



**Signal word: Danger**

##### **Hazard statement(s)**

H226	Flammable liquid and vapour
H302	Harmful if swallowed
H314	Causes severe skin burns and eye damage.
H318	Causes serious eye damage

##### **Prevention, Response, Storage and Disposal Precautionary Statement(s)**

P210	Keep away from heat/sparks/open flames/hot surfaces. No smoking.
P233	Keep container tightly closed.
P240	Ground/bond container and receiving equipment.
P241	Use explosion-proof electrical/ventilating/lighting/equipment
P242	Use only non-sparking tools.
P243	Take precautionary measures against static discharge
P260	Do not breathe dust or mist.
P264	Wash skin thoroughly after handling
P270	Do not eat, drink or smoke when using this product.



# SAFETY DATA SHEET

## Water Bath Algae Inhibitor

ITEM NO: CP-G4800, CPB-190700

ISSUE DATE: 12/15/14

REV.: 4.0

REV Date : 02/19/2025

DOCUMENT #:  
SDS 4800-4.0

P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
P301 + P312	IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.
P301 + P330 + P331	IF SWALLOWED: rinse mouth. Do NOT induce vomiting.
P303 + P361 + P353	IF ON SKIN (or hair): Remove/ Take off immediately all contaminated clothing. Rinse skin with water/ shower.
P304 + P340	IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
P305 + P351 + P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310	Immediately call a POISON CENTER or doctor/ physician.
P321	Specific treatment (see supplemental first aid instructions on this label).
P330	Rinse mouth.
P363	Wash contaminated clothing before reuse.
P370 + P378	In case of fire: Use dry sand, dry chemical or alcohol-resistant foam for extinction
P403 + P235	Store in a well-ventilated place. Keep cool.
P405	Store locked up.
P501	Dispose of contents/ container to an approved waste disposal plant.

### **SECTION 3 - COMPOSITION/INFORMATION ON INGREDIENTS**

Component	CAS #:	EINECS#:	Composition
Water	7732-18-5	231-791-2	Q.S
Benzalkonium Chloride	68424-85-1	270-325-2	40-60%
Ethanol	64-17-5	200-578-6	5-10%
Alkyl dimethyl amines ( C12-C16 )	68439-70-3	270-414-6	<1%

### **SECTION 4 - FIRST-AID MEASURES**

#### **4.1 Description of first aid measures**

**General advice:** Consult a doctor and show this safety data sheet.

- If inhaled:** Remove to fresh air and monitor breathing. If breathing becomes difficult, give oxygen. If breathing stops, give artificial respiration. Consult a doctor.
- In case of skin contact:** Immediately wash skin with copious amounts of soap and water for at least 15 minutes. Remove contaminated clothing and shoes and wash before reuse. Consult a doctor.
- In case of eye contact:** Flush with copious amounts of water for at least 15 minutes. Consult a doctor.
- If swallowed:** Rinse mouth with water. Do not induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Consult a doctor.

**4.2 Most important symptoms and effects, both acute and delayed:** To the best of our knowledge, the chemical, physical and toxicological properties have not been thoroughly investigated.

**4.3 Indication of immediate medical attention and special treatment needed:** Show this safety data sheet to



# SAFETY DATA SHEET

## Water Bath Algae Inhibitor

ITEM NO: CP-G4800, CPB-190700

ISSUE DATE: 12/15/14

REV.: 4.0

REV Date : 02/19/2025

DOCUMENT #:  
SDS 4800-4.0

the doctor in attendance. Immediate medical attention is required.

### **SECTION 5 - FIRE FIGHTING MEASURES**

#### **5.1 Extinguishing media**

**Suitable extinguishing media:** Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide

**5.2 Special hazards arising from the substance or mixture:** Carbon oxides, nitrogen oxides (NOx), Sulfur oxides, Hydrogen chloride gas, Hydrogen fluoride

**5.3 Precautions for fire-fighters:** Wear suitable protective clothing to prevent contact with skin and eyes and self-contained breathing apparatus.

**5.4 Further information:** no data available

### **SECTION 6 - ACCIDENTAL RELEASE MEASURES**

**6.1 Personal precautions, protective equipment and emergency procedures:** Do not take action without suitable protective clothing - see section 8 of SDS. Evacuate personnel to safe areas. Ensure adequate ventilation. Avoid breathing vapors, mist, dust or gas.

**6.2 Environmental precautions:** Do not let product enter drains.

**6.3 Methods and materials for containment and cleaning up:** Cover spillage with suitable absorbent material. Using non-spark tools, sweep up material and place in an appropriate container. Decontaminate spill site with 10% caustic solution and ventilate area until after disposal is complete. Hold all material for appropriate disposal as described under section 13 of SDS.

**6.4 Reference to other sections:** For required PPE see section 8. For disposal see section 13.

### **SECTION 7 - HANDLING AND STORAGE**

#### **7.1 Precautions for safe handling:**

Use in a chemical fume hood, with air supplied by an independent system. Avoid inhalation, contact with eyes, skin and clothing. Avoid the formation of dust and aerosols. Use in a well-ventilated area. Keep away from sources of ignition. Avoid prolonged or repeated exposure.

**7.2 Conditions for safe storage, including any incompatibilities:** Store in cool, well-ventilated area. Keep away from direct sunlight. Keep container tightly sealed until ready for use. . Moisture sensitive. Keep in dry place.

**7.3 Specific end use(s):** Use in a laboratory fume hood where possible. Refer to employer's COSHH risk assessment.

### **SECTION 8 - EXPOSURE CONTROLS/PERSONAL PROTECTION**

#### **8.1 Control parameters: Components with workplace control parameters:**

Component	ACGIH TLV	OSHA PEL	NIOSH IDLH
Ethyl alcohol	STEL: 1000 ppm	(Vacated) TWA: 1000 ppm (Vacated) TWA: 1900 mg/m <sup>3</sup> TWA: 1000 ppm TWA: 1900 mg/m <sup>3</sup>	IDLH: 3300 ppm TWA: 1000 ppm TWA: 1900 mg/m <sup>3</sup>
Component	Quebec	Mexico OEL (TWA)	Ontario TWA/EL



# SAFETY DATA SHEET

## Water Bath Algae Inhibitor

ITEM NO: CP-G4800, CPB-190700

ISSUE DATE: 12/15/14

REV.: 4.0

REV Date : 02/19/2025

DOCUMENT #:  
SDS 4800-4.0

Ethyl alcohol	TWA: 1000 ppm TWA: 1880 mg/m <sup>3</sup>	TWA: 1000 ppm TWA: 1900 mg/m <sup>3</sup>	STEL: 1000 ppm
---------------	--	--	----------------

- ❖ **ACGIH:** American Conference of Governmental Industrial Hygienists.
- ❖ **OSHA:** Occupational Safety and Health Administration.
- ❖ **NIOSH IDLH:** The National Institute for Occupational Safety and Health Immediately Dangerous to Life or Health.

### 8.2 Exposure Controls

- i. **Appropriate engineering controls:** Use in a fume hood where applicable. Ensure all engineering measures described under section 7 of SDS are in place. Ensure laboratory is equipped with a safety shower and eye wash station.

### 8.3 Personal protective equipment

- i. **Eye/face protection:** Face shield and safety glasses Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).
- ii. **Skin protection:** Wear appropriate protective eyeglasses or chemical safety goggles as described by OSHA's eye and face protection regulations in 29 CFR 1910.133 or European Standard EN166
- iii. **Body Protection:** Wear appropriate protective clothing. Complete suit protecting against chemicals, The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.
- iv. **Respiratory Protection:** Follow the OSHA respirator regulations found in 29 CFR 1910.134 or European Standard EN 149. Use a NIOSH/MSHA or European Standard EN149 approved respirator if exposure limits are exceeded or if irritation or other symptoms are experienced.
- v. **Control of environmental exposure:** Do not let product enter drains

## SECTION 9 - PHYSICAL AND CHEMICAL PROPERTIES

### Information on basic physical and chemical properties

<b>Appearance</b>	Liquid	<b>Vapor Pressure</b>	59 hPa at 20 °C
<b>Odor</b>	No Data Available	<b>Vapor Density</b>	No Data Available
<b>Odor Threshold</b>	No Data Available	<b>Relative Density</b>	0.960 g/cm <sup>3</sup>
<b>pH</b>	7.5 at 100 g/l	<b>Water Solubility</b>	No Data Available
<b>Melting / Freezing Point</b>	No Data Available	<b>Partition Coefficient</b>	No Data Available
<b>Initial Boiling Point Range</b>	78 - 100 °C	<b>Auto-Ignition Temperature</b>	No Data Available
<b>Flash Point</b>	40.5 °C - Closed Cup	<b>Decomposition Temperature</b>	No Data Available
<b>Evaporation Rate</b>	No Data Available	<b>Viscosity</b>	No Data Available
<b>Flammability (Solid, Gas)</b>	No Data Available	<b>Explosive Properties</b>	No Data Available
<b>Upper / Lower Flammability Or Explosive Limits</b>	No Data Available	<b>Oxidizing Properties</b>	No Data Available

## SECTION 10 - STABILITY AND REACTIVITY

**Stability:** Stable under recommended transport or storage conditions.

**Conditions to Avoid:** Heat, moisture air or water, incompatibles

**Incompatible Materials:** Strong oxidizing agents



# SAFETY DATA SHEET

## Water Bath Algae Inhibitor

ITEM NO: CP-G4800, CPB-190700

ISSUE DATE: 12/15/14

REV.: 4.0

REV Date : 02/19/2025

DOCUMENT #:  
SDS 4800-4.0

### SECTION 11 - TOXICOLOGICAL INFORMATION

#### ACUTE TOXICITY

Component	LD50 Oral	LD50Dermal	LC50 inhalation
Ethanol	3450 mg/kg ( Mouse )	Not Listed	20000 ppm/10H ( Rat )

**Carcinogenicity:** The table below indicates whether each agency has listed any ingredient as a carcinogen

Component	CAS Number	IARC	NTP	ACGIH	OSHA
Ethanol	64-17-5	Group 1	Known	A3	Listed

### SECTION 12 - ECOLOGICAL INFORMATION

**Eco Toxicity:** Do not empty into drains.

Component	Freshwater Algae	Freshwater Fish	Water Flea	Microtox
Ethanol	EC50 (72h) = 275 mg/l (Chlorella vulgaris)	Fathead minnow (Pimephales promelas) LC50 = 14200 mg/l/96h	EC50 = 9268mg/L/48h EC50 = 10800 mg/L/24h	Photobacterium phosphoreum:EC50 = 34634 mg/L/30 min Photobacterium phosphoreum:EC50 = 35470 mg/L/5 min

**Persistence and degradability:** No data available

**Bioaccumulative potential/ Mobility in soil:** No data available

**Results of PBT and vPvB assessment:** No data available

**Other adverse effects:** May be harmful to the aquatic environment.

**Chronic Toxicity:** There are no known carcinogenic chemicals in this product

### SECTION 13 - DISPOSAL CONSIDERATIONS

Dissolve or mix the material with a combustible solvent and burn in a chemical incinerator equipped with an afterburner and scrubber. Dispose of container and unused contents in accordance with federal, state and local requirements.

### SECTION 14 - TRANSPORT INFORMATION

	UN#	Proper Shipping Name	Hazard Class	Packaging Group
<b>DOT</b>	2920	Corrosive Liquids, Flammable, N.O.S. (Quaternary Ammonium Compounds, Benzyl-C12-16-Alkyldimethyl, Chlorides, Ethanol)	8 (3)	II
<b>TDG</b>	2920	Corrosive Liquids, Flammable, N.O.S. (Quaternary Ammonium Compounds, Benzyl-C12-16-Alkyldimethyl, Chlorides, Ethanol)	8 (3)	II
<b>IATA</b>	2920	Corrosive Liquids, Flammable, N.O.S. (Quaternary Ammonium Compounds, Benzyl-C12-16-Alkyldimethyl, Chlorides, Ethanol)	8 (3)	II
<b>IMDG/IMO</b>	2920	Corrosive Liquids, Flammable, N.O.S. (Quaternary Ammonium Compounds, Benzyl-C12-16-Alkyldimethyl, Chlorides, Ethanol)	8 (3)	II



# SAFETY DATA SHEET

## Water Bath Algae Inhibitor

ITEM NO: CP-G4800, CPB-190700

ISSUE DATE: 12/15/14

REV.: 4.0

REV Date : 02/19/2025

DOCUMENT #:  
SDS 4800-4.0**SECTION 15 - REGULATORY INFORMATION****USA FEDERAL REGULATION  
SARA 311/312 HAZARDOUS CATEGORIZATION**

Acute Health Hazard	Chronic Health Hazard	Fire Hazard	Sudden Release Of Pressure Hazard	Reactive Hazard
Yes	Yes	Yes	No	No

**TSCA 12(b):** Not Applicable**SARA 313:** Not Applicable**Clean Water Act/ Clean Air Act:** Not Applicable**OSHA:** Not Applicable**California Proposition 65:** This product does not contain any Proposition 65 chemicals.**CERCLA:** Not Applicable**STATE RIGHT TO KNOW**

Massachusetts	New Jersey	Pennsylvania	Illinois	Rhode Island
Listed	Listed	Listed	-	-

**HMIS Rating**

Health hazard	Chronic Health Hazard	Flammability	Physical Hazard
3	*	2	0

**NFPA Rating**

Health Hazard	Fire Hazard	Reactivity Hazard
3	2	0

**US DEPARTMENT OF TRANSPORTATION**

REPORTABLE QUANTITY (RQ)	DOT MARINE POLLUTANT	DOT SEVER MARINE POLLUTANT
No	No	No

**U.S. DEPARTMENT OF HOMELAND SECURITY:** This product does not contain any DHS chemicals.**Canada:** This product has been classified in accordance with the hazard criteria of the Controlled Products Regulations (CPR) and the SDS contains all the information required by the CPR.**This SDS complies with the requirements of Regulation (EC).****SECTION 16 - OTHER INFORMATION**

The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. ChemPure Brand Chemicals shall not be held liable for any damage resulting from handling or from contact with the above product. Refer to the website for additional terms and conditions of sale. This product is sold for laboratory research and development purposes use only.