



CERTIFICATE OF CONFORMANCE

Part No: B29894

Description: Cap-Closure, PPS Assy, Set of 6

P/N(s) B29894 – The component(s) / specific lot(s) used in the assembly conform to Beckman print and process control specifications and are manufactured under ISO 9001:2015 guidelines.

B29894:

Contains B16370, cap/closure assembly; which contains B16369, and uses a Glass-Filled Polyphenylsulfide material which meets the requirements of ISO 10993 and USP Class VI, and a stainless-steel wire handle 363656.

A handwritten signature in black ink, appearing to read 'William Swanson'.

William Swanson

Senior Quality Assurance Engineer