



October 27, 2020

## CERTIFICATE OF CONFORMANCE

Part No: C15124  
Lot No: P27705  
Date of Mfr: 10/21/2020  
Description: ASSEMBLY, 1000ML BOTTLES, PP+PPS, SET OF 2

P/N(s) C15124– The component(s) / specific lot(s) used in the assembly conform to Beckman print and process control specifications and are manufactured under ISO 9001 guidelines.

C15124:

Contains B16370, cap/closure assembly; which contains B16369, and uses a Glass-Filled Polyphenylsulfide material which meets the requirements of ISO 10993 and USP Class VI, and a stainless-steel wire handle 363656.

Contains C15129, autovent plug assembly; which contains 393744 autovent plug, polypropylene copolymer with glass fiber reinforcement, which meets the requirements of USP Class VI. Also contained

within this assembly are A98816, O-ring and C11932, Quad Ring, which meet the requirements of USP Class VI.

Contains 393895, bottle, 1000mL, polypropylene; which complies with the U.S. Pharmacopeia criteria for Class VI.

A handwritten signature in black ink, appearing to read 'William Swanson'.

William Swanson

Senior Quality Assurance Engineer