

## CERTIFICATE OF CONFORMANCE

Part No: 356989  
Lot No: P36306  
Date of Mfr: 1/20/2022  
Description: Bottle, with Screw-On Cap, 230mL, Polypropylene (PP), Conical, set of 6

Part No. 356989 – The component(s) / specific lot(s) used in the assembly conform to Beckman print and process control specifications and are manufactured under ISO 9001 guidelines.

The 230mL conical PP bottle with screw-on cap set contains p/n 344691, cap, bottle, 38mm, which is manufactured from a polypropylene material that complies with the U.S. Pharmacopoeia criteria for a Class VI plastics and the FDA Code of Federal Regulations, Title 21 §177.1520 (c)3.1a; p/n 356988, bottle, rotor, 230mL, conical, which is manufactured using a polypropylene material that complies with the U.S. Pharmacopoeia criteria for a Class VI plastics and the FDA Code of Federal Regulations, Title 21 §177.1520 (c) 3.1a & 3.2a.



William Swanson

Senior Quality Assurance Engineer