



May 25, 2022

CERTIFICATE OF CONFORMANCE

Part No: 355605
Lot No: P38460
Date of Mfr.: 5/16/2022
Description: Bottle Assembly, Polycarbonate, 500 mL, 69 x 160 mm (qty. 6)

P/N(s) 355605 – The component(s) / specific lot(s) used in the assembly, conform to Beckman print and process control specifications and are manufactured under ISO 9001 guidelines. The 500mL PC bottle assembly contains p/n 334419, 870411, 344692 and 339245.

P/N 334419, stopper, bottle; which uses a material that complies with the FDA Code of Federal Regulations, Title 21 § 177.1520. Also contained within is p/n 870411, O-ring, which is a general-purpose neoprene.

P/N 344692, cap, bottle, 48mm; which uses a material that complies with the FDA Code of Federal Regulations, Title 21 § 177.1520 (c) 3.1a and the U.S. Pharmacopoeia criteria for Class VI plastics.

P/n 339245, bottle, 500mL polycarbonate; which uses a material that complies with the USP 23 Plastics Class VI and ISO 10993-Part 1 testing requirements.

A handwritten signature in black ink, appearing to read 'William Swanson'.

William Swanson
Senior Quality Assurance Engineer