



July 25, 2023

CERTIFICATE OF CONFORMANCE

Part No: 355618

Lot No: P35875

Date of Mfr.: 7/17/2023

Description: Bottle Assembly, Polycarbonate (PC), 1 x 3-1/2, (Box-6)

P/N(s) 355618 – The component(s) / specific lot(s) used in the assembly conform to Beckman print and process control specifications and are manufactured under ISO 9001 guidelines.

355618:

Contains 340382, bottle, rotor, 1 x 3-1/2, polycarbonate; which uses a material that complies with U.S. Pharmacopeia criteria for Class IV plastics and ISO 10993-1.

Cap: 338824 aluminum cap; which uses an aluminum alloy material.

Plugs The polymer material used to fabricate the components within the above listed assemblies is a general purpose Noryl.

O-Ring: The elastomeric material used to fabricate the components within the above listed assemblies is a general-purpose neoprene material.

A handwritten signature in black ink, appearing to read 'William Swanson'.

William Swanson
Senior Quality Assurance Engineer