

Beckman Coulter
 MIAMI AND KENTUCKY REAGENTS MANUFACTURING
 CERTIFICATE OF ANALYSIS

| | | |
|------------------------|--------------------------------|---------------------------------|
| IDENTIFICATION: | Product Name: | CytoFLEX Daily QC Fluorospheres |
| | Part Number: | B53230 |
| | Lot Number (Finished Product): | 15AQAF |
| | Lot Number (Sub-assembly): | 15AQA |
| | Date of Manufacture: | 20 Sep 2023 |
| | Expiration Date: | 01 Aug 2026 |

| TEST | SPECIFICATION | RESULTS | Pass/Fail |
|---|---|---------------------|-----------|
| CV%: FITC, PE, PER-CP, APC, V450 and V500 | ≤ 6% | Meets Specification | Pass |
| Mean Diameter | 3.0-3.5µm | Meets Specification | Pass |
| Particle Count | 1.27 x 10 ⁷ to 1.73 x 10 ⁷ Particles/mL | Meets Specification | Pass |
| Coefficient of Variation | <3% | Meets Specification | Pass |
| Percent in Peak: FSC | >80% | Meets Specification | Pass |
| | | | |

The product was manufactured in accordance with the US FDA Quality System Regulation (GMP).

| APPROVAL | TITLE | DATE |
|--|---------------------------------|-------------|
| DocuSigned by: <i>Franklin Obando</i> | Sr. Quality Assurance Scientist | 27-Sep-2023 |

Signer Name: Franklin Obando
 Signing Reason: I am the author of this document
 Signing Time: 27-Sep-2023 | 7:57:17 AM EDT

CE676AF86C3949529F9BF8D00E2E94F8

PRDRELCAB53230-AF 08/2022

B51501 Rev AF