

Beckman Coulter
MIAMI AND KENTUCKY REAGENTS MANUFACTURING
CERTIFICATE OF ANALYSIS

IDENTIFICATION:	Product Name:	CytoFLEX Ready to Use Daily QC Fluorospheres
	Part Number:	C65719
	Lot Number (Finished Product):	15AQEF
	Lot Number (Sub-Assembly):	15AQE
	Date of Manufacture:	20-NOV-2023
	Expiration Date:	22-SEP-2025

TEST	SPECIFICATION	RESULTS	Pass/Fail
rCV: FITC, PE, APC, V450, NUV450 and UV405	≤5%	Meets Specification	Pass
Mean Diameter	3.0-3.4 μm	Meets Specification	Pass
Particle Count	1.2 x 10 ⁶ to 1.4 x 10 ⁶ Particles/mL	Meets Specification	Pass
Percent in Peak: FSC	>80%	Meets Specification	Pass

The product was manufactured in accordance with the US FDA Quality System Regulation (GMP).

DocuSigned by: APPROVAL	TITLE	DATE
<i>Esme (Esmeralda) Iglesias-Castillo</i>	QUALITY ASSURANCE SCIENTIST I	04-Dec-2023

Signer Name: Esme (Esmeralda) Iglesias-Castillo
Signing Reason: I am the author of this document
Signing Time: 04-Dec-2023 | 9:29:26 AM EST

B683536A50444E1CBFAF6254E42D42C3