



CERTIFICATE OF CONFORMANCE

Part No: 361690
Lot No: Non-Specific
Description: Bottle Assembly, 500mL, Polycarbonate (PC), set of 6

P/N(s) **361690** – The component(s) / specific lot(s) used in the assembly conform to Beckman print and process control specifications and are manufactured under ISO 9001:2008 guidelines.

361690 Bottle Assembly which contains: **339245** – bottle, 500mL, polycarbonate, which is made of a material that complies with the U.S. Pharmacopeia criteria for Class VI plastics, the FDA Code of Federal Regulations, Title 21 § 177.1580, and FDA modified ISO 10993-1 testing. **358679** – cap, 48mm, which uses a general purpose molding compound. **927860** – O-ring, which is a general purpose silicone.

**Valid for all current lot numbers, specific to these part number(s),
as of the date of issuance of this certificate.**

Issuance date: January 1, 2011