



CERTIFICATE OF CONFORMANCE

Part No: 355624
Lot No: Non-Specific
Description: Bottle & Cap Assembly, 100mL, Polypropylene (PP), set of 6

P/N(s) **355624**– The component(s) / specific lot(s) used in the assembly conform to Beckman print and process control specifications and are manufactured under ISO 9001:2008 guidelines.

355624 Bottle & Cap Assembly which contains: **356951** – bottle, rotor 100mL, which is made of a polypropylene material that complies with the U.S. Pharmacopeia criteria for Class VI plastics and the FDA Code of Federal Regulations, Title 21 § 177.1520 (c) 3.1a & 3.2a. **358632** – cap, bottle, which uses a general purpose molding compound. **889633** – O-ring, which is a general purpose nitrile material.

**Valid for all current lot numbers, specific to these part number(s),
as of the date of issuance of this certificate.**

Issuance date: January 1, 2011