



CERTIFICATE OF CONFORMANCE

Part No: 357003
Lot No: Non-Specific
Description: Bottle, with Screw-On Cap, 50mL, Polypropylene (PP), set of 25

P/N(s) **357003** – The component(s) / specific lot(s) used in the assembly conform to Beckman print and process control specifications and are manufactured under ISO 9001:2008 guidelines.

357003 Bottle with Screw-On Cap which contains: **356287** – bottle, 50mL, polypropylene, which is made of a material that complies with the U.S. Pharmacopeia criteria for Class VI plastics and the FDA Code of Federal Regulations, Title 21 § 177.1520 (c) 3.1a & 3.2a. **356284** – cap, bottle, 24mm, which uses a polypropylene material that complies with the U.S. Pharmacopeia criteria for Class VI plastics and the FDA Code of Federal Regulations, Title 21 §177.1520 (c) 3.1a.

**Valid for all current lot numbers, specific to these part number(s),
as of the date of issuance of this certificate.**

Issuance date: January 1, 2011