



## CERTIFICATE OF CONFORMANCE

Part No: 355618  
Lot No: Non-Specific  
Description: Bottle Assembly, 26.3mL, Polycarbonate (PC), set of 6

P/N(s) **355618** – The component(s) / specific lot(s) used in the assembly conform to Beckman print and process control specifications and are manufactured under ISO 9001:2008 guidelines.

**355618** Bottle Assembly which contains: **335258** – insert, bottle, which is made of a general purpose molding compound. **338824** – cap, bottle, rotor, uses an aluminum alloy material. **340382** – bottle, polycarbonate, uses a material that complies with the U.S. Pharmacopeia criteria for Class VI plastics, the FDA Code of Federal Regulations, Title 21 §177.1580, and FDA modified ISO 10993-1 testing. **870385** – O-ring, which is a general purpose neoprene.

**Valid for all current lot numbers, specific to these part number(s),  
as of the date of issuance of this certificate.**

Issuance date: January 1, 2011