

Compatibility Chart

UNIVERSAL PIPETTE TIPS (BARRIER TIPS)

Genesee Cat#				26-400	26-401	26-404	26-407	26-412	26-415	26-430
Brand	Model	Volume Range	Channels	10µl	10µl Reach	20µl	100µl	200µl	300 µl	1000µl Reach
Brand	Transferpette S	0.1-1µl, 0.1-2.5µl, 0.5-10µl	1	✓	✓					
Brand	Transferpette S	2-20µl, 10-100µl	1							
Brand	Transferpette S	20-200µl	1			✓	✓	✓		
Brand	Transferpette S	100-1000µl	1							✓
Brand	Transferpette S	0.5-10µl	8/12	✓	✓					
Brand	Transferpette S	20-200µl	8/12			✓	✓	✓		
Brand	Transferpette S	30-300µl	8/12				✓	✓	✓	
CAPP	Aero	0.2-2µl, 0.5-10µl	1	✓	✓					
CAPP	Aero	2-20µl	1			✓				
CAPP	Aero	10-100µl	1			✓	✓	✓		
CAPP	Aero	25-200µl	1			✓	✓	✓		
CAPP	Aero	50-300µl	1						✓	
CAPP	Aero	100-1000µl	1							✓
CAPP	Aero	0.5-10µl	8/12	✓	✓					
CAPP	Aero	5-50µl	8/12				✓	✓	✓	
CAPP	Bravo	0.1-2µl, 0.5-10µl	1	✓	✓					
CAPP	Bravo	2-20µl	1			✓				
CAPP	Bravo	5-50µl	1				✓	✓	✓	
CAPP	Bravo	10-100µl	1				✓	✓		
CAPP	Bravo	20-200µl	1				✓	✓	✓	
CAPP	Bravo	100-1000µl	1							✓
CAPP	Bravo	0.5-10µl	8/12	✓	✓					
Eppendorf	Research	0.1-2.5µl, 0.5-10µl	1	✓	✓					
Eppendorf	Research	2-20µl	1			✓				
Eppendorf	Research	10-100µl	1			✓	✓	✓		
Eppendorf	Research	20-200µl	1			✓	✓	✓	✓	
Eppendorf	Research	30-300µl	1			✓	✓	✓	✓	
Eppendorf	Research	100-1000µl	1							✓
Eppendorf	Research	0.5-10µl	8/12	✓	✓					
Eppendorf	Research	10-100µl	8/12			✓	✓	✓		
Eppendorf	Research	30-300µl	8/12			✓	✓	✓	✓	
Eppendorf	Research Plus	0.1-2.5µl, 0.5-10µl	1	✓	✓					
Eppendorf	Research Plus	2-20µl ultra micro	1		✓					
Eppendorf	Research Plus	2-20µl	1			✓				
Eppendorf	Research Plus	10-100µl	1			✓	✓	✓		
Eppendorf	Research Plus	20-200µl	1			✓	✓	✓	✓	
Eppendorf	Research Plus	30-300µl	1			✓	✓	✓	✓	
Eppendorf	Research Plus	100-1000µl	1							✓
Eppendorf	Research Plus	0.5-10µl	8/12	✓	✓					
Eppendorf	Research Plus	20-300µl	8/12			✓	✓	✓	✓	
Eppendorf	Reference	0.1-2.5µl, 0.5-10µl	1	✓	✓					
Eppendorf	Reference	2-20µl	1			✓				
Eppendorf	Reference	10-100µl	1			✓	✓	✓		
Eppendorf	Reference	100-1000µl	1							✓
Eppendorf	Reference 2	0.1-2.5µl, 0.5-10µl	1	✓	✓					
Eppendorf	Reference 2	2-20µl	1			✓	✓	✓		
Eppendorf	Reference 2	10-100µl	1			✓	✓	✓		
Eppendorf	Reference 2	20-200µl	1			✓	✓	✓		
Eppendorf	Reference 2	30-300µl	1			✓	✓	✓	✓	
Eppendorf	Reference 2	100-1000µl	1							✓
Gilson	Pipetman Classic	0.2-2µl, 1-10µl	1	✓	✓					
Gilson	Pipetman Classic	2-20µl	1			✓				
Gilson	Pipetman Classic	10-100µl	1			✓	✓	✓		
Gilson	Pipetman Classic	20-200µl	1			✓	✓	✓	✓	
Gilson	Pipetman Classic	100-1000µl	1							✓
Gilson	Pipetman Classic	2-20µl	8/12			✓				
Gilson	Pipetman Classic	20-200µl	8/12			✓	✓	✓	✓	
Gilson	Pipetman Classic	30-300µl	8/12			✓	✓	✓	✓	
Gilson	Pipetman L/M/G	0.2-2µl, 1-10µl	1	✓	✓					
Gilson	Pipetman L/M/G	2-20µl	1			✓				
Gilson	Pipetman L/M/G	10-100µl	1			✓	✓	✓		
Gilson	Pipetman L/M/G	20-200µl	1			✓	✓	✓	✓	
Gilson	Pipetman L/M/G	100-1000µl	1							✓
Gilson	Pipetman L/M/G	0.5-10µl	8/12	✓	✓					
Gilson	Pipetman L/M/G	2-20µl	8/12			✓				
Gilson	Pipetman L/M/G	20-200µl	8/12			✓	✓	✓	✓	
Gilson	Pipetman L/M/G	30-300µl	8/12			✓	✓	✓	✓	
Gilson	Pipetman Neo	0.2-2µl, 1-10µl	1	✓	✓					
Gilson	Pipetman Neo	2-20µl	1			✓				
Gilson	Pipetman Neo	10-100µl	1			✓	✓	✓		
Gilson	Pipetman Neo	20-200µl	1			✓	✓	✓	✓	
Gilson	Pipetman Neo	100-1000µl	1							✓
Hamilton	Softgrip	0.2-2µl	1	✓	✓					
Hamilton	Softgrip	1-10µl	1		✓					

[illegible]

[illegible]